

Orthotic Essentials Product Catalogue

As a key NHS partner, we've grown into a leading UK manufacturer of orthotic products, uniting bespoke and standard devices with innovation and traditional skills.

Contents

Introduction	4
Head and Neck	10
Spine and Torso	26
Arm and Shoulder	46
Wrist and Hand	54
Hip and Knee	74
Ankle and Foot	94
Adult Insoles	174
Paediatric	208
Index	256

Welcome to the Steeper Orthotic Essentials Product Catalogue

Our Orthotic Essentials range offers a selection of carefully handpicked stock products that enhance and improve patient mobility. A full head-to-toe range, these high-quality, off-the-shelf options offer clinicians value for money and are readily available within a short lead time.

For any queries regarding the products within this catalogue, or any others outside of this range, please contact orthotics@steepergroup.com.

Icons key

Throughout this catalogue you'll find icons to help you better understand each product. Here's what they mean:



This product is available to order to your patient's custom measurements.



This product can be ordered to arrive pre-assembled. Please specify when ordering.



This product isn't kept in stock but it is available to order. Delivery times may vary.



This product is sold in a pair.



This product is available in a universal size that fits most patients.



Click this button to bring you back to the original contents page.

Bespoke custom manufacturing service

To complement our comprehensive stock products range, we also manufacture a variety of fully bespoke products at our Steeper manufacturing sites in the UK. These products include:



Creating custom devices for orthotic application demands rare skills with tried-and-tested techniques. We proudly employ expert technicians with these skills, who are committed to crafting the perfect orthoses to provide superior stabilisation, support and comfort to patients.

With specialisms from paediatric all the way to geriatric, our technical teams produce products that provide exceptional fit and functionality.

For further information on our bespoke custom products and manufacturing capabilities, please contact orthotics@steepergroup.com.

ARM & SHOULDER


Ordering Information

How to order

Steeper accepts orders by the following methods:

- **Phone** 0113 207 0435
- **Fax** 0113 272 5689
- Email order@steepergroup.com

Please specify product code, quantity required and any special delivery information.



The following options are available as standard:

- First class post
- Next working day

If you require a timed delivery the following options are available:

- Next working day before 9:00*
- Next working day before 10:30
- Next working day before 12:00 noon
- Saturday delivery**
- Same day delivery**

For next day delivery, orders must be placed with Steeper by 16:00.

Additional delivery charges may apply for timed delivery, please contact Customer Services for further details.

*Subject to location, restrictions may apply.

**Additional charges and restrictions may apply, please contact Customer Services for further details.

We can offer order confirmation on all orders via e-mail or fax.

Delivery shortages/satisfaction

On receipt of your goods, please check that the delivery note matches what you have received. In the unlikely event of receiving the wrong product or quantity, please notify our Customer Services team within 7 days.



Returns procedure

All stock may be returned within a 30 day period unless otherwise specified. It must be in its original packaging and in a saleable condition. Please note, soiled/contaminated goods must not be returned.

To return a product, you must complete a Returns Form obtainable from our Customer Service team.

We do not charge restocking fees for our range of orthotic essentials products. Custom and made to order products are non-returnable as these have been made outside the standard size ranges. If in doubt please contact our Customer Services team for clarity on 0113 207 0435.

Faulty/damaged goods

In the unlikely event that your order is received in a damaged or faulty condition, we kindly request you contact our Customer Service department. You will then be required to complete a Complaints Form, which should be submitted along with the damaged product for further assessment and resolution.

Single patient use

Any product supplied by Steeper is for single patient use only and should be fitted by a suitably qualified professional.

Re-use of any products is potentially dangerous and will invalidate the CE mark or any warranty.

ADULT INSOLES


Meet Our Model Jacob Chambers


Jacob was born with **Cerebral Palsy Hemiplegia** affecting the right side of his body and has spent many hours with clinicians and other medical professionals throughout his life, much like many of our very own patients. Despite the challenges that come with such a condition, plus the addition of ADHD, he has never let his disability hold him back! A bright and energetic young man, he channels his energy into his love of sports including international-level competitive swimming.

A professional and experienced model, Jacob has worked with many industry leaders such as global image agency, Getty Images, and featured in the Paralympics Advert in 2021.

Alongside modelling, Jacob also loves acting and has attended the Manchester School of Acting every Saturday for many years, honing his craft. He played the part of a brother, Sammy, in Creeped Out 2 and was cast as the main character, Lemar, in a short film called Liccie Bit.

Collaborating with Jacob for our Orthotics Essentials Catalogue just made sense. He resembles a lot of our patients in their resilience and determination to not let their disability stand in their way.





PAEDIATRIC

HEAD & NECK

Starlight Standard

Features and Benefits

- Impact absorbing foam construction
- · Cotton jersey lining
- · Soft to the skin
- Wicks away moisture
- · Good moisture circulation
- Outer layer made from robust cotton for longevity
- Lightweight approximately 100g
- Machine washable
- Available in head circumferences 45cm to 62cm in 1cm increments

Indications for Use

 Simple yet effective protection for people at risk of head injury

Styles Available



Note: This product is non-returnable for hygiene reasons.

*Starlight Standard in Brown is not kept in stock, delivery times may vary.

Ordering Information

Colour	Product Code
Brown	5180-
Red	5180-1-
Blue	5180-2-

= Circumference of head (45cm to 62cm in 1cm increments)

Example: For a red Starlight Standard Helmet measuring 50cm - 5180-1-50





Starlight Quilted Soft


Features and Benefits

- The Starlight Quilted Soft is upholstered with highly efficient high-tech foam
- The foam has excellent shock absorbing properties
- The forces applied by a fall are redistributed over a larger area
- Adapts well to the shape of the head offering improved comfort
- Available in head circumferences 45cm to 62cm in 1cm increments
- Excellent moisture wicking properties

Indications for Use

 Advanced protection for people at risk of head injury

Styles Available



Note: This product is non-returnable for hygiene reasons.

Ordering Information

Colour	Product Code
Oxford Rose	5180SG-1-
Oxford Light Blue	5180SG-2-

= Circumference of head (45cm to 62cm in 1cm increments)

Example: For an Oxford Rose Starlight Quilted Soft Helmet measuring 50cm - 5180SG-1-50



PAEDIATRIC

Starlight Secure Fabric



Features and Benefits

- Impact absorbing soft/firm foam combination construction
- · High quality cotton
- · Soft to the skin
- Lightweight approximately 100g 130g (size dependent)
- Easy to wash

Indications for Use

 Advanced protection for people at risk of head injury

Note: Starlight Secure model shown with optional forehead protector, neck cushion and chin guard.

Note: This product is non-returnable for hygiene reasons.

Features and Benefits

- Impact absorbing soft/firm foam combination construction
- · Leather lined
- Soft to the skin
- Lightweight approximately 150g 240g (size dependent)

Starlight Secure Leather

- · Heat and moisture regulating
- Breathable design
- Easy to clean

Indications for Use

 Advanced protection for people at risk of head injury

Note: This product is non-returnable for hygiene reasons.

Ordering Information

Size	Head Circ. (cm)	Product Code	Accessories	Product Code
1	47 - 49	C401101B-1	Chin Guard	C4012B-1
2	50 - 52	C401102B-1	Forehead Protector	C4013B-1
3	53 - 55	C401103B-1	Neck Cushion	C4014B-1
4	56 - 59	C401104B-1		
5	60 - 62	C401105B-1		

Is grey not your colour?

Made to Measure service is available for this product!

Other colours are available when ordering made to measure.





Ordering Information

Size	Head Circ. (cm)	Product Code	Accessories	Product Code
1	47 - 49	C401101	Chin Guard	C4012
2	50 - 52	C401102	Forehead Protector	C4013
3	53 - 55	C401103	Neck Cushion	C4014
4	56 - 59	C401104		
5	60 - 62	C401105		

Is beige not your colour?

Made to Measure service is available for this product!

Other colours are available when ordering made to measure.







Basic Foam Collar

Features and Benefits

- Symmetrical design provides comfortable support
- Suitable for night use
- Touch and close fastening
- Made from medium density foam covered with stockinette

Indications for Use

- Spondylitis
- Spondylosis
- · Rheumatoid arthritis
- Osteoarthritis

Note: Made to measure service available, please provide neck circumference and front depth.

Ordering Information

Size	Neck Circ. (cm)	Front Depth (cm)	Product Code
X-Small	33 - 37	7	BFC/XS
Small	38 - 42	8	BFC/S
Medium	43 - 46	9	BFC/M
Large	47 - 51	10	BFC/L
X-Large	52 - 56	11	BFC/XL

Need a different size to fit your patient?

Made to Measure service is available for this product!

To order, please provide neck circumference and front depth.





More InformationWear and Care Instructions

Shaped Foam Collar

Features and Benefits

- Symmetrical design provides comfortable support, suitable for night use
- Increased shaping to reflect underlying anatomy and provide an improved fit
- Touch and close fastening
- Made from medium density foam covered with stockinette

Indications for Use

- Spondylitis
- Spondylosis
- · Rheumatoid arthritis
- Osteoarthritis

Note: Made to measure service available, please provide neck circumference and front depth.

Ordering Information

Size	Neck Circ. (cm)	Front Depth (cm)	Product Code
Small	38 - 42	6.5	SFC/S
Medium	43 - 46	7	SFC/M
Large	47 - 51	8	SFC/L

Need a different size to fit your patient?

Made to Measure service is available for this product!

To order, please provide neck circumference and front depth.



HEAD & NECK



Heads First Collar

Features and Benefits

- Open design makes it breathable and helps keep the area cool
- Attractive, minimal appearance
- Quick and easy fitting
- · Comfort for daily use and increased patient compliance
- Access to a tracheal tube
- Supplied pre-assembled and ready to
- · Can be bent to the desired fit by hand
- Constructed from surgical rubber tubing, around a malleable core. The tubing is covered with a brushed nylon material that allows the hook and loop attachments to adhere anywhere

Indications for Use

Originally designed for people suffering from arthritis, but has been found to be appropriate and of great benefit to people suffering from many other disabilities including:

- · Progressive neurological conditions, such as motor neurone disease
- · Congenital conditions such as cerebral palsy
- Whiplash
- · Head and spinal cord injury
- Burns

Note: For an accurate fit, measure from the chin to the corner of the jaw.

Ordering Information

Size	Jaw to Chin Length (cm)	Product Code
Junior	9	HFC4
Small	10	HFC3
Medium	11.5	HFC2
Large	13	HFC1

Accessories	Product Code		
Accessories	Large/Medium	Small/Junior	
Front Support	FS01	FS02	
Chin Pad	CP1	CP2	
Neck Strap	NS1	NS2	

Universal PROGLIDE 172



Features and Benefits

The PROGLIDE 172 cervical collar's easyglide system has 12 adjustments available, allowing for improved fit and aims to be a one-size-fits-all solution. It maintains posture and with improved padding, allows more effective pressure distribution to decrease the risk of redness and skin breakdown.

- Ergonomically adjustable control system
- 12 individual adjustments
- · Pulley system cervical collar
- · Soft cloth inner liner
- Tracheal opening
- MRI, CT and X-Ray compatible
- Effective between C2 C5

Indications for Use

- Cervical spine precaution for trauma patients
- Immobilisation for pre and post cervical spine surgery
- Degenerative disorders
- Spinal stenosis
- Spondylolisthesis
- Anterior cervical fusion
- Halo removal
- Dens type 1 cervical fractures of C2 vertebrate
- Anterior discectomy
- Teardrop fracture of the vertebral body
- Cervical strain
- Neural involvement
- Atlanto-occipital dislocation


Description	Size	Product Code	
172 Collar	Universal	PG172	
Spare Liner Set	Universal	PGPAD172	

HEAD & NECK

Universal PROGLIDE 174

Features and Benefits

The PROGLIDE 174 cervical collar is designed to offer maximum support with the addition of an extended sternal pad for increased occipital support. Featuring the same easy-glide system as the PROGLIDE 172. it aims to be a one-size-fits-all solution.



- Ergonomic adjustable control system
- 12 individual adjustments
- Pulley system cervical collar
- · Soft cloth inner liner
- Tracheal opening
- MRI, CT and X-Ray compatible

Indications for Use

- · Cervical spine precaution for trauma patients
- Immobilisation for pre and post cervical spine surgery
- Degenerative disorders
- Spinal stenosis
- Spondylolisthesis
- Anterior cervical fusion
- Halo removal
- Dens type 1 cervical fractures of C2 vertebrate
- · Anterior discectomy
- Teardrop fracture of the vertebral body
- Cervical strain
- Neural involvement
- Atlanto-occipital dislocation

Ordering Information

Description	Size	Product Code
174 Collar	Universal	PG174
Spare Liner Set	Universal	PGPAD174



The Miami J Select is designed to support successful patient outcomes through its wide range of height adjustability.

Featuring a patented, first-tomarket locking mechanism to protect against tampering and intuitive blue Patient Touch-Points, this collar encourages patient compliance and confidence.

The anti-microbial Sorbatex padding prevents bacteria buildup while being easy to clean and replace when needed.

Miami J® Select Collar

Features and Benefits

- · Easy height adjustability with a tamperproof locking mechanism - 'Left to Lock'
- Anti-microbial and biocompatible Sorbatex[™] padding prevents bacterial build-up
- Universal design
- 7 precise fitting sizes with opportunity to fine-tune
- Miami J[®] inspired chin tray
- One-piece, anatomically inspired collar back with Flex Edge® designed to reduce occipital pressure
- Easy to don and doff with patented Reproducible Fit option
- Intuitive blue Patient Touch-Points to send patients home with confidence -Blue is You™
- Simplified comfort adjustments with integrated and patented Sternal-Relief Dial
- Large tracheal opening improved function for ICU patients
- Easy to clean and replace multi-coloured
- · X-ray and CT lucent, MR-safe
- Compatible with Occian™ Back & Miami JTO® accessories

Indications for Use

- Pre and post c-spine surgery
- C-spine precaution for trauma patients

Description	Product Code
Collar Set	MJSR-101
Miami J Select Collar	MJS-101

Accessories	Product Code
Miami J Select Replacement Pads	MJSP-101
Miami J Select Tabs Reproducible Fit Tabs (10 pack)	MJST-101

HEAD & NECK



Miami J[®] Cervical Collar

Features and Benefits

- Sorbatex[™] padding is antibacterial and clinically proven to inhibit microbial growth, which enhances comfort and helps ensure patient compliance
- MRI, CT and X-ray lucent
- · Bioengineered to minimise pressure points in key known areas: chin, occiput, trapezius and clavicle
- Comfort adjustment buttons for occiput and sternal pad removable when lying and swallowing

Indications for Use

- Post trauma
- Stable cervical fractures
- Cervical spondylitis
- · Herniated cervical disc
- RA and OA cervical spine
- Cervical spondylosis
- Post operative

Ordering Information

Size	Product Code			
3126	Collar Set*	Collar	Replacement Pads	
Super Short**	MJR-200S	MJ-200S	MJP-100	
Stout	MJR-200L	MJ-200L	MJP-200L	
X-Small	MJR-250	MJ-250	MJP-250	
Short	MJR-300	MJ-300	MJP-100	
Regular	MJR-400	MJ-400	MJP-100	
Tall	MJR-500	MJ-500	MJP-100	


^{*}Collar sets come with one additional set of replacement pads.

Sizing Chart

Osteoporitic and/

or ankylosing spondylitic

Speciality sizes: 30% of the population




Standard sizes: 60% of the population

shoulders

(MJ-200B)

*For a neck circumference >25", try the Miami J Bariatric Collar



 Short-necked women

Phenotype based sizing system makes immobilisation.

The Miami J front and back panels may be interchanged to accommodate patient anatomy, e.g. a short front may be used with a tall back. Choose collar size based on patient neck height, not patient body size. Super short consist of a regular back panel and a super short front panel.

500 Does your patient have a long, tall neck?

size selection simple whilst drastically reducing margin of error in fitting a patient. Each collar has bioengineered size specific front and back sections designed to work together for superior

22

^{**}Super Short collars are packaged with Regular Back Panels.

Aspen Vista® **Universal Collar** Features and Benefits

- Innovative height adjustment technology
- Adjustable front depth using a unique dial adjustment mechanism
- The Vista collar is six collars in one, allowing for improved fit
- Reducing storage and inventory costs
- 100% cotton latex free

Indications for Use

- Post-trauma
- Stable cervical fractures
- Cervical spondylitis
- Cervical spondylosis
- R/A & O/A cervical spine
- · Herniated cervical disc
- Postoperative
- Motor Neurone Disease

Note: *The Collar Set comes with one collar and an additional set of replacement pads.

Ordering Information

Description	Size	Product Code
Aspen Vista Collar	Universal	984000
Aspen Vista Collar Set*	Universal	984002

Accessories	Product Code
Replacement Pads (Sets)	984020
Back Panel - Standard	984006
Back Panel - Large	984008

Michael's Craftsmanship **Creates Positive Turning Points**

For adults and children with self-injurious behaviours and sensory disorders, orthotic helmets can provide life-changing treatment. At Steeper, one of our company's core values is Patient Focus and Michael embodies this value every day when he crafts top quality bespoke helmets for our patients.

Case Study

Michael's skills and attention to detail have been instrumental in the development of new helmet designs that are more comfortable and effective for patients. For example, a patient with a sensory disorder who wore protective earmuffs came to our orthotic clinic needing an orthotic helmet. He could have crafted a helmet big enough to go over the patients' earmuffs, but if they were to be removed, the helmet wouldn't fit anymore. Michael knew there had to be a better way of giving the patient the noise reduction they required from the earmuffs, and the protection they needed from the helmet, so created the first-ever Steeper helmet with built-in earmuffs.

He's also pioneered helmets with tinted visors for patients who are intolerant to daylight.

He's worked with many clinicians over the years, and he says, "My tip for Orthotists to get the best result first time would be to focus on the circumference and earto-ear measurements, if these are correct the fit will be too! However. if it's not suitable to measure your patient or if you can't, I'll happily work with you to get the best result."





SPINE & TORSO

Universal Abdominal Binder

Features and Benefits

- Neoprene binder with touch and close fastening and bound edges
- Universal sizing to reduce stock holding
- · Adjustable to allow for stoma care

Indications for Use

Post-op support/ventral hernia

Ordering Information

Size	Front Depth (cm)	Back Depth (cm)	Product Code
Universal	22	22	NB



Universal Deep Abdominal Binder

Features and Benefits

- Deep abdominal binder with simple touch and close fastening
- Adjustable to allow for stoma care
- Universal sizing to reduce stock holding

Indications for Use

Post-op support/ventral hernia

Ordering Information

Size	Front Depth (cm)	Back Depth (cm)	Product Code
Universal	22	30	AB



Universal Temporary Hernia Support

Features and Benefits

- A separate 150mm diameter, 18mm thick soft hernia pad can be located as
- Supplied with 1x hernia pad, additional pads can be purchased separately

Indications for Use

- · Abdominal hernia
- General posture support
- Lumbar strains
- Sciatica
- Lumbosacral strains

Size	Front Depth (cm)	Product Code
Universal	23	1043
Universal	30	1044

Accessories	Product Code
Hernia Pad	1041/P



Standard Maternity Support

Features and Benefits

- · Providing gentle, yet firm support to alleviate lumbar pain
- · Lift and back support can be controlled independently
- · Soft elasticated insets and shaping provide comfort
- Easy and inconspicuous to wear under most clothing
- Support expands to accommodate growth during various stages of pregnancy



- Abdominal support
- Lower back pain



Dual Maternity Belt

Features and Benefits

- Provides support to the sacroiliac and symphysis pubis to ease sacroiliac, supra pubic or low back pain
- The 3-touch-and-close strap configuration allows the patient to obtain optimum fit through a range of strap positions
- Strong, lightweight elastic gives excellent support without being restrictive

Indications for Use

- Sacroiliac support
- Low back pain
- Suprapubic support

Ordering Information

Size	Abdomen Circ. (cm)*	Product Code
1	63 - 81	SMS/S
2	81 - 112	SMS/M
3	112 - 142	SMS/L
4	142 - 172	SMS/XL

*Abdominal circumference should be measured at the widest point.



More Information Wear and Care Instructions

Ordering Information

Size	Hip Circ. (cm)*	Product Code
XX-Small	90 - 96	331/MS
X-Small	96 - 102	332/MS
Small	102 - 108	333/MS
Medium	108 - 116	334/MS
Large	116 - 126	335/MS
X-Large	126 - 136	336/MS
XX-Large	136 - 146	337/MS
XXX-Large	146 - 156	338/MS

^{*}Circumference around hips and under bump.

30

PAEDIATRIC



PRIM MAX Medium Lumbosacral Back Support

Features and Benefits

- Front parts made with seamless technology
- · Support systems that facilitate placement and fitting
- Micro hook and loop fastening for secure closure and fastening
- Paravertebral steel stays ensure maximum containment
- Latex free



- Lumbago
- Lumbosciatica
- Lumbocruralgia
- Discarthrosis
- · Disc disease
- Osteoporosis
- Osteomalacia
- Postoperative disk herniation





Ordering Information

Size	Hip Circ. (cm)	Front Height (cm)	Back Height (cm)	Product Code
X-Small	75 - 85	20	33	6033G-XS
Small	80 - 95	20	33	6033G-S
Medium	90 - 105	20	33	6033G-M
Large	105 - 115	20	33	6033G-L
X-Large	115 - 125	20	33	6033G-XL
XX-Large	125 - 135	20	33	6033G-XXL





PRIMSPINE Semi-Rigid Lumbosacral Back Support

Features and Benefits

- Semi-rigid back support belt
- Made of "BoostFit" stretch fabric
- Highly breathable
- Double channel structure that enhances the adaptation to the body, providing better compression and comfort
- Low profile option has a front height of 17cm and back height of 26cm
- · High profile option has a front height of 20cm and back height of 30cm
- · Fully breathable material with moisture wicking properties
- · Sleek slimline design

Indications for Use

- Lower back pain
- Lumbosciatica
- Osteoarthritis
- Incipient osteoporosis
- Osteomalacia
- Disc protrusion
- Postoperative recovery
- Muscle weakness
- Sprain

		Product Code		
Size	Hip Circ. (cm)	Low Profile F: 17cm / B: 26cm	High Profile F: 20cm / B: 30cm	
X-Small	70 - 80	PRS630G-XS	PRS631G-XS	
Small	77 - 90	PRS630G-S	PRS631G-S	
Medium	87 - 100	PRS630G-M	PRS631G-M	
Large	95 - 110	PRS630G-L	PRS631G-L	
X-Large	107 - 120	PRS630G XL	PRS631G-XL	
XX-Large	117 - 130	PRS630G-XXL	PRS631G-XXL	



Oasis LSO

Features and Benefits

Oasis is our most popular bivalve spinal system and ideal for patients who require maximum support, but cannot tolerate most rigid LSO braces.

- The traction closure system provides easy compression adjustments
- Revolutionary spacer fabric offers increased breathability, heat and moisture removal, less bulk and longer life span than trilaminates
- Flexibility of closing mechanism ensures an excellent fit across various body shapes
- Effective between L2 L5



- Burst fracture
- Post-op lumbar spine
- Stabilisation
- · Chronic back pain
- Compression fracture
- Degenerative, herniated or bulging disc
- Spinal instability
- Post-op disc surgery
- Spondylolysis
- Spondylolisthesis
- Spinal stenosis
- Facet syndrome
- Tumours or inflammation
- Post-op laminectomy

Ordering Information

Size	Waist Circ. (cm)	Hip Circ. (cm)	Product Code
X-Small	66 - 82.3	78.7 - 94	OALS-XS
Small	76.2 - 91.4	88.9 - 104.1	OALS-SM
Medium	86.4 - 101.6	99 - 114.3	OALS-MD
Large	96.5 - 111.8	109.2 - 124.5	OALS-LG
X-Large	106.7 - 121.9	119.4 - 134.6	OALS-XL
XX-Large	116.8 - 132	129.5 - 144.8	OALS-XXL
XXX-Large	127 - 142.2	139.7 - 154.9	OALS-XXXL





Terrapin LSO Bivalve Brace

Features and Benefits

- Bivalve design allows for easy donning and doffing
- Provides abdominal compression
- Adjustable A/P control
- Removable and washable liner
- Made from high density polyethylene
- Standard option has a front height of 10" (25.4cm)
- Short option has a front height of 81/2" (21.6cm)

Indications for Use

- Post operative spinal immobilisation
- Compression fracture
- Spinal instability
- Chronic back pain
- Degenerative, herniated or bulging disc

Size	Hip Circ.	Hip Circ.	Product Code		
Size	(cm)	(inch)	Standard	Short	Pendulous Abdomen
1	76.2 - 86.4	30 - 34	LSO-1	-	-
2	86.4 - 94	34 - 37	LSO-2	-	-
3	94 - 104.1	37 - 41	LSO-3	LSO-3S	-
4	104.1 - 114.3	41 - 45	LSO-4	LSO-4S	LSO-4P
5	114.3 - 124.5	45 - 49	LSO-5	-	LSO-5P



PRIM Dorsolumbar Back Support

Features and Benefits

- Thoraco-lumbar soft spinal system designed to aid numerous spinal conditions
- · Made of soft, breathable fabric
- Washable
- Shoulder straps promote scapular retraction
- Fully mouldable stays

Indications for Use

- Dorsalgia
- Osteoporosis
- Postoperative period
- Herniated disc
- Spondyloarthrosis
- Spondylolisthesis
- Spondylolysis
- Other degenerative disc diseases
- Hyperkyphosis



PRIM SPINAIR SP500 TLSO

Features and Benefits

- A complete and functional orthosis for the progressive treatment of muscular hypotonia and osteoporosis
- Breathable
- Washable
- Padded straps to protect underarms
- Fully adaptable "Elastic on/off" zone at the base for maximum fastening
- Easy to use level 1-3 tensioning system
- Fully mouldable lightweight aluminium support

Indications for Use

- Osteoporosis
- Hyperkyphosis
- Vertebral arthritis
- Low back pain
- Spondyloarthrosis

Ordering Information

Size	Hip Circ. (cm)	Front Height (cm)	Back Height (cm)	Product Code
X-Small	75 - 85	23	55	6037G-XS
Small	80 - 95	23	55	6037G-S
Medium	90 - 105	23	55	6037G-M
Large	105 - 115	23	55	6037G-L
X-Large	115 - 125	23	55	6037G-XL
XX-Large	125 - 135	23	55	6037G-XXL

Ordering Information

Size	Hip Circ. (cm)	Height from S1 to C7 (cm)	Product Code
S1	78 - 93	51 - 59	SP500-S1
M1	90 - 105	51 - 59	SP500-M1
M2	90 - 105	59 - 68	SP500-M2
L2	103 - 120	59 - 68	SP500-L2
L3	103 - 120	68 - 75	SP500-L3

Accessories	Product Code
SPINAIR Accessory (Extender Band For SP500)	SP500ACC

36

Hyperextension Brace



Features and Benefits

- · Lightweight aluminium construction
- Padded pelvic band can be locked or free to adjust as patient changes from standing to sitting, avoiding brace shift
- Simple height and width adjustment of sternal component
- Three step auto-lock system allows for easier access
- · Lightweight and hard-wearing polymer hinge with no metal springs guarantees the right amount of movement between the frames and pelvic band



- Immobilisation following inter-vertebral disc surgery
- Stable vertebral compression fractures of the lower thoracic and lumbar spine





Ordering Information

Description	Pelvis Circ. (cm)	Height of Brace (cm)*	Product Code
Small	60 - 75	34 - 42	MP100S
Medium Short	75 - 90	34 - 41	MP100MS
Medium	75 - 90	36 - 46	MP100M
Medium Long	75 - 90	40 - 52	MP100ML
Large Short	90 - 105	36 - 46	MP100LS
Large	90 - 105	40 - 52	MP100L
X-Large Short	105 - 115	39 - 49	MP100XLS
X-Large	105 - 115	42 - 54	MP100XL

*Measured from the proximal side of the body of the sternum to the upper edge of the pubis (pubic tubercle)







Note: To fit the VENUM Universal HYBRID TLSO correctly, please choose the desired size on the slide panel from the charts on Page 41.

VENUM Universal HYBRID TLSO

Features and Benefits

- Reduced stock holding
- Cincher technology
- Cost effective
- Step-down design
- Streamline design
- Mechanically advanced single or dual compression closure system provides a 6:1 pull ratio
- Posterior panel is available in independent or central support panel designs which creates a natural support system for the spine that can be customised to provide the perfect fit
- Ergonomic grips can be customised to patient size by extending or retracting cable length
- Made from durable, breathable, soft. hypo-allergenic materials that wick away perspiration
- Fits waist circumferences from 60cm -157.5cm; optional extender panel adds 25cm to circumference
- Effective between T7 L5

Indications for Use

- Chronic back pain
- **Kyphosis**
- Osteoporosis or osteoporotic fractures
- Degenerative, herniated or bulging disc
- Spondylolysis
 - Spondylolisthesis
- Spinal stenosis
- Facet syndrome
- Tumours
- · Spinal instability
- · Burst fracture
- Pre or post-op spinal lumbar stabilisation
- Pre or post-op disc surgery

Description	Size	Product Code
VENUM Universal HYBRID TLSO	Universal	VIPHYBTLSO3DLK
Extender Panel	Universal	STEXT



VENUM Universal TLSO PPX

Features and Benefits

The VENUM Universal TLSO PPX is the ultimate one-size-fits-all semi-rigid comfort brace for the thoraco lumbar-sacral region. This brace reduces pain, whilst providing support and immobilisation to the thoraco lumbar-sacral spine. Featuring unique anterior pads and a highly effective cincher system, the VENUM Universal TLSO PPX is designed to produce maximal intra-abdominal pressure to distract the vertebral bodies of the spine. As the brace is universal in nature, it benefits from fitting every patient type and can be stored easily where limited storage is available.

- Reduced stock holding
- Cincher technology
- Cost effective
- Step-down, streamline design
- Motion restriction of the thoracic and lumbar regions, T7 L5
- Mechanically advanced single or dual compression closure system provides a 6:1 pull ratio
- Posterior panel is available in independent or central support panel designs which creates a natural support system for the spine that can be customised to provide the perfect fit
- Ergonomic grips can be customised to patient size by extending or retracting cable length
- Made from durable, breathable, soft, hypo-allergenic materials that wick away perspiration
- Fits waist circumferences from 60cm 157.5cm; optional extender panel adds 25cm to circumference

Indications for Use

- · Chronic back pain
- Degenerative, herniated or bulging disc
- Spondylolysis
- Spondylolisthesis
- Spinal stenosis
- Facet syndrome
- Tumours
- Spinal instability
- Burst fracture
- Pre or post-op spinal thoracolumbar stabilisation
- Pre or post-op disc surgery

Note: To fit the VENUM Universal TLSO PPX correctly, please choose the desired size on the slide panel from the charts below.

Ordering Information

Description	Size	Product Code
VENUM Universal TLSO PPX	Universal	VIPHYBTLSO4PPX
Extender Panel	Universal	STEXT

VENUM Universal TLSO and Hybrid TLSOSize Chart

Size	Waist Measurement (cm)	Size	Waist Measurement (cm)
X-Small	60 - 65	X-Small/Small	65 - 70
Small	70 - 75	Small/Medium	75 - 80
Medium	80 - 85	Medium/Large	85 - 90
Large	90 - 95	Large/X-Large	95 - 100
X-Large	100 - 105	X-Large/2X-Large	105 - 110
2X-Large	110 - 115	2X-Large/3X-Large	115 - 120
3X-Large	120 - 125	3X-Large/4X-Large	125 - 130
4X-Large	130 - 135	4X-Large/5X-Large	135 - 140
5X-Large	140 - 145	Add extender panel	157.5 - 182.5
Panel fully open	142.5 - 157.5	-	-





Oasis TLSO


Features and Benefits

Oasis is our most popular bivalve spinal system and ideal for patients who require maximum support, but cannot tolerate most rigid TLSO braces. Revolutionary spacer fabric offers increased breathability, heat and moisture removal, less bulk and longer life span than trilaminates. The Oasis TLSO includes a Dorsal Lumbar Kit (DLK) with Sternal Panel to prevent forward flexion. Effective between T7 - L5.



- · Burst fracture
- Chronic back pain
- Compression fracture
- Degenerative, herniated or bulging disc
- Kvphosis
- Laminectomy
- Osteoporosis
- Paraplegia
- Postoperative
- Spinal stabilisation
- Spinal stenosis
- Spondylolysis
- Spondylolisthesis





Ordering Information

Size	Waist (cm)	Hip (cm)	Product Code
X-Small	66 - 82.3	78.7 - 94	OATL4-XS-PPK
Small	76.2 - 91.4	88.9 - 104.1	OATL4-S-PPK
Medium	86.4 - 101.6	99 - 114.3	OATL4-M-PPK
Large	96.5 - 111.8	109.2 - 124.5	OATL4-L-PPK
X-Large	106.7 - 121.9	119.4 - 134.6	OATL4-XL-PPK
XX-Large	116.8 - 132	129.5 - 144.8	OATL4-XXL-PPK
XXX-Large	127 - 142.2	139.7 - 154.9	OATL4-XXXL-PPK



More Information

Video: How to fit the Oasis TLSO brace

Terrapin Bivalve TLSO Brace

Features and Benefits

- The Bivalve technology allows for easy donnina
- The Ilio-loc provides abdominal compression, positioning, and prevents migration
- Adjustable A/P control
- Removable and washable liner
- Made from a high density polyethylene
- Standard option has a front height of 10" (25.4cm)
- Short option has a front height of 81/2" (21.6cm)
- Effective between T7 L5



- Burst fracture
- Chronic back pain
- Compression fracture
- Degenerative, herniated or bulging disc
- **Kyphosis**
- Laminectomy
- Osteoporosis
- Paraplegia
- Postoperative
- Spinal stabilisation
- Spinal stenosis
- Spondylolysis
- Spondylolisthesis

Ordering Information

Size Hip Circ.		Hip Circ.	Product Code		
Size	(cm)	(inch)	Standard	Short	Pendulous Abdomen
1	76.2 - 86.4	30 - 34	TLSO-1	-	-
2	86.4 - 94	34 - 37	TLSO-2	-	-
3	94 - 104.1	37 - 41	TLSO-3	TLSO-3S	-
4	104.1 - 114.3	41 - 45	TLSO-4	TLSO-4S	TLSO-4P
5	114.3 - 124.5	45 - 49	TLSO-5	-	-



More Information

Video: How to fit the Terrapin TLSO brace

PAEDIATRIC



CTLSO Bivalve

Features and Benefits

A rigid modular spinal system designed to limit spinal movement to promote the healing of spinal fractures. Each brace is padded with breathable foam to improve comfort for the patient.

- A rigid spinal system designed to provide patients with maximum stability and support to the cervical-lumbar spine from C3 - L5
- · Restricts and controls unwanted motion
- Durable and easy to clean and can be easily modified
- Combines rigid anterior and posterior polymer shells
- Chape and hook and loop fastening system

Indications for Use

- · Postoperative internal fixation
- Multiple level fractures
- · Burst fracture
- Compression fracture
- · Spinal stenosis

Ordering Information

Female				
Size	Xyphoid (cm)	Waist (cm)	Hip (cm)	Product Code
X-Small	57.5 - 72.5	47.5 - 62.5	62.5 - 77.5	STKCTLSOCF-XS
Small	67.5 - 82.5	57.5 - 72.5	72.5 - 87.5	STKCTLSOCF-SM
Medium	77.5 - 92.5	67.5 - 82.5	82.5 - 97.5	STKCTLSOCF-MD
Large	87.5 - 102.5	77.5 - 92.5	92.5 - 107.5	STKCTLSOCF-LG
X-Large	97.5 - 112.5	87.5 - 102.5	102.5 - 117.5	STKCTLSOCF-XL
XX-Large	107.5 - 122.5	97.5 - 112.5	112.5 - 127.5	STKCTLSOCF-2X
XXX-Large	117.5 - 132.5	107.5 - 122.5	122.5 - 137.5	STKCTLSOCF-3X
Male				
Size	Xyphoid (cm)	Waist (cm)	Hip (cm)	Product Code
X-Small	72.5 - 87.5	60 - 75	70 - 85	STKCTLSOCM-XS
Small	82.5 - 97.5	70 - 85	80 - 95	STKCTLSOCM-SM
Medium	92.5 - 107.5	80 - 95	90 - 105	STKCTLSOCM-MD
Large	102.5 - 117.5	90 - 105	100 - 115	STKCTLSOCM-LG
X-Large	112.5 - 127.5	100 - 115	110 - 125	STKCTLSOCM-XL
XX-Large	122.5 - 137.5	110 - 125	120 - 135	STKCTLSOCM-2X
XXX-Large	132.5 - 147.5	120 - 135	130 - 145	STKCTLSOCM-3X

Made to Measure

Features and Benefits

Our custom-built manufacturing facility enables us to rapidly manufacture and control every stage of the production process, ensuring quality and on-time delivery of any orthosis required.

We can manufacture your spinal orthosis utilising the latest CAD/CAM technology which allows us to carve and mould each spinal orthosis to your patients unique measurement or scanned image, and your specification. This ensures the very best fit, support, and comfort.

Please contact customer service at **orthotics@steepergroup.com** for more information.



Made to Measure Cast
TLSO Soft Body
Anterior




Made to Measure TLSO Bivalve (ATE version)



Made to Measure Scoliosis Brace





ARM & SHOULDER

HEAD & NECK



Epicondylitis Strap

Features and Benefits

- Epicondylitis strap on a breathable fabric support, with seamless "high-frequency" technology for better skin protection
- Anatomical and re-positionable condyle pad that unloads the tendon and produces a massage effect
- Adjustable elastic strap for dynamic compression adjustment
- Anatomical structure for enhanced forearm fit
- Hypoallergenic non-slip support to prevent the strap from rotating

Indications for Use

- Sport or work-related lateral epicondylitis (tennis elbow)
- Medial epicondylitis (golfer's elbow)



Ordering Information

Size	Product Code
Maximum elbow circumference 30cm	T29



Universal High Arm Sling

Features and Benefits

- · Lightweight breathable mesh
- Touch and close strap fastening

Indications for Use

• Forearm, elbow and arm immobilisation

Ordering Information

Size	Product Code
Universal	0860



Tennis Elbow Sleeve

Features and Benefits

- Fully breathable easy to use pull sleeve wicks moisture away creating a cooling effect
- The strap allows control of compression for maximum effect and pain relief
- Pre-formed shaping for improved intimacy of fit, offering effective pain relief

Indications for Use

Medial and lateral epicondylitis

Size	Elbow Joint Circ. (cm)	Product Code
Small	20.5 - 23	TES/CX/BE/S
Medium	23 - 25.5	TES/CX/BE/M
Large	25.5 - 28	TES/CX/BE/L
X-Large	28 - 31	TES/CX/BE/XL



ARM & SHOULDER



Formfit® Shoulder Brace

Features and Benefits

- Streamlined fit and sporty look
- Breathable and comfortable
- Provides abduction of the shoulder (when the abduction pillow is attached), or external rotation of the shoulder (when the external rotation pillow is attached)
- Available with a customisable aeroplane kit for patients who require precise positioning following shoulder surgery
- Optional axilla pillow available for purchase to help distribute weight over a larger surface area and decrease pressure on the ulnar nerve



- Immobilisation for rotator cuff repairs
- · Bankart repairs
- Capsular shifts
- Glenohumeral dislocations/subluxations
- Soft tissue repairs/strains
- Shoulder instabilities/injuries
- Post-surgical rehabilitation

Ordering Information

Size	Length (cm)	Product Code			
312e	Length (cm)	Regular	External Rotation	With Abduction	
Small	Up to 29	SSS-SO	SSS-ER	SSS	
Medium	29 - 40	SSR-SO	SSR-ER	SSR	
Large	40 - 44.5	-	-	SSL	



Neuro-Lux Shoulder II® (Sporlastic)

Features and Benefits

- Freedom of movement allows virtually unrestricted use of affected arm
- Orthosis can be taken on and off with one hand
- Different coloured plastic buckles make it easy to put on correctly
- Improved repositioning of joint (humerus) whilst correcting subluxation and inner rotation
- Stabilisation of shoulder joint
- Prevents secondary injuries of crown cap, tendons, muscles and shoulder nerves
- Improvement of posture through movement of arm pendulum
- Avoidance of shoulder-hand syndrome

Indications for Use

- Flaccid and painful movement of the shoulder and arm musculature after hemiplegia
- Plexus injuries
- · Craniocerebral trauma
- Cervical spine trauma
- Peripheral nerve injuries


Sizo	Size Lower Arm Chest Cir		Product Code	
Size	Circ (cm)	(cm)	Left	Right
1	21 - 25	75 - 90	07261/1L	07261/1R
2	21 - 25	90 - 105	07261/2L	07261/2R
3	21 - 25	105 - 125	07261/3L	07261/3R
4	25 - 29	75 - 90	07261/4L	07261/4R
5	25 - 29	90 - 105	07261/5L	07261/5R
6	25 - 29	105 - 125	07261/6L	07261/6R
7	25 - 35	75 - 90	07261/7L	07261/7R
8	25 - 35	90 - 105	07261/8L	07261/8R
9	25 - 35	105 - 125	07261/9L	07261/9R


Rebound® Post Op Elbow

Features and Benefits

- Includes a drop-lock hinge and four telescoping paddles
- Lightweight
- Promotes user comfort with straps that stay in place
- Anatomical left and right fit with mouldable stability cuffs
- Blue patient touch-points for easy education
- humerus, proximal radius or ulna
- Tendon and ligament injuries or repairs
- Triceps tendon repair

Indications for Use





- Elbow hyperextension
- · Chronic elbow injuries



Description	Product Code
Right	B-8391580
Left	B-8391680





WRIST & HAND

HEAD & NECK



Universal Palm Protector

Features and Benefits

- Fleece-lined, lightweight palm cover
- Designed to allow full finger movement whilst reducing finger contractures and skin damage/breakdown in the palmar region
- Easily laundered
- Easily applied using hook and loop fastening straps

Indications for Use

- · Post surgery/injury positioning
- Hypertonic contractures
- Rheumatoid arthritis

Ordering Information

Description	Size	Product Code
Right	Universal	PPR
Left	Universal	PPL



CMC Brace

Features and Benefits

- Support for the thumb CMC joint in the treatment of osteoarthritis and medial subluxation. The shaped strap wraps around the thumb CMC joint and the first metacarpal for targeted CMC joint support
- Made from Neoprene

Indications for Use

- · CMC joint arthritis
- CMC instability
- Tendinitis
- Pain associated with repetitive motion

Ordering Information

56

Sina Dalm Sina (am)	Product Code		
Size	Palm Circ. (cm)	Left	Right
Small	16.5 - 17.5	CMC/SL	CMC/SR
Medium	17.5 - 20	CMC/ML	CMC/MR
Large	20 - 24	CMC/LL	CMC/LR



Thumb Spica - Neoprene

Features and Benefits

- Neoprene thumb spica featuring a removable flexible stay which restricts movement at the CMC and MCP joints, leaving the IP joint free
- Warmth and compression
- Removable flexible stay

Indications for Use

- CMC joint arthritis
- CMC instability
- Tendinitis
- Pain associated with repetitive motion
- Osteoarthritis
- Painful/unstable 1st MCP joint
- Rheumatoid arthritis
- Thumb strains
- Sprains and post-cast removal

Size Wrist Circ. (cm)		Product Code	
Size	Wrist Circ. (cm)	Left	Right
Small	14 - 15	TS/BE/S/L	TS/BE/S/R
Medium	15 - 16.5	TS/BE/M/L	TS/BE/M/R
Large	16.5 - 18	TS/BE/L/L	TS/BE/L/R
X-Large	18 - 19	TS/BE/XL/L	TS/BE/XL/R



Push OrthoThumb Brace CMC

Features and Benefits

- Hand remains stable and allows good grip function
- Provides pain relief when experiencing thumb stress
- Surrounding joints remain free and doesn't limit their movement
- Pre-formed flexible plastic enclosure
- Aluminium splint can be manipulated for individual thumb contours and optimal stabilisation of the joint
- Two hook and loop fastening straps fasten around the back of the hand which allows the brace to be easily fitted by the individual
- The material and rounded edges allows for maximum comfort
- Antimicrobial
- Can be used in water and can be dried easily

Indications for Use

- Osteoarthritis of the CMC-1 joint
- Postoperative after care of the CMC-1 joint
- Status after arthroplasty of the CMC-1 ioint
- Instability of the CMC-1 joint

Measure the circumference of the affected hand and use the chart to pick the correct size. If you are in between sizes, select the larger size.



Ordering Information

Size	MCD Cine (om)	Product Code	
Size	MCP Circ. (cm)	Left	Right
0	15 - 17.5	PUS122B	PUS121B
1	17.5 - 19.5	PUS067B	PUS064B
2	19.5 - 22.5	PUS068B	PUS065B
3	22.5 - 27	PUS069B	PUS066B

Push Care Wrist Brace



Features and Benefits

- Provides pressure around the wrist joint to improve stability
- Allows for adjustable pressure for individual needs, the wrist joint is not immobilised and functionality of the hand is retained
- The palm of the hand maintains freedom, ensuring a good grasp can be achieved
- Padded material around the ulnar area of the wrist
- Sympress material inner offers maximum comfort and keeps the skin dry under the brace



- · Residual instability
- Moderate to serious distortions
- Strain complaints in the wrist joint
- Therapy resistant strain complaints to the wrist extensors
- · Wrist fractures without dislocations
- · Aftercare for wrist fractures
- Arthrosis with distension
- RSI (repetitive strain injury)

Size	Muich Cine (cm)		Product Code
Size	Wrist Circ. (cm)	Left	Right
1	13 - 15	PUS074B	PUS070B
2	15 - 17	PUS075B	PUSO71B
3	17 - 19	PUS076B	PUS072B
4	19 - 21	PUS077B	PUS073B

Push Med Wrist Brace

Features and Benefits

- Provides high levels of support and protection to the wrist joint without restriction of the functionality of the hand
- Anatomically shaped stiffeners limit palmary and dorsal flexion without inhibiting it completely
- Elastic strap provides adjustable compression
- Easy to fit
- Can be applied using one hand meaning the individual can fit it themselves
- Padded material around vulnerable areas of the wrist



Indications for Use

- Residual instability
- Moderate to serious distortions
- Strain complaints in the wrist
- Therapy-resistant strain complaints of the wrist extensors
- · Wrist fractures without dislocations
- Follow-up treatment of wrist fractures
- Arthrosis with hydrops
- · Repetitive strain injury



Short Wrist Brace

Features and Benefits

- Standard single compression one-piece wrist brace provides good conformity around the wrist and hand
- Removable, mouldable metal stay
- Made from Neoprene

Indications for Use

- · Flaccid drop wrist
- Repetitive strain injury
- Rheumatoid arthritis
- · Carpal tunnel syndrome

Ordering Information

Size	Muich Cine (cm)	Product Code	
Size	Wrist Circ. (cm)	Left	Right
1	13 - 15	PUS082B	PUS078B
2	15 - 17	PUS083B	PUS079B
3	17 - 19	PUS084B	PUS080B
4	19 - 21	PUSO85B	PUS081B

Size Wrist Circ. (cm	Muist Cius (sus)	Product Code	
	Wrist Circ. (cm)	Left	Right
X-Small	10 - 12.5	SWB/E/XS/L	SWB/E/XS/R
Small	12.5 - 15	SWB/E/S/L	SWB/E/S/R
Medium	15 - 17.5	SWB/E/M/L	SWB/E/M/R
Large	17.5 - 20	SWB/E/L/L	SWB/E/L/R
X-Large	20 - 22.5	SWB/E/XL/L	SWB/E/XL/R

Universal Wrist Brace

Features and Benefits

- Breathable felt and neoprene mix
- Aluminium palmer splint and stays can be shaped to fit the wrist with additional hook and loop strap for enhanced fit
- Designed for easy on and off application
- Circumferential seam for improved fit around waist

Indications for Use

- Arthritis

- Carpal Tunnel

Note: Paediatric sizes also available, see Page 214.



- Sprain
- Swelling

Ordering Information

Description	Size	Product Code
Left	Adult	SWB/L/UNIVERSAL
Right	Adult	SWB/R/UNIVERSAL





Wrist and Thumb Spica

Features and Benefits

- Two-piece dual compression with malleable aluminium thumb splint to hold thumb in abduction
- 5cm wide elasticated strap
- Made from Neoprene

Indications for Use

- Osteoarthritis
- Painful/unstable 1st MCP joint
- Rheumatoid arthritis

Size	Wrist Circ (cm)	Product Code		
Size	Wrist Circ. (cm)	Left	Right	
Small	12.5 - 15	WTB/CX/BE/S/L	WTB/CX/BE/S/R	
Medium	15 - 17.5	WTB/CX/BE/M/L	WTB/CX/BE/M/R	
Large	17.5 - 20	WTB/CX/BE/L/L	WTB/CX/BE/L/R	
X-Large	20 - 22.5	WTB/CX/BE/XL/L	WTB/CX/BE/XL/R	



Wrist and Thumb Spica - Neoprene

Features and Benefits

- Restricted wrist and thumb motion without limiting finger motion
- Easy fit and comfortable
- Lightweight and padded
- Contours to the anatomy of the wrist and thumb
- Breathable Neoprene material
- Hook and loop fastening



- Gamekeeper's thumb
- Post-cast healing
- Soft tissue injuries
- Selective scaphoid injuries



Comfyprene® Contour Hand Orthosis

Features and Benefits

- The Comfyprene® Contour Hand Orthosis has unique side flaps that hold the wrist in place allowing for full hand contouring
- C-bar and thumb support with included strap
- 4 straps come as standard
- Unique butterfly designed side straps hold wrist and forearm in place
- Contours the whole hand, wrist and forearm
- Instead of removable terrycloth covers, the Comfyprene splint inserts are covered with soft neoprene and headliner material, which creates a slim profile for everyday comfort
- Splints are adjustable and re-adjustable using the Comfy[™] patented reinforced spine
- Straps may be positioned anywhere on the splint for optimal fit
- · Contracture management

Indications for Use

- Use after sports injuries
- · Removal of casts
- Sports medicine
- Arthritic patient's support of the weakened extremity



Ordering Information

Size	Muist Cine (sue)	Product Code		
Size	Wrist Circ. (cm)	Left	Right	
Small	15 - 18	WTS/S/L	WTS/S/R	
Medium	18 - 20	WTS/M/L	WTS/M/R	
Large	20 - 23	WTS/L/L	WTS/L/R	
X-Large	23 - 25	WTS/XL/L	WTS/XL/R	

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code	
			Left	Right
Adult	17.8+	27.9	HCO-101-CP-LB-LT	HCO-101-CP-LB-RT

WRIST & HAND

Universal

Comfy Air Hand Orthosis

Features and Benefits

- The Comfy Air Hand Orthosis is an exceptionally useful splint for the severely contracted hand
- It has an adjustable air bladder that allows gradual extension of the MP, DIP and PIP joints, as well as an adjustable wrist for customised fit
- Ingeniously designed by a therapist, one simply inserts the unfitted splint with the bladder under the fingers and then pumps air to inflate the bladder and increase finger extension. (Headliner material only)
- Ambidextrous

Indications for Use

- · Ideal for the severely contracted hand
- Inflatable air bladder allows for gradual MP, DIP and PIP joint extension
- Adjustable wrist
- Provides low pressure stretch therapy
- Hand air bladder can be used separately without full orthosis as needed
- Can fit into even the most tightly clenched hand





Comfy Hand-WristFinger Orthosis

Features and Benefits

- Provides support and positioning for the weak or deformed wrist, hand or fingers
- Serves as an excellent resting splint to prevent trauma to joints and positions to increase ROM
- Side wings prevent ulnar and radial deviation
- Arm cuffs easily adjust to width of arm
- Terrycloth
- Headliner
- Ambidextrous

Indications for Use

- Flaccid drop wrist
- Repetitive strain injury
- · Rheumatoid arthritis
- Carpal tunnel syndrome

Note: Paediatric sizes also available, see Page 216.

Ordering Information

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Adult	17.8 - 20.3	31.8	HA101-A

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Adult Small	15.2 - 17.8	27.9	H101-AS
Adult Large	20.3+	34.3	H101-AL

SoftPro® Functional Resting WHFO

Features and Benefits

- "Bend-to-fit" endoskeleton frame can be hand shaped in the sagittal plane as well as medially, laterally, and rotationally to accommodate most angles of the wrist/ hand/and fingers without the use of tools
- The endoskeleton frame is a closed cell foam, covered by a laminated foam/ cloth cover that can be laundered as necessarv
- Extra long wrist straps hold the wrist in the proper position and applies a gentle
- · Maintains and supports the wrist, hand and fingers in a comfortable functional position



• Support and protection of a weak or injured wrist, hand, and fingers, including thumb




ComfySplints™ Deviation Hand Thumb Orthosis

Features and Benefits

- Unique, below the wrist, swivel that allows for lateral and medial positions to accommodate for ulnar and radial
- Progressively change wrist and finger positions and amount of wrist deviation by adjusting the internal frame and locking the frame's swivel with included
- Correct or control wrist deviation motion even when changes in the position of the wrist and hand are made
- Ambidextrous



- · Arthritic deformities
- Neuromuscular deformities
- Flaccid, weak or spastic hand conditions



Ordering Information

Size	Width of MP	Length	Product Code	
Size	Joints (cm)	(cm)	Left	Right
Small	6.5 - 7.5	23	SP31FFR-S-L	SP31FFR-S-R
Medium	7.5 - 9	25	SP31FFR-M-L	SP31FFR-M-R
Large	9 - 10	28	SP31FFR-L-L	SP31FFR-L-R

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Adult	17.8 - 20.3	31.8	DHT-101

68



ComfySplints™ Rest Hand

Features and Benefits

- Supports the weak or deformed hand with ultra-comfort
- Removable finger separator included
- Super soft and removable navy headliner cover is machine washable (laundry bag
- Treated with antimicrobial barrier for skin breakdown prevention
- · Patented reinforced spine can bend and adjust as needed (no tools or heat required)



- Arthritis
- CVA



- Sprain/strains

Ordering Information

Size	Wrist to Tip Total of Middle Length		Product Code	
Size	Finger (cm)		Left	Right
Adult Small	13.3 - 15.9	26.7	RH-101-AS-LEFT	RH-101-AS-RIGHT
Adult	15.9 - 18.4	29.2	RH-101-LEFT	RH-101-RIGHT
Adult Large	18.4+	34.3	RH-101-AL-LEFT	RH-101-AL-RIGHT





ComfySplints™ Large Pan Hand Thumb Orthosis

Features and Benefits

- Designed to support the thumb without stressing the web space between the index finger and the thumb known as the "Thenar Eminence"
- Soft straps accommodate sensitive skin, oedema and fragile skin
- Removable Terrycloth material covers are fully machine washable
- Fully and easily adjustable at all joints to make a custom fit without the use of tools or heat guns
- Side wings prevent ulnar/radial deviation
- Ambidextrous


Indications for Use

- Sprains/Strains/Fractures
- Sports injuries
- · Neurological impairments
- Contractures
- CP
- · Wound care

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Adult	17.8 - 20.3	31.8	LPHT-101

ComfySplints™ Spring Loaded Goniometer Hand Orthosis






- Features a spring that gently pushes extremity to the selected degree of
- Use the goniometer dial to set the range of flexion and extension from 0 to 180° in 10° increments
- · Allen key supplied
- Blocks flexion or extension while allowing motion in opposite direction
- Can be set in a desired static position
- Re-adjustable custom fit of wrist, MP, PIP and DIP joints
- Ambidextrous



- Contractures of the wrist and hand
- Flaccid drop wrist
- Repetitive strain injury
- Rheumatoid arthritis







Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Adult	17.8 - 20.3	31.8	SGH-101





Ortel Hernia Truss



Features and Benefits

- Removable and repositionable silicone pads increase comfort and effectiveness
- Buckle system makes unilateral version left/right adjustable

Indications for Use

• Reducible unilateral inguinal hernia

Note: This product is non-returnable for hygiene reasons.





Hernia Trusswith Bilateral Pads

Features and Benefits

 Bilateral hernia belt is a reusable medical device intended to support treatment and rehabilitation according to the relevant indications

Indications for Use

- Patients with minor inguinal hernias who cannot be subjected to surgery for various reasons
- Following double-sided inguinal hernia surgery in the case of weak muscletendon system to prevent recurrence

Note: This product is non-returnable for hygiene reasons.

Ordering Information

Size	Maist Cive (sm)	Product Code		
	Waist Circ. (cm)	Unilateral	Bilateral	
1	75 - 85	27300120100199	27200120100199	
2	86 - 97	27300120200199	27200120200199	
3	98 - 110	27300120300199	27200120300199	
4	111 - 125	27300120400199	27200120400199	

Size	Hip Circ. (cm)	Product Code
Small	70 - 79	PT0110-S
Medium	80 - 90	PT0110-M
Large	91 - 102	PT0110-L
X-Large	103 - 115	PT0110-XL
XX-Large	116 - 130	PT0110-XXL

ARM & SHOULDER



Hernia Support Briefs Men's

Features and Benefits

· Reusable medical device intended for men to support treatment and rehabilitation according to the relevant indications

Indications for Use

- Following surgery of large inguinal or femoral hernias
- In the case of postoperative oedemas of testes or epididymis
- Following surgery or puncture of hydrocele testis
- Prevention of inguinal or femoral hernia recurrence
- In patients who cannot be subjected to surgery for various reasons

Note: This product is non-returnable for hygiene reasons.



Hernia Support Briefs Women's

Features and Benefits

· Reusable medical device intended for women to support treatment and rehabilitation according to the relevant indications

Indications for Use

- Following surgery of large inquinal or femoral hernias
- In patients who cannot be subjected to surgery for various reasons
- Prevents recurrence of inguinal and femoral hernias
- Prevents oedemas and pain after inguinal and femoral hernia surgery

Note: This product is non-returnable for hygiene reasons.

Ordering Information

Size	Waist Circ. (cm)	Hip Circ. (cm)	Product Code
X-Small	<70	<85	PT0111-XS
Small	70 - 80	85 - 95	PT0111-S
Medium	81 - 90	96 - 105	PT0111-M
Large	91 - 100	106 - 115	PT0111-L
X-Large	101 - 115	116 - 125	PT0111-XL
XX-Large	116 - 130	126 - 135	PT0111-XXL

Size	Waist Circ. (cm)	Hip Circ. (cm)	Product Code
Small	70 - 80	85 - 95	PT0112-S
Medium	81 - 90	96 - 105	PT0112-M
Large	91 - 100	106 - 115	PT0112-L
X-Large	101 - 115	116 - 125	PT0112-XL
XX-Large	116 - 130	126 - 135	PT0112-XXL



VENUM UniversalSoft Hip Abduction Brace

Features and Benefits

The one-size-fits-all VENUM Soft Hip Abduction brace incorporates the most advanced technology and materials. A revolutionary soft hip management system for patients of all shapes and sizes providing maximum comfort and support. The LSO part of the brace operates using a cincher system to maximise control. The belt thigh cuff and hinge adjust to fit the dimension of most patients. The brace applies a dynamic abduction force, so that the femoral head remains in the socket.

- Universal design
- Cincher technology
- Size indicator displays belt size to ensure correct and comfortable fitting
- Joint with adjustable abduction/ adduction, flexion/extension
- Single or double-pull cable pulley system
- Stock reduction, saves money and space
- Materials are comfortable, breathable, soft, wick-away perspiration and are durable and resilient
- Waist circumference range: 60cm -157.5cm
- Extender panel adds an additional 25cm (sold separately)
- Ergonomic arthritic grips can be customised to patient size

Indications for Use

• For patients at risk of hip dislocation, suitable for all levels of instability

Material: 30% Polyurethane 70% Polyamide

Ordering Information

Description	Size	Patient Height (cm)	Product Code
VENUM Soft Hip Abduction Brace	Universal	<188	VHHIPABD
Long Joint*	Universal	>188	TLCHIPJNTLG
Extender Panel	Universal	-	STEXT

*Note: The Venum Universal Soft Hip Aduction Brace comes with a standard joint for patients of height below 188cm, for patients over 188cm a long joint will be required.

To fit the VENUM Universal Hip Brace correctly, please choose the desired size on the slide panel from the chart below.

VENUM Universal Braces Size Chart

Size	Waist Measurement (cm)	Size	Waist Measurement (cm)
X-Small	60 - 65	X-Small/Small	65 - 70
Small	70 - 75	Small/Medium	75 - 80
Medium	80 - 85	Medium/Large	85 - 90
Large	90 - 95	Large/X-Large	95 - 100
X-Large	100 - 105	X-Large/2X-Large	105 - 110
2X-Large	110 - 115	2X-Large/3X-Large	115 - 120
3X-Large	120 - 125	3X-Large/4X-Large	125 - 130
4X-Large	130 - 135	4X-Large/5X-Large	135 - 140
5X-Large	140 - 145	Add extender panel	157.5 - 182.5
Panel fully open	142.5 - 157.5	-	-



Arthritic Knee Sleeve

Features and Benefits

- Simple knee sleeve
- Pull on design, elasticated for easy use of donning
- Provides compression and warmth for painful joints
- Posterior shaping eradicates bunching in popliteal area

Indications for Use

Arthritis



Stabilised Knee Sleeve

Features and Benefits

- Incorporates medial and lateral flexible spiral stays for enhanced medial/lateral stability
- Provides compression and support
- Touch and close straps for improved fit
- Universal left/right application
- Made from Airprene

Indications for Use

- Medio lateral stability and patella control knee instability
- Post-op rehabilitation
- Strains/sprains
- Bursitis
- Tendonitis osteoarthritis
- · Rheumatoid arthritis
- Osteoarthritis

Note: When measuring the knee, ensure it is fully extended. For thigh circumference, measure 15cm above the centre of the patella. For calf circumference, measure 15cm below the centre of the patella.

Ordering Information

Size	Knee Circ. (cm)	Product Code
Small	31 - 34	803/AKS
Medium	34 - 37	805/AKS
Large	37 - 42	807/AKS
X-Large	42 - 46	808/AKS
XX-Large	49 - 54	809/AKS

Size	Thigh Circ. (cm)	Knee Circ. (cm)	Calf Circ. (cm)	Product Code
Small	34 - 36	25 - 28	26 - 28	SKS/CX/BE/S
Medium	36 - 39	28 - 30	28 - 32	SKS/CX/BE/M
Large	39 - 41	30 - 33	32 - 35	SKS/CX/BE/L
X-Large	41 - 44	33 - 36	35 - 37	SKS/CX/BE/XL



ARM & SHOULDER

PAEDIATRIC



Standard Knee Brace

Features and Benefits

- The front opening design allows easy application for the arthritic or elderly
- Metal side hinges provide medial/lateral stability while allowing flexion
- Opposing fastenings for upper and lower touch and close straps counter rotation and allow correct fitting
- Removable hinges
- Universal left/right application
- Made from Airprene

Indications for Use

- Medio lateral stability
- Osteoarthritis
- Rheumatoid arthritis

Note: When measuring the knee, ensure it is fully extended. For thigh circumference, measure 15cm above the centre of the patella. For calf circumference, measure 15cm below the centre of the patella.



Note: Measure 15cm above the centre of the patella.

Stabimed® Knee Brace

Features and Benefits

- Short soft brace with the possibility of extension and flexion limitation
- Easy to handle thanks to wrap design
- Silicone dot coating supports a secure fit on the leg
- Medi Airtex+ material with intelligent thermoregulation provides a high level of comfort
- Special knit technology at the hollow of the knee allows ease of flexion
- The physioglide hinge mimics the roll-glide movement of the knee and supports a secure fit on the leg
- The compressive spacer-knit can improve proprioception
- Comprehensive set of wedges for flexible extension and flexion limitation including 5° hyperextension wedge. The hinge design enables the wedges to be changed without using a tool

Indications for Use

- Injuries of the collateral ligaments
- Instability of the knee
- Meniscus injuries
- · Osteoarthritis of the knee

Ordering Information

Size	Thigh Circ. (cm)	Knee Circ. (cm)	Calf Circ. (cm)	Product Code
X-Small	29 - 34	18 - 22	26 - 31	SKB/CX/BE/XS
Small	34 - 39	22 - 27	31 - 37	SKB/CX/BE/S
Medium	39 - 47	27 - 29	37 - 42	SKB/CX/BE/M
Large	47 - 53	29 - 34	42 - 46	SKB/CX/BE/L
X-Large	53 - 59	34 - 40	46 - 53	SKB/CX/BE/XL

Ordering Information

Size	Thigh Circ. (cm)	Product Code
X-Small	31 - 37	826/1
Small	37 - 45	826/2
Medium	45 - 51	826/3
Large	51 - 57	826/4
X-Large	57 - 65	826/5
XX-Large	65 - 74	826/6



More Information

Wear and Care Instructions

Collamed® Knee Brace

Features and Benefits

- Long soft brace with the possibility of extension and flexion limitation
- The compressive spacer-knit can improve proprioception
- Easy to handle thanks to wrap design
- Silicone dot coating supports a secure fit on the leg
- Medi Airtex+ material with intelligent thermoregulation provides a high level of comfort
- · Special knit technology at the hollow of the knee allows ease of flexion
- The physioglide hinge mimics the roll-glide movement of the knee and supports a secure fit on the leg
- · Comprehensive set of wedges for flexible extension and flexion limitation including 5° hyperextension wedge. The hinge design enables the wedges to be changed without using a tool



- Injuries to the collateral ligament (post traumatic/postoperative)
- Injuries to the cruciate ligament (post traumatic/postoperative)
- Instability of the knee
- Meniscus injuries
- Osteoarthritis of the knee



Note: Measure 15cm above the centre of the patella.



Rebound® Post-Op Knee Brace

Features and Benefits

- Simple, tool-free design
- Multi-purpose Smart-Fit[™] packaging
- OneTouch[™] easy to adjust ROM hinge with lock-out capability
- Colour-coding and numbered frame for patient education
- Telescoping frame from 47cm 68cm (Adult) and 39.5cm - 54cm (paediatric)
- Quick-fit, gripping buckles for easy application
- The brace is lightweight and the edges are protected to improve patient comfort

Indications for Use

- ACL, PCL, MCL, LCL repairs/injuries
- Osteochondral repairs
- Meniscal repairs
- Condylar fractures
- Tibial plateau fractures
- Sprains/strains of the knee
- Patellar tendon repairs
- HTOs (High Tibial Osteotomy)

Note: Measure 15cm above the centre of the patella.

Ordering Information

Size	Thigh Circ. (cm)	Product Code
X-Small	31 - 37	824/1
Small	37 - 45	824/2
Medium	45 - 51	824/3
Large	51 - 57	824/4
X-Large	57 - 65	824/5
XX-Large	65 - 74	824/6
XXX-Large	74 - 84	824/7

Size	Description	Thigh Circ. (cm)	Product Code
Universal	Cool	41 - 81	B-8381480
Universal	Full Foam	41 - 81	B-8382480



X-ROM™ Post-Op Knee Brace

Features and Benefits

- Fast, simple application and adjustment through innovative strap, tab and buckle designs
- Precisely control ROM with settings in 10° increments from -10° to 90° extension and -10° to 120° flexion
- Four independently telescoping sliders
- Ability to contour to the patient's leg to aid varus and valgus support
- Comfortably accommodates a wide range of patient shapes and sizes (from 4'11" to 6'5")
- Enhanced ankle pads and strategically placed contact points between the leg and brace help reduce slippage
- Anti-microbial pads help neutralise odour while increasing comfort and durability

Indications for Use

 To aid immobilization and protected range of motion associated with ACL, PCL, LCL, and MCL surgeries and meniscal repairs

Size	Colour	Product Code
Universal	Black	11-2181-9



ARM & SHOULDER



OA Ease Knee Brace

Features and Benefits

- Wraparound design
- Simple rubber pull-tabs, colour-coded guide dots and a convenient hinge alignment guide make application and removal simple
- 3-Point Leverage System with Dynamic Force System™ strap and shells unload the affected knee compartment
- Lightweight and breathable fabric
- Machine washable
- Low-profile brace with a slim hinge can be worn under clothes
- Fine-tune strap tension for individual pain relief, with marked dosing indicator lines
- Extension kits and liners are available to further enhance fit and comfort

Indications for Use

- Mild-to-moderate OA pain
- Mild-to-moderate unicompartmental knee osteoarthritis and/or degenerative meniscal tears
- Other conditions that may benefit from unicompartmental unloading





Note: For thigh circumference, measure 15cm above the centre of the patella. For calf circumference, measure 15cm below the centre of the patella.

Medial				
Size	Calf Circ.	Thigh Circ.	Product Code	
Size	(cm)	(cm)	Right Medial	Left Medial
X-Small	28 - 32	36 - 41	FOA0518711	FOA0618711
Small	32 - 36	41 - 46	FOA0518712	FOA0618712
Medium	36 - 41	46 - 53	FOA0518713	FOA0618713
Large	41 - 51	53 - 60	FOA0518714	FOA0618714
X-Large	46 - 51	60 - 67	FOA0518715	FOA0618715
XX-Large	51 - 56	67 - 75	FOA0518716	FOA0618716
XXX-Large	56 - 61	75 - 81	FOA0518717	FOA0618717

Lateral				
Size	Calf Circ. Thigh Circ.		Product Code	
Size	(cm)	(cm)	Right Lateral	Left Lateral
X-Small	28 - 32	36 - 41	FOA0528711	FOA0628711
Small	32 - 36	41 - 46	FOA0528712	FOA0628712
Medium	36 - 41	46 - 53	FOA0528713	FOA0628713
Large	41 - 51	53 - 60	FOA0528714	FOA0628714
X-Large	46 - 51	60 - 67	FOA0528715	FOA0628715
XX-Large	51 - 56	67 - 75	FOA0528716	FOA0628716
XXX-Large	56 - 61	75 - 81	FOA0528717	FOA0628717

Unloader® One X


Features and Benefits

- 3-Point Leverage System consisting of two Dynamic Force System (DFS) straps, rigid shells and a single upright hinge
- No condyle pressure during sitting or swing phase
- A slender, perforated frame with soft edges
- Discreet whether worn over or under clothing
- Skin-friendly liners and pads as well as a smooth inner surface of shells reduce the likelihood of pressure points
- The gastroc strap is positioned over the thickest part of calf to help prevent slipping
- Tool-free SmartDosing® dials enable the patient to fine-tune the strap tension
- Colour coded touchpoints and Quick-Fit buckles make it easy to use
- Alignment guide on hinge
- Machine washable (washbag included with product)

Indications for Use

- Moderate to severe single compartment knee osteoarthritis (OA)
- Degenerative meniscal tears
- Other unicompartmental knee conditions that may benefit from unloading, such as:
 - Articular cartilage defect repair
 - Meniscal repair
 - Avascular necrosis
 - Condylar bone marrow lesions
 - Tibial plateau fractures

Note: Recommended for patients with inside leg measure over 85cm.




Ordering Information

Standard	Standard (Leg length more than 85cm)					
Size	Calf Circ.	Product Code				
Size	(cm)	Right Medial	Left Medial	Right Lateral	Left Lateral	
X-Small - Medium	25 - 43	UOX0518713	UOX0618713	UOX0528713	UOX0628713	
Large - X-Large	43 - 62	UOX0518715	UOX0618715	UOX0528715	UOX0628715	

Note: Recommended for patients with inside leg measure under 85cm.

For calf circumference, measure 15cm below the centre of the patella.





Short (Le	Short (Leg length less than 85cm)					
Size	Calf Circ.	Product Code				
Size	(cm)	Right Medial	Left Medial	Right Lateral	Left Lateral	
X-Small - Medium	25 - 43	UOX1518713	UOX1618713	UOX1528713	UOX1628713	
Large - X-Large	43 - 62	UOX1518715	UOX1618715	UOX1528715	UOX1628715	

92





Ankle Sleeve Plus

Features and Benefits

- Two-strap fastening for ease of application allows support to strap on like a sock
- Universal left/right application

Indications for Use

• Mild sprains of the ankle requiring minimal support and compression

Ordering Information

Size	Shoe Size (UK)	Product Code
Small	5 - 6	ASP/CX/BE/S
Medium	7 - 8	ASP/CX/BE/M
Large	9 - 10	ASP/CX/BE/L
X-Large	11 - 12	ASP/CX/BE/XL



Airtex TranspirableAnkle Brace

Features and Benefits

- Figure-of-eight straps
- Front opening with hook and loop fasteners
- Fixed to the sole of the foot with a seam
- 100% Latex free

Indications for Use

- Minor sprains
- Prevention of sprains while playing sport
- Tendonitis

Ordering Information

Size	Ankle Circ. (cm)	Product Code
Small	20 - 22	T114-S
Medium	22 - 24	T114-M
Large	24 - 26	T114-L
X-Large	26 - 28	T114-XL





Features and Benefits

- Provides compression around the ankle
- Improves the sense of stability of the joint
- Compression can be adjusted as and when required
- Can be worn with any shoe
- Padding around vulnerable areas of the ankle for added comfort
- Easy to fit

Indications for Use

- Minor distortions
- Slight capsule irritation
- Slight ankle instability
- Arthrosis without instability

Size	Amble Cine (am)	Product Code		
Size	Ankle Circ. (cm)	Left	Right	
1	26 - 29	PUS015B	PUS010B	
2	29 - 32	PUS016B	PUS011B	
3	32 - 35	PUS017B	PUS012B	
4	35 - 38	PUS018B	PUS013B	
5	38 - 41	PUS019B	PUS014B	



Features and Benefits

- Provides long-term support of the ankle joint
- Based on the functions of tape bandage
- The degree of pressure and pronation/ supination limits is adjustable
- Silicone application ensuring good position retention
- Comfort zones located on the most vulnerable areas of the ankle such as the Achilles tendon
- Made of Sympress™ which is a high quality material made for ultimate comfort
- Microfibres keep the skin dry and ensure optimal comfortability



- · Moderate residual instability
- Treatment of acute lateral ankle ligament lesions (alternative to taping)
- Treatment of stable ankle fractures
- Aftercare for ankle fractures
- Postoperative treatment of ligament injuries and fractures
 - Capsular irritation and tendonitis following ankle fracture or ligament lesions
- Arthrosis with instability
- Secondary prevention of ankle ligament lesions (under moderate stress with accumulation of fluid)
- Impingement

Ordering Information

Size	Ankla Civa (ama)	Product Code	
512e	Ankle Circ. (cm)	Left	Right
1	26 - 29	PUS022B	PUS0017B
2	29 - 32	PUS023B	PUSO018B
3	32 - 35	PUS024B	PUSO019B
4	35 - 38	PUS025B	PUS020B
5	38 - 41	PUS026B	PUS021B



Push Aequi Ankle Brace

Features and Benefits

- An effective ankle orthosis which strikes a balance between efficient mechanical support of a disabled ankle and good user comfort
- The rigid support is situated exclusively on the medial side of the ankle and has a recess for the malleolus
- The lateral side is made from preformed foam, which easily adapts to individual anatomy

Indications for Use

- Serious (residual) instability (under high stress)
- Treatment of acute lateral ankle ligament lesions
- After-care of conservative or operated ankle fractures
- Arthrosis with instability without accumulation of fluid
- Secondary prevention of ankle ligament lesions

Note: Paediatric sizes also available, see Page 228.

Size	Showh Hool Cive (am)	Product Code	
Size	Short Heel Circ. (cm)	Left	Right
1	27 - 31	320111	320121
2	31 - 34	320112	320122
3	34 - 40	320113	320123



Malleo Dynastab[®] Boa[®] **Ankle Brace**

Features and Benefits

- Offers a high level of support
- Lightweight
- Low profile meaning it can be worn in any type of shoe
- Additional straps and a metallic splint for maximum stabilisation
- Anatomic design
- Micro-aerated 3D fabric ensure optimum comfort
- Quick and easy to fit
- Boa® Lacing System

Indications for Use

- Mild to moderate sprain
- · Resumption of sporting activity e.g. basketball, volleyball, handball, etc.



Formfit® Ankle Brace with Figure 8

Features and Benefits

- · Provides stability and protection for the
- Low profile design fits easily into shoes
- Quick lace up feature makes is easy to don and doff
- Breathable material
- Inversion/eversion control

Indications for Use

- Ankle injury rehabilitation/ immobilisation (Grade I/II sprains)
- Support chronic instability
- · Prophylactic use

Ordering Information

Size	Ankle Circ. (cm)	Product Code
1	19 - 22	23510120100199
2	22 - 26	23510120200199
3	26 - 30	23510120300199

Size	Ankle Circ. (cm)	Product Code	
X-Small	28 - 30.5	B-212000001	
Small	30.5 - 33	B-212000002	
Medium	33 - 35.5	B-212000003	
Large	35.5 - 38	B-212000004	
X-Large	38 - 40.5	B-212000005	

A60™ Ankle Support

Features and Benefits

- Incorporating a stabilizer moulded at a 60° angle and breathable material
- Lightweight anatomical design
- Fits easily into athletic footwear to provide cool, comfortable protection
- · Helps prevent sprains and roll-over
- Anatomically designed for prophylactic support
- Single strap for simplified application
- Replaces time-consuming lacing
- Available in black

Indications for Use

- Chronic mild instability
- Prophylaxis



Airsport™

Features and Benefits

- Provides Air-Stirrup® support for inversion and eversion protection
- Foam-filled Aircell technology
- Cushions and supports the ankle
- Anterior Talofibular (ATF) cross-strap and integral forefoot and shin wraps
- Provide additional ligament stabilisation and compression
- Strong, coated fabric
- Durable and long-lasting
- Step-in style and automatic heel adjustment
- Easy on/off application and fast adjustments
- Universal fit
- Fits left or right foot

Indications for Use

• Grade I or II ankle sprain, chronic instability

Ordering Information

Size	Men's Shoe	Women's Shoe Size (UK)	Product Code	
Size	Size (UK)		Left	Right
Small	<7	<8	DJO-02TSL	DJO-02TSR
Medium	7.5 - 11.5	9 - 13	DJO-02TML	DJO-02TMR
Large	12+	13.5+	DJO-02TLL	DJO-02TLR

Size	Men's Shoe	Women's Shoe	Product Code	
Size	Size (UK)	Size (UK)	Left	Right
X-Small	<4	<3	02MXSL	02MXSR
Small	4.5 - 6	3.5 - 6.5	02MSL	02MSR
Medium	6.5 - 10.5	7 - 10	02MML	02MMR
Large	11 - 12.5	10.5 - 12	02MLL	02MLR
X-Large	13+	12.5+	02MXLL	02MXLR

104



AirLift™ PTTD

Features and Benefits

- · Versatile stirrup design
- Anatomically correct design
- New intuitive Integrated Inflation System
- Aircast® technology
- Durable cross strapping
- Easy step-in design
- Off-the-shelf customised fit for comfortable relief from PTTD symptoms

Indications for Use

- Posterior tibial tendon dysfunction (PTTD)
- Adult acquired flatfoot (AAF)



Rebound Foot-Up®

Features and Benefits

- Lightweight
- Contoured to the shape of the leg
- Comfortable to wear
- Easy to use
- Designed to provide visible improvement in gait by supporting the foot the moment it is raised

Indications for Use

• Foot drop and/or other similar conditions



Ordering Information

Rebound Foot-Up				
Size	Ankle Circ. (cm)	Product Code		
Small/Medium	18 - 25	RFU50001		
Large/X-Large	26 - 33	RFU50002		

Foot Wrap				
Size	Foot Circ. (cm)	Product Code		
Universal	20 - 28	RFU50101		

Rebound Foot-Up & Foot Wrap Bundle				
Size	Ankle Circ. (cm)	Product Code		
Small/Medium	18 - 25	RFU500011		
Large/X-Large	26 - 33	RFU500021		

Accessories			
Description	Contents	Product Code	
Attachment Kit	Plastic inlay, elastic strap, 1x nut & screw	RFU50201	
Nuts & Bolts	5x Nuts, 2 xShort & 3x Long Screws	RFU50202	

Size	Men's Shoe	Women's Shoe	Product Code	
3126	Size (UK)	Size (UK)	Left	Right
Small	<6	<6.5	02PSL-1	02PSR-1
Medium	6 - 10.5	7 - 10	02PML-1	02PMR-1
Large	11+	10+	02PLL-1	02PLR-1



Foot-Up shown with shoeless wrap

Features and Benefits

- Foot-Up is a lightweight ankle orthosis that offers dynamic support for drop foot or similar complaints
- The revolutionary orthosis provides visible improvement in gait by supporting the foot the moment it is raised
- Lightweight and difficult to see under
- Constructed from a padded breathable material, it is comfortable to wear with or without footwear (shoeless wrap required for use without footwear)

Indications for Use

- · Peroneal neuropathies
- Drop foot
- Dysfunctional gait

Styles Available





Foot-Up shown with footwear



Ordering Information

Siz		Ankla Civa (ana)	Product Code	
512	:e	Ankle Circ. (cm)	Black	Beige
1	Medium	18 - 21	07810-1	07811-1
2	Large	22 - 26	07810-2	07811-2
3	X-Large	27 - 33	07810-3	07811-3

Shoeless Wrap*				
Size		Fred Circ (core)	F	Product Code
312	:e	Foot Circ. (cm)	Black	Beige
1	Small	<20	27810-1	27811-1
2	Medium	20 - 24	27810-2	27811-2
3	Large	24 - 28	27810-3	27811-3
4	X-Large	>28	27810-4	27811-4

To order contact order@steepergroup.com

Plastic Insert	
Description	Product Code
Black	17810
Beige	17811

Foot-Up with Shoeless Wrap Bundle			
Description Product Code			
Black	07810-		
Beige	07811-		

= Size of Foot-Up (1	- 3)
= Size of Shoeless V	Vrap (1 - 4)

Example: To order a black size 1 (medium) Foot-Up with a size 2 (medium) Shoeless Wrap, the code would be **07810-12.**

Note: *Shoeless Wrap comprises of foot wrap only and must be ordered with complete Foot-Up kit unless the patient already has one.

106





PRIM Dorsi Assist

Features and Benefits

- · Consists of two separated pieces, a supramalleolar strap over the ankle and a non-slip support union to fix between the tongue and shoelaces by means of a central anchorage hook type. Both parts are joined by an elastic strap, which allows variable amounts of dorsiflexion
- On the back of the ankle band, there is a gel pad which is cushioned for increased comfort on the achilles tendon
- Discrete and flexible, it takes up no space in the shoe making it easy to fit
- To be fitted when lacing the shoe
- Has the function to produce dorsiflexion on the take-off phase of gait, indicated in flaccid palsy

Indications for Use

- Drop foot
- Peroneal neuropathies
- Dysfuntional gait

Textile Drop Foot Orthosis				
Size Measurement (cm) Product Code				
X-Small	13 - 17	ATX01-XS		
Small	17 - 21	ATX01-S		
Medium	21 - 25	ATX01-M		
Large	25 - 29	ATX01-L		



Plantar Band				
Size	Measurement	F	Product Code	
512 e	(cm)	Left	Right	
Small	17 - 21	ATX11-SL	ATX11-SR	
Medium	21 - 25	ATX11-ML	ATX11-MR	
Large	25 - 29	ATX11-LL	ATX11-LR	







Accessories	
Description	Product Code
Shoe Support	ATX12
Elastics Small	ATX13-S
Elastics Medium	ATX13-M
Elastics Large	ATX13-L
Metal Hooks (2 Pack)	ATX14

HEAD & NECK



Neurodyn Foot Lift Orthosis

Features and Benefits

- Fabric orthosis for dynamic compensation of flaccid paralysis
- · Orthopaedic rehabilitation should compensate for lost functionality without hindering remaining functions
- · Dynamic instead of rigid splints comfort instead of pressure
- Free plantar flexion

Indications for Use

- Flaccid paralysis
- Reduced dorsiflexion

Styles Available



Ordering Information

Black			
Size	Circ. of Leg Above Ankle (cm)	Left	Right
1	<19	070771/L	070771/R
2	19 - 21	070772/L	070772/R
3	21 - 23	070773/L	070773/R
4	23 - 25	070774/L	070774/R
5	25 - 27	070775/L	070775/R

Beige	Beige			
Size	Circ. of Leg Above Ankle (cm)	Left	Right	
1	<19	070771/L/BG	070771/R/BG	
2	19 - 21	070772/L/BG	070772/R/BG	
3	21 - 23	070773/L/BG	070773/R/BG	
4	23 - 25	070774/L/BG	070774/R/BG	
5	25 - 27	070775/L/BG	070775/R/BG	





Push AFO

Features and Benefits

- Allows flexion and extension of ankle joint during stance phase
- Adjustable elastic strap allowing for adaptable dorsiflexion assist during swing phase
- Can be used in conjunction with insoles inside brace
- Allows for a natural gait
- Lifts the foot during the swing phase, controls the plantar flexion post heel strike and then allows unimpeded rolling of the foot

Indications for Use

- Loss of control of the ankle dorsiflexors during swing phase and during loading response (heel strike/foot flat), possibly with ankle instability
- Disorders possibly involving loss of the function of the ankle dorsiflexors are: CVA, multiple sclerosis, postpolio syndrome, muscular dystrophy, spinal cord injury, polyneuropathy



Size	Shoe Size	Short Heel	Produc	ct Code
Size (E	(EU)	Size (cm)	Left	Right
1	36 - 41	27 - 31	PUS115B	PUS118B
2	39 - 44	31 - 34	PUS116B	PUS119B
3	42 - 47	34 - 40	PUS117B	PUS120B



Swedish Type Lightweight AFO

Features and Benefits

- Heat mouldable splint providing dorsiflexion assistance and lateral stability for ankle/foot area
- Easily trimmed and low profile with open heel to fit into a shoe easily
- Wide calf strap with touch/close
- Removable/washable padding for comfort and hygiene

Indications for Use

· Flaccid foot drop after neurological impairment



Swedish Type Lightweight AFO - Carbon Reinforced

Features and Benefits

- Heat mouldable splint providing dorsiflexion assistance
- · Easily trimmed with low profile and open heel to fit into a shoe easily
- Wide calf strap with touch and close fastening and removable/washable padding for comfort and hygiene

Indications for Use

 Flaccid foot drop after neurological impairment

Ordering Information

Size	Heel - Calf	Footplate	Product Code	
3126	Height (cm)	Length (cm)	Left	Right
Small	24	26	AFO-1S-L	AFO-1S-R
Medium	29	27.5	AFO-2M-L	AFO-2M-R
Large	33	28	AFO-3L-L	AFO-3L-R
X-Large	38	30	AFO-4XL-L	AFO-4XL-R

More Information Wear and Care Instructions

Ordering Information

Size	Heel - Calf Footp	Footplate	Produ	ct Code
Size	Height (cm)	Length (cm)	Left	Right
Small	24	26	AFO-1S-LC	AFO-1S-RC
Medium	29	27.5	AFO-2M-LC	AFO-2M-RC
Large	33	28	AFO-3L-LC	AFO-3L-RC
X-Large	38	30	AFO-4XL-LC	AFO-4XL-RC



PAEDIATRIC

ARM & SHOULDER

ANKLE & FOOT

PAEDIATRIC



Leaf Spring Orthosis

Features and Benefits

- Heat mouldable splint providing dorsiflexion assistance and limited plantar flexion
- Easily trimmed to fit into a shoe and with padded touch and close fastening
- ³/₄ and full foot versions are available

Indications for Use

· Flaccid foot drop

Ordering Information

3/4 Foot						
Size	Shoe Size	Heel - Calf Height			t Code	
3126	(UK)	(cm)	(cm)	Left	Right	
Small	2.5 - 5.5	33	14.5	PLSC-100-L	PLSC-100-R	
Medium	5.5 - 8.5	35.5	16	PLSC-200-L	PLSC-200-R	
Large	8.5 - 10.5	35.5	18.5	PLSC-300-L	PLSC-300-R	
X-Large	10.5 - 12.5	38	19.5	PLSC-400-L	PLSC-400-R	

Full Foot						
C:	Shoe Size	Heel - Calf Height	Footplate	Product Code		
Size	(UK)	(cm)	Length (cm)	Left	Right	
Small	2.5 - 5.5	33	21.5	PLS-FF-100-L	PLS-FF-100-R	
Medium	5.5 - 8.5	34.5	25	PLS-FF-200-L	PLS-FF-200-R	
Large	8.5 - 10.5	34.5	28	PLS-FF-300-L	PLS-FF-300-R	
X-Large	10.5 - 12.5	38	29	PLS-FF-400-L	PLS-FF-400-R	







AFO Light

Features and Benefits

- Cosmetically-appealing look
- Uses Flex-Foot® technology
- Made from carbon fibre
- Strong and lightweight solution

Indications for Use

- Drop foot secondary to:
 - Cerebro vascular accident (CVA)
 - Other neurological conditions
 - Multiple sclerosis
 - Mild knee instability
 - Charcot Marie Tooth disease (HMSN)
 - Lower extremity overuse injuries

Note: Please ensure you are satisfied with this product before applying the padding kit as it can only be returned in its original condition.

Size	UK Shoe	UK Shoe Size	Calf Height	Product Code	
Size	Size (Mens)	(Womens)	(cm)	Left	Right
Small	<7	4.5 - 6.5	32 - 36	O-121022	O-121012
Medium	7.5 - 10	7 - 9.5	34 - 38	O-121024	O-121014
Large	10.5 - 13.5	-	36 - 40	O-121026	O-121016
X-Large	11.5 - 14.5	-	38 - 42	O-121028	O-121018

Accessories	Product Code
AFO Light Refurb Kit	O-122101



Features and Benefits

- Limits the plantar flexion in the upper ankle joint
- · Magnetic closure simplifies opening and closing with one hand
- Great energy storage and return by leveraging carbon stiffness
- Can be trimmed to size with sturdy scissors
- Highly cosmetic



- Peroneal paresis
- · Weakness of dorsiflexion of the foot
- Drop foot
- Hemiplegia



Features and Benefits

PeroSupport.pro

- Dynamic re-adjustment of the foot
- Restriction of plantar flexion and support of dorsal flexion
- · Prevents lateral ankle twisting (supination position)
- Prevents side-slips (lateralisation)
- Great energy storage and return by leveraging carbon stiffness
- Can be trimmed to size with sturdy scissors



- Weakness of dorsiflexion of the foot
- Malalignment of the ankle joint as a result of weakness of dorsiflexion of the foot
- Hemiplegia



Ordering Information

Size Shoe Size		Shoe Size	Height	Produ	Product Code	
Size	(UK) (EU) (cm)		Left	Right		
X-Small	2 - 3.5	34 - 36	32.5	197R4/L0	197R4/R0	
Small	3.5 - 55	36 - 39	33.5	197R4/L1	197R4/R1	
Medium	5.5 - 8	39 - 42	34.5	197R4/L2	197R4/R2	
Large	8 - 10.5	42 - 45	35.5	197R4/L3	197R4/R3	
X-Large	10.5 - 13	45 - 48	37.5	197R4/L4	197R4/R4	

Spare Accessories		Product Code
3x pads, 1x magnetic fastener and 10x hook and loop dots	Universal	197R10

Ordering Information

Size	Shoe Size Shoe Size Height		Product Code		
Size	(UK)	(EU) (cm)	Left	Right	
Small	3.5 - 55	36 - 39	34	197R6/L1	197R6/R1
Medium	5.5 - 8	39 - 42	35	197R6/L2	197R6/R2
Large	8 - 10.5	42 - 45	35	197R6/L3	197R6/R3
X-Large	10.5 - 13	45 - 48	36	197R6/L4	197R6/R4

Spare Accessories	Shoe Size (UK)	Shoe Size (EU)	Product Code
Calf pad with fastener, 12 hook and	3.5 - 8	36 - 42	197R7/1
loop dots	8 - 13	42 - 48	197R7/2

116

ADULT INSOLES

ADULT INSOLES

HEAD & NECK



ToeOFF® Short

Features and Benefits

- Specifically developed to accommodate individuals with shorter stature, slightly broader-shaped tibia, and wider foot at the metatarsal area
- ToeOFF® Short has a heel height that is 7mm in all sizes

Indications for Use

- Designed to manage footdrop in conditions such as
 - Stroke
 - Multiple Sclerosis
 - Post Polio Syndrome
 - Muscular Dystrophy
 - Spinal Cord Injuries
 - Traumatic Brain Injuries
 - Guillian-Barre Syndrome
 - Charcot-Marie-Tooth

 - Myelomeningocele
 - Neuropathy
 - Cerebral Palsy
- Supports gait in conditions such as Posterior Tibialis Tendon Dysfunction (PTTD) and toe amoutations

Ordering Information

www.steepergroup.com

Size	Shoe Size	Shoe Size	Footplate Length	Height	Produc	ct Code
Size	(UK)	(EU)	(cm)	(cm)	Left	Right
Small	1 - 2.5	32 - 35	21.5	29.5	289901011	289902011
Medium	2.5 - 5	35 - 38	23.5	32	289901012	289902012
Large	5 - 8	38 - 42	25.5	34	289901013	289902013





ToeOFF® 2.0

Features and Benefits

- The ToeOFF® 2.0 is the updated version of the ToeOFF® Original
- Lightweight, carbon reinforced rehabilitation orthosis that absorbs force on heel strike
- · Worn at the anterior side of the leg, bridging the ankle and preventing plantar flexion
- SoftKIT™ 2.0 comprises padding shaped to precisely fit the orthosis and comes supplied with two tibia pads
- ComfortKIT™ 2.0 comprises thicker padding, designed for the more sensitive skin, fits the orthosis precisely
- Durable and strong straps
- · Available in black

Indications for Use

- Mild to severe footdrop
- · Mild to moderate ankle instability
- Impaired knee control
- Limb proprioception deficit
- Mild to moderate spasticity
- Partial foot amoutation

Note: To order ToeOFF® Original please contact Customer Services.

Size Shoe Size		Footplate	Height	Product Code		
Size	(UK) Length (cm) (cm)		Left	Right		
X-Small	1 - 3	21	36	289201010	289202010	
Small	3 - 4	23	38	289201011	289202011	
Medium	5 - 9	24.5	40.5	289201012	289202012	
Large	10 - 11	27	43	289201013	289202013	
X-Large	12 - 13	28.5	45	289201014	289202014	

Size		Product Code						
	SoftKIT™ 2.0*	ComfortKIT™ 2.0*	Replacement Straps					
X-Small	287500010	287510010	287600010					
Small	287500011	287510011	287600011					
Medium	287500012	287510012	287600012					
Large	287500013	287510013	287600013					
X-Large	287500014	287510014	287600014					

^{*}This accessory comes in a pair.

SPINE & TORSO

ToeOFF® 2½



Features and Benefits

- Features the same dynamic functions and stability as ToeOFF® 2.0 but has a lower heel height of 7mm, to adapt to modern shoes
- Lower toe lift which results in more space in the toe box
- MikroFIX pre-applied to simplify the padding application
- Durable and strong straps
- The straps can be applied on the interior or external side of the AFO

Indications for Use

- Designed to manage footdrop in conditions such as:
 - Stroke
 - Multiple Sclerosis
- Post Polio Syndrome
- Muscular Dystrophy
- Spinal Cord Injuries
- Traumatic Brain Injuries
- Guillian-Barre Syndrome
- Charcot-Marie-Tooth
- Myelomeningocele
- Neuropathy
- Cerebral Palsy
- Supports gait in conditions such as Posterior Tibialis Tendon Dysfunction (PTTD) and toe amputations
- · Can also be used for partial foot amputations, most proximal level is Chopart

Ordering Information

Size	Shoe Size	Shoe Size	Footplate	Height	Produc	t Code
Size	(UK)	(EU)	Length (cm)	(cm)	Left	Right
X-Small	1 - 2.5	32 - 35	21.5	36	289801010	289802010
Small	2.5 - 5	35 - 38	23	38	289801011	289802011
Medium	5 - 8	38 - 42	24.5	40.5	289801012	289802012
Large	8 - 10.5	42 - 45	27	43	289801013	289802013
X-Large	> 0.5	>45	28.5	45	289801014	289802014





Features and Benefits

- ½ heel height of 7mm in all sizes
- Smoother transition (flow throughout the gait cycle
- Increased ROM in sagittal plane
- Accommodates lower shoe heel heights
- Accommodates more gentle contours of shoe insoles
- More clearance for forefoot in the shoe toe box
- Easier to customise to meet alignment and biomechanical needs
- MikroFIX pre-applied which makes it quicker to attach a SoftKIT™ 2.0, ComfortKIT™ 2.0 or CoverKIT™ 2.0
- Closure straps are the same as for ToeOFF® 2.0
- · Lightweight and thin
- Open heel

Indications for Use

- · Mild proximal deficits
- Mild spasticity
- Limb proprioception deficit
- · Functions well for higher activity levels

Ordering Information

Size	Shoe Size	Shoe Size	Footplate Height		Produc	t Code
Size	(UK)	(EU)		Left	Right	
X-Small	1 - 2.5	32 - 35	21.5	36	289601010	289602010
Small	2.5 - 5	35 - 38	23	38	289601011	289602011
Medium	5 - 8	38 - 42	24.5	40.5	289601012	289602012
Large	8 - 10.5	42 - 45	27	43	289601013	289602013
X-Large	>10.5	>45	28.5	45	289601014	289602014

ANKLE & FOOT



BlueROCKER® Short

Features and Benefits

- Specifically developed to accommodate individuals with a shorter stature, slightly broader-shaped tibia, and wider foot at the metatarsal area
- Tibia padding should always be worn for patient comfort
- Both models have a heel height of 7mm in all sizes

Indications for Use

- Foot drop in combination with weak plantar flexors, and in some cases limited knee control
- May also be used for bilateral problems and partial foot amputations, most proximal level is Chopart





BlueROCKER® 2.0

Features and Benefits

- The BlueROCKER® 2.0 is the updated version of BlueROCKER® Original
- Identical in shape to the ToeOFF® 2.0 but offers more rigid orthotic control
- Developed primarily for bilateral patients and those with more involved pathologies
- The extra stability will improve both balance and posture and give the user greater security, especially individuals with weak quadricep muscles
- Usually the preferred orthosis to be used in conjunction with a socket and toe filler for management of partial foot amputations
- SoftKIT™ 2.0 comprises padding shaped to precisely fit the orthosis and comes supplied with two tibia pads
- ComfortKIT™ 2.0 comprises thicker padding, designed for the more sensitive skin, fits the orthosis precisely

Indications for Use

- Footdrop
- Weakness or impairment in multiple leg muscle groups
- · Impaired balance
- · Limb proprioception
- · Partial foot amputation

Ordering Information

Size	Shoe Size	Shoe Size	Footplate Length	Height	Product Code	
Size	(UK)	(EU)	(cm)	(cm)	Left	Right
Small	1 - 2.5	32 - 35	21.5	29.5	289951011	289952011
Medium	2.5 - 5	35 - 38	23.5	32	289951012	289952012
Large	5 - 8	38 - 42	25.5	34	289951013	289952013

Ordering Information

Size	Footplate	Hoight (cm)	Product Code		
Size	Length (cm)	Height (cm)	Left	Right	
X-Small	21	36	289401010	289402010	
Small	23	38	289401011	289402011	
Medium	24.5	40.5	289401012	289402012	
Large	27	43	289401013	289402013	
X-Large	28.5	45	289401014	289402014	

To order contact order@steepergroup.com



BlueROCKER® 2½

Features and Benefits

- The BlueROCKER® 2 ½ has the same dynamic functions and stability as a BlueROCKER® 2.0 with the exception that the 2½ has a lower heel height. 7mm, to adapt to modern shoes
- The toe lift is also slightly lower
- Spacious toe box
- Has all MikroFIX[™] pre-applied to simplify the fitting of the padding
- Durable and strong straps
- The straps can be applied for either inside or outside closure on the orthosis
- · Adaptable to the heel height of the shoe

Indications for Use

- BlueROCKER® is designed to manage footdrop in conditions such as:
 - Stroke
 - Multiple Sclerosis
 - Post Polio Syndrome
 - Muscular Dystrophy
 - Spinal Cord Injuries

 - Traumatic Brain Injuries
 - Guillian-Barre Syndrome
 - Charcot-Marie-Tooth

 - Myelomeningocele
 - Neuropathy
 - Cerebral Palsy
- Supports gait in conditions such as Posterior Tibialis Tendon Dysfunction (PTTD) and toe amputations
- Partial foot amoutations, most proximal level is Chopart

Ordering Information

Size	Shoe Shoe Footplate Height		Product Code			
Size	(UK)	(EU)		Left	Right	
X-Small	1 - 2.5	32 - 35	21.5	36	289701010	289702010
Small	2.5 - 5	35 - 38	23	38	289701011	289702011
Medium	5 - 8	38 - 42	24.5	40.5	289701012	289702012
Large	8 - 10.5	42 - 45	27	43	289701013	289702013
X-Large	>10.5	>45	28.5	45	289701014	289702014



Comfort Night Splint

Features and Benefits

- Adjustable flexion straps provide a gentle stretching of the plantar fascia and achilles tendon
- Soft elastane lined foam interior is soft, washable and breathable to eliminate heat build-up
- Elastane lined foam padded straps maintain comfort across the flexor tendons
- The orthosis helps to prevent pain and further injury by limiting contractions of the soft tissues at night
- Lightweight, low profile shell is sturdy and breathable
- Toe wedge to increase the plantar fascia stretch
- Ambidextrous fit left/right

Indications for Use

- Plantar fasciitis
- · Achilles tendonitis
- Plantar flexion contractures

Ordering Information

Size	Shoe Size (UK)	Product Code
Small	<6	UFO/S/E
Medium	6 - 9	UFO/M/E
Large	9+	UFO/L/E

To order contact order@steepergroup.com



More Information

PediSupport.night Features and Benefits lower lea

- Provides optimal immobilisation of the
- · Reduces a plantar flexed position of the foot and helps prevent contracture of the achillies tendon
- Maintains a plantar grade position
- Toe wedge to increase the plantar fascia
- Additional calf pad to offload the heel for patients in a supine position



- · Plantar fasciitis
- Postoperative immobilisation of the lower lea
- Dorsal flexor resp. peroneal paralysis
- For night postural support
- · Pes equinus position





DynaPro® AFO

Features and Benefits

- The semi-rigid insert can be positioned to accommodate the angle of the ankle/ foot combines with a removable pull strap to initiate a therapeutic dorsiflexion stretch to help combat mild plantar flexion contracture with or without inversion/eversion
- The cover has a non-skid surface sewn into the bottom to improve grip during transfers
- Toe post to provide protection from shear forces during in bed usage

Indications for Use

- Plantar flexion contractures of the ankle
- To protect the heel from unwanted pressure while in bed

Note: Paediatric sizes also available, see Page 232.



126

Size	Show Size (UK)	Shoe Size (EU)	Product Code
Small	3.5 - 5.5	36 - 39	190R1/S
Medium	5.5 - 7	39 - 41	190R1/M
Large	7 - 9.5	41 - 44	190R1/L
X-Large	9.5 - 11	44 - 46	190R1/XL

Ordering Information

Size	Length of Foot (cm)	Max Calf Circ. (cm)	Product Code
Small	15 - 20	38	DPAFO21S
Medium	20 - 25	42	DPAFO22M
Large	25 - 30	48	DPAFO23L

To order contact order@steepergroup.com



Features and Benefits

 Utilises Flex-E-Core technology for adjustment without the need for tools or heat

Indications for Use

- · Plantar fasciitis
- Achilles tendonitis
- · Night splinting
- Foot drop

Note: Paediatric sizes also available, see Page 233.

Ordering Information

Size	Max Foot Length (cm)	Brace Height (cm)	Product Code
Adult	25	35	LMU-A

Repose Wedge

Features and Benefits

• The Repose Wedge has been designed primarily for off-loading pressure from the heels and can be used wherever support or elevation of the lower legs is required. The Repose Wedge offers the same pressure reducing properties as other Repose products. Particularly suited to leg elevation to aid vascular conditions and a useful aid in posture support

Indications for Use

• The Repose Wedge is a versatile product suitable for a wide range of applications within pressure area

Ordering Information

Size	Description	Product Code
Universal	Wedge and pump	6711100
Universal	Wedge with PU strap and pump	6711111



Podus Boot

Features and Benefits

- A lightweight, quality plastic shell with poly/pile liner that helps treat and prevent lower extremity disorders associated with trauma or immobility
- An adjustable toe extension helps protect toes from bedding and the anti-rotation bar may be repositioned to assist with internal or external rotation
- Removable transfer pad on the dorsal side for limited patient ambulating
- Fits right or left foot
- Latex free

Indications for Use

- Pressure neurosis
- · Ankle contractures
- Post-op hip abduction/rotation

Ordering Information

Size	Mens Shoe Size (UK)	Womens Shoe Size (UK)	Product Code
Regular	<10	<11	POD-R
X-Large	>10.5	>11.5	POD-XL

Liners (3 Pack)				
Size	Product Code			
Regular	79-90560			
X-Large	79-90561			

To order contact order@steepergroup.com



PAEDIATRIC




A-Boot Universal


- Ambulation version (A-Boot) comes with removable, soft, non-slip sole
- shield against mould and mildew, stains, bacteria, fungus and algae
- removable and machine washable
- accommodate all foot sizes
- Rotator bar controls hip and leg rotation
- to prevent pressure areas through pressure reduction
- · Comes with extender kit to accommodate for wider or oedematous


Indications for Use

Note: Paediatric sizes also available, see

Features and Benefits







- · A toe extension is available to
- Boots position heel in floating position
- Universal sizing to reduce stock holding

- Wider or oedematous feet
- Ankle contractures
- Heel ulcerations

Page 231.

Ordering Information

Size	Product Code
Universal	A-BOOT





Universal PRAFO 650 Kodel

Features and Benefits

- Designed for recumbent or ambulant protection of the heel and foot
- Zero pressure or shear force at the heel
- Universal sizing to reduce stock holding
- Resists plantarflexion contractures developing
- Custom modifications are possible
- One version fits both left and right
- Calibrated adjustment of plantar/ dorsiflexion
- Integral walking base designed for ambulation
- Washable and replaceable kodel liner
- Adjustable foot, toe and calf sections
- · Aluminium heel connector for strength
- · Rotation bar resists internal/external rotation when used in bed
- · Suitable for day and night use

Indications for Use

- Pressure ulcers
- Can be used to resist the development of contracture

Note: Paediatric sizes also available, see Page 230.

Size	Product Code	Accessories	Product Code
Universal	PRFO	Universal liner	PRAFO-LINER


Universal PRAFO APU (Adjustable Posterior Upright) Features and Benefits

- Designed for recumbent or ambulant protection of the heel and foot
- Resists plantarflexion contractures developing
- Custom modifications are possible
- · One version fits both left and right
- Calibrated adjustment of plantar/ dorsiflexion
- Integral walking base designed for ambulation
- Washable and replaceable kodel liner
- Aluminium heel connector for strength
- Rotation bar resists internal/external
- · Universal sizing to reduce stock holding

Indications for Use

- · Pressure ulcers
- · Can be used to resist the development

Note: Paediatric sizes also available, see Page 230.



- Zero pressure or shear force at the heel

- · Adjustable foot, toe and calf sections
- rotation when used in bed

- of contracture

Ordering Information

Size	Product Code	Accessories	Product Code
Universal	PRFO/APU	Universal liner	PRAFO-LINER





- The pad and strap version allows bulky wound dressings to be accommodated and can allow superior foot control when the PRAFO® is to be used in walking training and other circumstances when mobility is the primary focus
- Universal sizing to reduce stock holding
- Bi-planar adjustments with a screwdriver
- Opportunity to evaluate and train gait
- Ambulatory or recumbent usage
- All the benefits of the 650 and 650 APU - with increased application possibilities
- Calibrated adjustment of plantar/ dorsiflexion
- Integral walking base designed for ambulation
- Washable and replaceable kodel liner svstem
- Adjustable foot, toe and calf sections
- Aluminium heel connector for strength
- · Rotation bar resists internal/external rotation when used in bed
- Universal liners are also available to order - use code PRAFO-LINER

Indications for Use

• The Dual Action accommodates and helps control varus or valgus conditions of the ankle/foot complex, as well as allowing the practitioner a more precise and simplified method of setting the ankle/foot complex in the sagital plane

Note: Paediatric sizes also available, see Page 230.

Size	Product Code	Accessories	Product Code
Universal	PRFO/DA	Universal liner	PRAFO-LINER

ARM & SHOULDER

ANKLE & FOOT



Formfit® Walker

Features and Benefits

- Low, wide rocker bottom promotes a stable, natural gait
- Shock-absorbing sole for reduced impact at heel strike
- · Breathable, seamless liner
- · Available in black

Indications for Use

- For foot and ankle conditions that may benefit from motion restriction such as:
 - Grade II or III ankle sprains
- Stable fractures
- Soft tissue injuries
- During postoperative or rehab phase



Size	Shoe Size (UK)	Product Code
X-Small	1.5 - 3.5	W0200BLK
Small	4 - 6.5	W0400BLK
Medium	7 - 10	W0600BLK
Large	10.5 - 12	W0800BLK
X-Large	12.5 - 14.5	W1000BLK





Features and Benefits

- The Rebound Air Walker features linked medial and lateral air bladders with a single inflation point to provide balanced compression. These air bladders are shaped to cradle the Achilles, holding the heel in place for extra stability
- Moisture wicking soft goods and ventilated panels prevent microbial build up and helps to keep the skin dry and comfortable
- Straps fasten in front for easy application, adjustment and removal
- EVA material provides greater shock absorption in the heel
- Fully-integrated air pump and release valve provides easy, tool-free inflation and deflation for fine-tuned, consistent compression
- Flex Edge[™] over-moulding conforms to each patient's unique anatomy to help limit the formation of pressure points
- Rocker bottom promotes a natural, stable gait, whilst a semi-soft toe bumper that can be trimmed as required provides added protection

Indications for Use

- Grade 2 and 3 sprains
- Stable fractures
- · Postoperative use
- Trauma
- · Rehabilitation use

Ordering Information

Size	Men's Shoe Size (UK)	Women's Shoe Size (UK)	Product Code
X-Small	2 - 4	3.5 - 5.5	B-242900001E
Small	4 - 7	5.5 - 8.5	B-242900002E
Medium	7 - 10.5	8.5 - 12	B-242900003E
Large	10.5 - 12.5	12 - 14	B-242900004E
X-Large	12.5 - 15	14 - 16.5	B-242900005E

134

ARM & SHOULDER

ANKLE & FOOT



Rebound® Diabetic Walker

Features and Benefits

- Uniquely designed for effective treatment of plantar foot ulcers, the Rebound Diabetic Walker enables easy access and re-fitting during wound care. whilst improving the patient's mobility during the healing process by including a rocker bottom outsole for smooth roll-over
- Effective tri-laminate insole with removable shock-absorbing hexagons that promote ulcer offloading and a Plastazote top layer that evenly distributes the pressure
- Anatomically designed shell that minimises shear forces by immobilising the foot and ankle at 90° whilst providing protection
- Flex Edge technology for relief of pressure points, along with a patentpending pressure relief valve to prevent over-inflation
- Positioning portals allow the clinician to check the fit of the walker at the heel bottom, in order to ensure appropriate positioning of the foot in the walker
- Soft and breathable machine washable liner for hygiene promotion and prevention of infection

Indications for Use

· Plantar foot ulcers

Ordering Information

Size	Men's Shoe Size (UK)	Women's Shoe Size (UK)	Product Code
Small	4.5 - 7	5.5 - 8.5	B-242500002
Medium	7 - 10.5	8.5 - 12	B-242500003
Large	10.5 - 12.5	12 - 14	B-242500004
X-Large	12.5 - 15	14 - 16.5	B-242500005





Features and Benefits

- Inspired by ski boots, the VACOped features a unique design to enable patients to weight-bear earlier, enabling a faster return to daily activities
- Inner vacuum cushion provides equal and firm support by adjusting itself precisely to the shape of every foot and lower leg and allowing adjustments to be made for swelling
- Adjustable ROM bar allows the boot to be fixed or permit movement dependent on the indication
- Durable, yet lightweight Germanengineered plastic shell provides the injury site with the necessary stability and allows an increased freedom of movement in comparison to a plaster
- Multiple washable and exchangeable liners promote hygiene

Indications for Use

- Achilles tendon ruptures
- Ankle fractures
- Metatarsal fractures
- Calcaneus fractures
- Ankle replacements
- Arthrodesis

Size	Shoe Size (UK)	Calf Circ. (cm)	Ankle Circ. (cm)	Product Code
Small	2.5 - 6.5	23 - 39	14 - 29	VPED-01-4-S
Medium	6.5 - 10.5	35 - 51	22 - 40	VPED-01-4-M
Large	10.5 - 14	40 - 58	30 - 45	VPED-01-4-L

ARM & SHOULDER

VACOcast Diabetic

Features and Benefits

- The VACOcast Diabetic provides the same stability and offloading as a TCC, but is easy to apply and comfortable for the patient
- VACOcast Diabetic will help you achieve cost savings and better outcomes

Indications for Use

- Diabetic Foot Syndrome
- Ulcerations on the foot
- Postoperative care after toe and forefoot amputations
- Diabetic neuro-osteoarthropathy
- DNOAP and Charcot foot

Note: If patient is in-between two sizes, we recommend choosing the bigger size.





EVENup Shoe Lift

Features and Benefits

- Designed to offset the imbalance often experienced whilst wearing a walker boot type orthosis
- Stretches over the shoe
- For inside and outside the house
- Easy to apply and easy to adjust

Indications for Use

• Severe hip or knee pain, e.g. when wearing a walker boot or cast

Ordering Information

Size	Shoe Size (UK)	Calf Circ. (cm)		Foot Width From Ball (cm)	Product Code
Small	2.5 - 5.5	23 - 39	14 - 29	10.5	VPCD-01-1-S
Medium	6 - 9	35 - 51	22 - 40	12	VPCD-01-1-M
Large	9.5 - 14	40 - 58	30 - 45	13.5	VPCD-01-1-L

Ordering Information

Size	Shoe Size (UK)	Product Code
Small	2 - 5	AUSG-01-2-S
Medium	5.5 - 9	AUSG-01-2-M
Large	9.5 - 11	AUSG-01-2-L

ANKLE & FOOT

138

Neoprene Boots

Features and Benefits

- 4mm neoprene rubber faced and lined with jersey nylon
- Touch and close fastening
- EVA unit with 8mm sole and 25mm heel
- Comes with heel stiffeners
- Comfort foam insole
- Available in black
- Wide fitting
- Only available in pairs

Indications for Use

· Bulky dressing and fluctuating oedema



Felt Boots

Features and Benefits

- Supplied with heel and toe stiffeners
- Touch and close fastening
- EVA unit with 8mm sole and 25mm heel
- Comfort foam insole
- Available in black
- Wide fitting
- Only available in pairs

Indications for Use

· Bulky dressing and fluctuating oedema

Ordering Information

Shoe Size (UK)	Product Code	Shoe Size (UK)	Product Code
2	SHNB/02	9	SHNB/09
3	SHNB/03	10	SHNB/10
4	SHNB/04	11	SHNB/11
5	SHNB/05	12	SHNB/12
6	SHNB/06	13	SHNB/13
7	SHNB/07	14	SHNB/14
8	SHNB/08		

Shoe Size (UK)	Product Code	Shoe Size (UK)	Product Code
2	SHFB/02	9	SHFB/09
3	SHFB/03	10	SHFB/10
4	SHFB/04	11	SHFB/11
5	SHFB/05	12	SHFB/12
6	SHFB/06	13	SHFB/13
7	SHFB/07	14	SHFB/14
8	SHFB/08		

ARM & SHOULDER

142

PAEDIATRIC

ANKLE & FOOT

All Purpose Boot



Features and Benefits

- Wider sole and deep opening provides increased room to accommodate large dressings, bandages and casts
- Wide ankle strap eliminates heel slippage as well as securing the foot or cast in the shoe
- A rocker bottom sole encourages smooth roll through gait pattern and reduces plantar pressure
- · Cubical cushioning sole for selective pressure relief

Indications for Use

- · Postoperative protection for soft tissue trauma or surgery
- Plantar ulcers/wounds
- Closed toe allows dressing to keep clean
- Wider sole allows more room for dressings and casts

DH Active Offloading Post-Op Shoe

Features and Benefits

- An excellent alternative to post-op style shoes with static pads
- Accommodates the ambulatory and compliant patient for short-term treatment (1-4 months)
- Can accommodate either the left or right
 - Designed to accommodate bulky
- Super soft nylon upper with reinforced heel counter provides extra comfort and stability
- Patented removable toe pieces aides in wound accessibility for daily wound site inspection

Indications for Use

- · Plantar wounds
- Postoperative use
- · Plantar foot ulcers

Ordering Information

Size	Shoe Size (UK)	Insole Width (cm)	Insole Length (cm)	Product Code
X-Small	Up to 4	8.5	21.5	NV0321
Small	4.5 - 6.5	9.5	23.5	NV0322
Medium	7 - 8	10	25.5	NV0323
Large	8.5 - 10	10.5	27.5	NV0324
X-Large	10.5 - 12.5	11.5	29.5	NV0326

More Information Wear and Care Instructions

Ordering Information

Size	Men's Shoe Size (UK)	Women's Shoe Size (UK)	Product Code
Small	3.5 - 5.5	4 - 6	10340
Medium	6 - 8	6.5 - 8.5	10341
Large	8.5 - 10.5	9 - 10	10342
X-Large	11 - 13.5	10.5+	10343

To order contact order@steepergroup.com

Komet Therapeutic Shoes

Features and Benefits

- Pressure redistribution and relief in the forefoot area through a rigid roll-off sole
- Benefits the healing process of wounds in the forefoot area
- Cubical cushioning sole for selective pressure relief
- Reinforced heel area and straps in the forefoot area give the foot support
- Adjustable hook and loop fasteners on the flaps accommodate large volume bandages
- Wide angular sole shape in the forefoot area when using surgical techniques with Kirschner wire
- Profiled outsole for high slip resistance
- Anti-bacterial inner shoe lining and cover fabric of the padded insole
- Weight limit: 120kg



Indications for Use

- Wounds in the forefoot area caused by diabetes mellitus
- Peripheral artery disease (PAD)
- Other indication independent wound formation in the foot and ankle area
- Post-surgical for corrections and forefoot injuries

Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Width (cm)	Product Code
2.5 - 3.5	35 - 36	24.2	10.3	217R1/36
4 - 5	37 - 38	25.5	10.6	217R1/38
5.5 - 6.5	39 - 40	26.8	10.9	217R1/40
7 - 8	41 - 42	28.2	11.2	217R1/42
9 - 9.5	43 - 44	29.5	11.5	217R1/44
10.5 - 11	45 - 46	30.8	11.8	217R1/46
12 - 13	47 - 48	32.2	12.1	217R1/48
13.5 - 14.5	49 - 50	33.5	12.4	217R1/50



Cubical Cushioning Insole					
Shoe Size Shoe Size Product (UK) (EU) Code					
2.5 - 3.5	35 - 36	217R11/36			
4 - 5	37 - 38	217R11/38			
5.5 - 6.5	39 - 40	217R11/40			
7 - 8	41 - 42	217R11/42			
9 - 9.5	43 - 44	217R11/44			
10.5 - 11	45 - 46	217R11/46			
12 - 13	47 - 48	217R11/48			
13.5 - 14.5	49 - 50	217R11/50			

Forefoot Protective Cap				
Shoe Size (UK)	Shoe Size (EU)	Product Code		
2.5 - 3.5	35 - 36	217R12/36		
4 - 5	37 - 38	217R12/38		
5.5 - 6.5	39 - 40	217R12/40		
7 - 8	41 - 42	217R12/42		
9 - 9.5	43 - 44	217R12/44		
10.5 - 11	45 - 46	217R12/46		
12 - 13	47 - 48	217R12/48		
13.5 - 14.5	49 - 50	217R12/50		



Foot pressure measurement -Left: Daytime shoe Right: Komet



144

Merkur Therapeutic Shoes



Features and Benefits

- Features a forward rising and shortened outsole (negative heel)
- Pressure relief and pressure distribution during walking
- Benefits the healing process of wounds in the forefoot area
- Encompasses the heel and tarsus with a softly padded socket
- Heel cap, hook and loop straps and fasteners prevent slipping in the shoe
- Cubical cushioning sole for selective pressure relief
- Improved slip protection by integrated anti-slip nubs
- Anti-bacterial inner lining and cover fabric of the padded insole



Indications for Use

- Wounds in the forefoot area caused by diabetes mellitus
- Peripheral artery disease (PAD)
- Other indication independent wound formation in the foot and ankle area
- Post-surgical for corrections and forefoot injuries

Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Width (cm)	Product Code
X-Small	12 - 1	31 - 33	22	8.7	201R1/XS
Small	2 - 3.5	34 - 36	24	9.1	201R1/S
Medium	4 - 5.5	37 - 39	26	9.6	201R1/M
Large	6.5 - 8	40 - 42	28	10	201R1/L
X-Large	9 - 10.5	43 - 45	30	10.5	201R1/XL
XX-Large	11 - 13	46 - 48	32	10.9	201R1/XXL



Cubical Cushioning Insole					
Size	Shoe Size (UK)	Shoe Size (EU)	Product Code		
X-Small	12 - 1	31 - 33	201R12/XS		
Small	2 - 3.5	34 - 36	201R12/S		
Medium	4 - 5.5	37 - 39	201R12/M		
Large	6.5 - 8	40 - 42	201R12/L		
X-Large	9 - 10.5	43 - 45	201R12/XL		
XX-Large	11 - 13	46 - 48	201R12/XXL		






Foot pressure measurement -Left: Balance shoe Right: Merkur

Forefoot Protective Cap					
Size	Shoe Size (UK)	Shoe Size (EU)	Product Code		
X-Small	12 - 1	31 - 33	201R10/XS		
Small	2 - 3.5	34 - 36	201R10/S		
Medium	4 - 5.5	37 - 39	201R10/M		
Large	6.5 - 8	40 - 42	201R10/L		
X-Large	9 - 10.5	43 - 45	201R10/XL		
XX-Large	11 - 13	46 - 48	201R10/XXL		

Venus Therapeutic Shoes

Features and Benefits

- Features a forward rising and shortened outsole (negative heel)
- Pressure relief and pressure distribution during walking
- Benefits the healing process of wounds in the forefoot area
- Encompasses the heel and tarsus with a softly padded socket
- Heel cap, hook and loop straps and fasteners prevent slipping in the shoe
- Anti-bacterial inner lining and cover fabric of the padded insole



Post-surgical for corrections and forefoot injuries



Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Product Code
Small	3.5 - 5	36 - 38	16	200R1/S
Medium	5.5 - 8	39 - 42	17.5	200R1/M
Large	9 - 11	43 - 46	19	200R1/L
X-Large	12 - 15	47 - 51	20.5	200R1/XL





Foot pressure measurement -Left: Daytime shoe Right: Venus

Adjustable Protection Bracket					
Size	Shoe Size (UK)	Shoe Size (EU)	Product Code		
Small	3.5 - 5	36 - 38	200R14/S		
Medium	5.5 - 8	39 - 42	200R14/M		
Large	9 - 11	43 - 46	200R14/L		
X-Large	12 - 15	47 - 51	200R14/XL		



Forefoot Protective Cap				
Size	Shoe Size (UK)	Shoe Size (EU)	Product Code	
Small	3.5 - 5	36 - 38	200R10/S	
Medium	5.5 - 8	39 - 42	200R10/M	
Large	9 - 11	43 - 46	200R10/L	
X-Large	12 - 15	47 - 51	200R10/XL	

148

Merkur A Therapeutic Shoes

Features and Benefits

- Features a rigid, slightly ascending outsole towards the front (negative heel)
- Pressure relief and pressure distribution during walking
- Benefits the healing process of wounds in the forefoot area
- Encompasses the entire foot and ankle area with a softly padded upper
- Hook and loop fasteners at the wide flap prevent slipping within the shoe
- Cubical cushioning sole for selective pressure relief
- The flap can be opened widely and together with the adjustable hook and loop fasteners, it provides sufficient space for bandaged feet
- Improved slip protection by integrated anti-slip nubs
- Anti-bacterial inner lining and cover fabric of the padded insole



Indications for Use

- Wounds in the foot and ankle area caused by diabetes mellitus
- Peripheral artery disease (PAD)
- Other indication independent wound formation in the foot and ankle area
- Postoperative and trauma
- Gait unsteadiness (dizziness, visual impairment)

Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Width (cm)	Product Code
2.5 - 3.5	35 - 36	23.8	9	202R1/36
4 - 5	37 - 38	25.1	9.2	202R1/38
5.5 - 6.5	39 - 40	26.4	9.5	202R1/40
7 - 8	41 - 42	27.7	9.7	202R1/42
9 - 9.5	43 - 44	29	10	202R1/44
10.5 - 11	45 - 46	30.3	10.2	202R1/46
12 - 13	47 - 48	31.6	10.5	202R1/48
13.5 - 14.5	49 - 50	33	11.2	202R1/50



Cubical Cushioning Insole				
Shoe Size (UK)	Shoe Size (EU)	Product Code		
2.5 - 3.5	35 - 36	202R12/36		
4 - 5	37 - 38	202R12/38		
5.5 - 6.5	39 - 40	202R12/40		
7 - 8	41 - 42	202R12/42		
9 - 9.5	43 - 44	202R12/44		
10.5 - 11	45 - 46	202R12/46		
12 - 13	47 - 48	202R12/48		
13.5 - 14.5	49 - 50	202R12/50		



Foot pressure measurement -Left: Balance shoe Right: Merkur A

150

ANKLE & FOOT

Mars Therapeutic Shoes



Features and Benefits

- Features a padded midfoot roll and a shortened sole construction towards the heel area
- Pressure relief and pressure distribution during walking
- Benefits the healing process of wounds in the heel area
- Encompasses the ankle and tarsus with a softly padded socket
- Hook and loop fasteners hold the foot in a predefined position and thus prevent slipping in the shoe
- The removable padded insole prevents pressure peaks in the forefoot area
- Anti-bacterial inner lining and cover fabric of the padded insole
- Weight limit: 120kg



Indications for Use

- Wounds in the heel area caused by diabetes mellitus
- Peripheral artery disease (PAD)
- Other indication independent wound formation in the foot and ankle area
- Postoperative and trauma

Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Width (cm)	Product Code
Small	2.5 - 4	35 - 37	22.5	8.7	203R1/S
Medium	5 - 6.5	38 - 40	23.5	9	203R1/M
Large	7 - 9	41 - 43	25.5	9.5	203R1/L
X-Large	9.5 - 11	44 - 46	27.5	10	203R1/XL



Forefoot Prote	ection Cap		
Size	Shoe Size (UK)	Shoe Size (EU)	Product Code
Small	2.5 - 4	35 - 37	203R10/S
Medium	5 - 6.5	38 - 40	203R10/M
Large	7 - 9	41 - 43	203R10/L
X-Large	9.5 - 11	44 - 46	203R10/XL



Heel Protection Cap					
Size	Shoe Size (UK)	Shoe Size (EU)	Product Code		
Small	2.5 - 4	35 - 37	203R11/S		
Medium	5 - 6.5	38 - 40	203R11/M		
Large	7 - 9	41 - 43	203R11/L		
X-Large	9.5 - 11	44 - 46	203R11/XL		



Foot pressure measurement -Left: Daytime shoe Right: Mars

152

ANKLE & FOOT

Uranus Therapeutic Shoes

Features and Benefits



- Features a rigid, slightly ascending outsole towards the front (negative heel)
- Pressure relief and pressure distribution during walking
- Benefits the healing process of wounds in the forefoot area
- Encompasses the entire foot and ankle area with a softly padded upper
- · Hook and loop fasteners at the wide flap prevent slipping within the shoe
- Cubical cushioning sole for selective pressure relief
- The flap can be opened widely and together with the individually adjustable hook and loop fasteners, it provides sufficient space for bandaged feet
- Improved slip protection by integrated anti-slip nubs
- · Anti-bacterial inner lining and cover fabric of the padded insole
- Weight limit: 120kg



- · Wounds in the foot and ankle area caused by diabetes mellitus
- Peripheral artery disease (PAD)
- · Other indication independent wound formation in the foot and ankle area
- Postoperative and trauma
- Gait unsteadiness (dizziness, visual impairment)



Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Width (cm)	Product Code
2.5 - 3.5	35 - 36	23.8	9	204R1/36
4 - 5	37 - 38	25.1	9.2	204R1/38
5.5 - 6.5	39 - 40	26.4	9.5	204R1/40
7 - 8	41 - 42	27.7	9.7	204R1/42
9 - 9.5	43 - 44	29	10	204R1/44
10.5 - 11	45 - 46	30.3	10.2	204R1/46
12 - 13	47 - 48	31.6	10.5	204R1/48
13.5 - 14.5	49 - 50	33	11.2	204R1/50



Cubical Cushioning Insole				
Shoe Size (UK)	Shoe Size (EU)	Product Code		
2.5 - 3.5	35 - 36	204R12/36		
4 - 5	37 - 38	204R12/38		
5.5 - 6.5	39 - 40	204R12/40		
7 - 8	41 - 42	204R12/42		
9 - 9.5	43 - 44	204R12/44		
10.5 - 11	45 - 46	204R12/46		
12 - 13	47 - 48	204R12/48		
13.5 - 14.5	49 - 50	204R12/50		



Foot pressure measurement -Left: Balance shoe Right: Uranus

Jupiter Therapeutic Shoes



Features and Benefits

- Features a forward rising and shortened outsole (negative heel)
- Pressure relief and pressure distribution during walking
- Benefits the healing process of wounds in the forefoot area
- Encompasses the heel and tarsus with a softly padded socket
- Heel cap, hook and loop straps and fasteners prevent slipping in the shoe
- Cubical cushioning sole for selective pressure relief
- Integrated anti-slip nubs
- Anti-bacterial inner lining and cover fabric of the padded insole
- Weight limit: 120kg



Indications for Use

- Wounds in the forefoot area caused by diabetes mellitus
- Peripheral artery disease (PAD)
- Other indication independent wound formation in the foot and ankle area
- Post-surgical for corrections and forefoot injuries

Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Width (cm)	Product Code
X-Small	12 - 1	31 - 33	22	8.7	205R1/XS
Small	2 - 3.5	34 - 36	24	9.1	205R1/S
Medium	4 - 5.5	37 - 39	26	9.6	205R1/M
Large	6.5 - 8	40 - 42	28	10	205R1/L
X-Large	9 - 10.5	43 - 45	30	10.5	205R1/XL
XX-Large	11 - 13	46 - 48	32	10.9	205R1/XXL
XXX-Large	13.5 - 14.5	49 - 50	34	11.5	205R1/XXXL





Foot pressure measurement -Left: Balance shoe Right: Jupiter

Forefoot Protective Cap					
Size	Shoe Size (UK)	Shoe Size (EU)	Product Code		
X-Small	12 - 1	31 - 33	205R10/XS		
Small	2 - 3.5	34 - 36	205R10/S		
Medium	4 - 5.5	37 - 39	205R10/M		
Large	6.5 - 8	40 - 42	205R10/L		
X-Large	9 - 10.5	43 - 45	205R10/XL		
XX-Large	11 - 13	46 - 48	205R10/XXL		
XXX-Large	13.5 - 14.5	49 - 50	205R10/XXXL		



Cubical Cushioning Insole					
Size	Shoe Size (UK)	Shoe Size (EU)	Product Code		
X-Small	12 - 1	31 - 33	205R12/XS		
Small	2 - 3.5	34 - 36	205R12/S		
Medium	4 - 5.5	37 - 39	205R12/M		
Large	6.5 - 8	40 - 42	205R12/L		
X-Large	9 - 10.5	43 - 45	205R12/XL		
XX-Large	11 - 13	46 - 48	205R12/XXL		
XXX-Large	13.5 - 14.5	49 - 50	205R12/XXXL		

Saturn Therapeutic Shoes

Features and Benefits

- Features an ascending outsole towards the front (negative heel)
- Pressure relief and pressure distribution during walking
- Benefits the healing process of wounds in the forefoot area
- Embraces heel and tarsus with a softly padded shoe upper
- Heel cap, hook and loop fasteners prevent slipping within the shoe
- Anti-bacterial inner lining and cover fabric of the padded insole
- Weight limit: 120kg



- Wounds in the forefoot area caused by diabetes mellitus
- Peripheral artery occlusions
- Other wounds in the forefoot area
- Post-surgical after corrections and forefoot injuries



Size	Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Width (cm)	Product Code
X-Small	3.5 - 4	36 - 37	24.5	9.6	207R1/XS
Small	5 - 5.5	38 - 39	25.8	9.8	207R1/S
Medium	6.5 - 7	40 - 41	27.1	10	207R1/M
Large	8 - 9	42 - 43	28.4	10.2	207R1/L
X-Large	9.5 - 10.5	44 - 45	29.7	10.4	207R1/XL



Forefoot Protective Cap					
Size	Shoe Size (UK)	Shoe Size (EU)	Product Code		
X-Small	3.5 - 4	36 - 37	207R10/XS		
Small	5 - 5.5	38 - 39	207R10/S		
Medium	6.5 - 7	40 - 41	207R10/M		
Large	8 - 9	42 - 43	207R10/L		
X-Large	9.5 - 10.5	44 - 45	207R10/XL		

Product Code

209R11/36

209R11/38

209R11/40 209R11/42 209R11/44 209R11/46 209R11/50

Apollo Therapeutic Shoes

Features and Benefits

- Follows the concept of zero heel technology with its pronounced roll-off sole and thus relieves the entire foot during the swing and stance phase
- Benefits the healing process of wounds in the forefoot area
- The removable padded insole prevents pressure peaks in the foot sole area
- Anti-bacterial inner lining and cover fabric of the padded insole
- Weight limit: 120kg



Indications for Use

 Postoperative for healing of surgical interventions and injuries to the forefoot, e.g. hallux valgus surgery

Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Width (cm)	Product Code
2.5 - 3.5	35 - 36	23.9	8.2	209R1/36
4 - 5	37 - 38	25.2	8.5	209R1/38
5.5 - 6.5	39 - 40	26.5	8.8	209R1/40
7 - 8	41 - 42	27.9	9.1	209R1/42
9 - 9.5	43 - 44	29.2	9.4	209R1/44
10.5 - 11	45 - 46	30.5	9.7	209R1/46
12 - 13	47 - 48	31.9	10	209R1/48
13.5 - 14.5	49 - 50	33.1	10.3	209R1/50



Forefoot Pr	refoot Protective Cap Cushioning Insole					
Shoe Size (UK)	Shoe Size (EU)	Product Code	Shoe Size (UK)	Shoe Size		
2.5 - 3.5	35 - 36	209R12/36	2.5 - 3.5	35 - 36		
4 - 5	37 - 38	209R12/38	4 - 5	37 - 38		
5.5 - 6.5	39 - 40	209R12/40	5.5 - 6.5	39 - 40		
7 - 8	41 - 42	209R12/42	7 - 8	41 - 42		
9 - 9.5	43 - 44	209R12/44	9 - 9.5	43 - 44		
10.5 - 11	45 - 46	209R12/46	10.5 - 11	45 - 46		
12 - 13	47 - 48	209R12/48	12 - 13	47 - 48		
13.5 - 14.5	49 - 50	209R12/50	13.5 - 14.5	49 - 50		



Foot pressure measurement -Left: Daytime shoe Right: Apollo

160

Luna Therapeutic Shoes

Features and Benefits

- Sandal-like, wide dressing shoe with raised toe edge
- Designed to accommodate large-volume dressings and protect them against slipping
- Cubical cushioning sole for selective pressure relief
- Higher socket and sturdily reinforced heel area
- Integrated anti-slip nubs
- Anti-bacterial inner lining and cover fabric of the padded insole
- Weight limit: 120kg



- · Wounds in the foot and ankle area caused by diabetes mellitus
- Peripheral artery disease (PAD)
- · Other indication independent wound formation in the foot and ankle area
- · Postoperative and trauma





Ordering Information

Shoe Size	Shoe Size	e Size Inner Shoe Width		idth Product Code	
(UK)	(EU)	Length (cm)	Length (cm) (cm)	Left	Right
2.5 - 3.5	35 - 36	23.9	9.9	210R1/36L	210R1/36R
4 - 5	37 - 38	25.2	10.2	210R1/38L	210R1/38R
5.5 - 6.5	39 - 40	26.5	10.5	210R1/40L	210R1/40R
7 - 8	41 - 42	27.9	10.8	210R1/42L	210R1/42R
9 - 9.5	43 - 44	29.2	11.1	210R1/44L	210R1/44R
10.5 - 11	45 - 46	30.5	11.4	210R1/46L	210R1/46R
12 - 13	47 - 48	31.9	11.7	210R1/48L	210R1/48R
13.5 - 14.5	49 - 50	33.1	12	210R1/50L	210R1/50R





Foot pressure measurement -Left: Daytime shoe Right: Luna

Forefoot Protective Cap					
Shoe Size (UK)	Shoe Size (EU)	Produc	ct Code		
Silve Size (UK)	Silve Size (EU)	Left	Right		
2.5 - 3.5	35 - 36	210R12/36L	210R12/36R		
4 - 5	37 - 38	210R12/38L	210R12/38R		
5.5 - 6.5	39 - 40	210R12/40L	210R12/40R		
7 - 8	41 - 42	210R12/42L	210R12/42R		
9 - 9.5	43 - 44	210R12/44L	210R12/44R		
10.5 - 11	45 - 46	210R12/46L	210R12/46R		
12 - 13	47 - 48	210R12/48L	210R12/48R		
13.5 - 14.5	49 - 50	210R12/50L	210R12/50R		





Shoe Size (UK)	Shoe Size (EU)	Product Code		
Silve Size (UK)	Silve Size (EU)	Left	Right	
2.5 - 3.5	35 - 36	210R13/36L	210R13/36R	
4 - 5	37 - 38	210R13/38L	210R13/38R	
5.5 - 6.5	39 - 40	210R13/40L	210R13/40R	
7 - 8	41 - 42	210R13/42L	210R13/42R	
9 - 9.5	43 - 44	210R13/44L	210R13/44R	
10.5 - 11	45 - 46	210R13/46L	210R13/46R	
12 - 13	47 - 48	210R13/48L	210R13/48R	
13.5 - 14.5	49 - 50	210R13/50L	210R13/50R	

ANKLE & FOOT

Mars S Therapeutic Shoes

Features and Benefits

- · Heel relief shoe
- New improved sole construction effects targeted shifting of the body centre of gravity onto the forefoot
- · Comfortable foot roll-over
- The shoe can be worn on the left or right side thanks to the hook and loop pad
- The pad can be customised for better medial support and relief of the heel
- Long, padded insole prevents plantar pressure peaks
- · Individually adjust to the ankle, back of the foot and forefoot area of the shaft
- Anti-bacterial shoe lining and cover fabric of the insole

Indications for Use

- Wounds in the heel area caused by diabetes mellitus
- Peripheral artery occlusions
- Other wounds in the heel area
- Post-surgical and traumas





Foot pressure measurement -Left: Daytime shoe Right: Mars S

Forefoot Cap					
Size	Shoe Size (UK)	Shoe Size (EU)	Product Code		
Small	2 - 3.5	34 - 36	213R10/S		
Medium	4 - 5.5	37 - 39	213R10/M		
Large	6.5 - 8	40 - 42	213R10/L		
X-Large	9 - 10.5	43 - 45	213R10/XL		
XX-Large	11 - 13	46 - 48	213R10/XXL		





Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Width (cm)	Product Code
Small	2 - 3.5	34 - 36	24	8.7	213R1/S
Medium	4 - 5.5	37 - 39	26	9.1	213R1/M
Large	6.5 - 8	40 - 42	28	10	213R1/L
X-Large	9 - 10.5	43 - 45	30	10.5	213R1/XL
XX-Large	11 - 13	46 - 48	32	10.9	213R1/XXL





Heel Protection Cap		Cushioning Insole			
Size	Shoe Size (UK)	Product Code	Size	Shoe Size (UK)	Product Code
Small	2 - 3.5	213R11/S	Small	2 - 3.5	213R12/S
Medium	4 - 5.5	213R11/M	Medium	4 - 5.5	213R12/M
Large	6.5 - 8	213R11/L	Large	6.5 - 8	213R12/L
X-Large	9 - 10.5	213R11/XL	X-Large	9 - 10.5	213R12/XL
XX-Large	11 - 13	213R11/XXL	XX-Large	11 - 13	213R12/XXL

164

ANKLE & FOOT



Features and Benefits

- Long term surgical shoe
- · Accommodates existing, large-volume dressings and bandages
- Keeps feet warm and dry (light splashes
- Encompasses the entire foot and ankle area with a softly padded socket
- Hook and loop fasteners at the wide flap prevent slipping within the shoe
- Removable padded insole prevents pressure peaks
- Wide opening
- · Anti-bacterial inner lining and cover fabric of the padded insole



Indications for Use

- Wounds in the foot and ankle area caused by diabetes mellitus
- Peripheral artery disease (PAD)
- Other indication independent wound formation in the foot and ankle area
- Postoperative and trauma

Ordering Information

Shoe Size	e Size Shoe Size Inner Shoe Width		Product Code		
(UK)	(EU)	Length (cm)	(cm)	Left	Right
3.5	36	24.3	9.3	214R1/36L	214R1/36R
4	37	24.9	9.4	214R1/37L	214R1/37R
5	38	25.6	9.6	214R1/38L	214R1/38R
5.5	39	26.3	9.7	214R1/39L	214R1/39R
6.5	40	26.9	9.9	214R1/40L	214R1/40R
7	41	27.6	10	214R1/41L	214R1/41R
8	42	28.3	10.2	214R1/42L	214R1/42R
9	43	29	10.3	214R1/43L	214R1/43R
9.5	44	29.6	10.5	214R1/44L	214R1/44R
10.5	45	30.3	10.6	214R1/45L	214R1/45R
11	46	30.9	10.8	214R1/46L	214R1/46R
12	47	31.5	10.9	214R1/47L	214R1/47R
13	48	32.2	11.1	214R1/48L	214R1/48R
13.5	49	32.7	11.3	214R1/49L	214R1/49R
14.5	50	33.3	11.4	214R1/50L	214R1/50R



Cushioning Insole				
Chan Cina (LIIV)	Chan Cina (FII)	Product Code		
Shoe Size (UK)	Shoe Size (EU)	Left	Right	
3.5	36	214R12/36L	214R12/36R	
4	37	214R12/37L	214R12/37R	
5	38	214R12/38L	214R12/38R	
5.5	39	214R12/39L	214R12/39R	
6.5	40	214R12/40L	214R12/40R	
7	41	214R12/41L	214R12/41R	
8	42	214R12/42L	214R12/42R	
9	43	214R12/43L	214R12/43R	
9.5	44	214R12/44L	214R12/44R	
10.5	45	214R12/45L	214R12/45R	
11	46	214R12/46L	214R12/46R	
12	47	214R12/47L	214R12/47R	
13	48	214R12/48L	214R12/48R	
13.5	49	214R12/49L	214R12/49R	
14.5	50	214R12/50L	214R12/50R	



Foot pressure measurement -Left: Daytime shoe Right: Neptun

To order contact order@steepergroup.com

Neptun B Therapeutic Shoes

Features and Benefits

- Long term surgical shoe
- Accommodates existing, large-volume dressings and bandages
- Keeps feet warm and dry (light splashes only)
- Encompasses the entire foot and ankle area with a softly padded socket
- Hook and loop fasteners at the wide flap prevent slipping within the shoe
- Removable padded insole prevents pressure peaks in the foot sole area
- Particularly wide outsole, wide opening options in the flap and heel area providing generous space for bandaged feet
- High shaft for sufficient support
- Integrated anti-slip nubs
- Stiffened outsole for extra relief of the injured foot
- Anti-bacterial shoe lining and cover fabric of the insole



Indications for Use

- Wounds in the foot and ankle area caused by diabetes mellitus
- Peripheral artery disease (PAD)
- Other indication independent wound formation in the foot and ankle area
- Postoperative and trauma
- Rehabilitation of lymphoedema, venous leg ulcer, elephantiasis, obesity

Shoe Size Inner Shoe Width		Width	Product Code		
(UK)	(EU)	Length (cm)	(cm)	Left	Right
2.5 - 3.5	35 - 36	23.9	9.9	215R1/36L	215R1/36R
4 - 5	37 - 38	25.2	10.2	215R1/38L	215R1/38R
5.5 - 6.5	39 - 40	26.5	10.5	215R1/40L	215R1/40R
7 - 8	41 - 42	27.9	10.8	215R1/42L	215R1/42R
9 - 9.5	43 - 44	29.2	11.1	215R1/44L	215R1/44R
10.5 - 11	45 - 46	30.5	11.4	215R1/46L	215R1/46R
12 - 13	47 - 48	31.9	11.7	215R1/48L	215R1/48R
13.5 - 14.5	49 - 50	33.1	12	215R1/50L	215R1/50R



Cushioning Insole					
Chan Cina (UIV)	Shoo Sino (FII)	Produc	t Code		
Shoe Size (UK)	Shoe Size (EU)	Left	Right		
2.5 - 3.5	35 - 36	212R11/36L	212R11/36R		
4 - 5	37 - 38	212R11/38L	212R11/38R		
5.5 - 6.5	39 - 40	212R11/40L	212R11/40R		
7 - 8	41 - 42	212R11/42L	212R11/42R		
9 - 9.5	43 - 44	212R11/44L	212R11/44R		
10.5 - 11	45 - 46	212R11/46L	212R11/46R		
12 - 13	47 - 48	212R11/48L	212R11/48R		
13.5 - 14.5	49 - 50	212R11/50L	212R11/50R		




Pluto S Therapeutic Shoes

Features and Benefits

- Features a rigid, slightly ascending outsole towards the front (negative heel)
- Pressure relief and pressure distribution during the gait cycle
- Hook and loop fasteners at the wide flap prevent slipping within the shoe
- Front flap prevents lateralisation of the forefoot; the foot feels secure
- If it is not required, the front flap can be cut off
- Cubical cushion sole for selective pressure relief
- Flap can be opened widely providing sufficient space for bandaged feet
- Integrated anti-slip nubs
- Anti-bacterial inner lining and cover fabric of the padded insole
- Weight limit: 120kg



- Wounds in the foot and ankle area caused by diabetes mellitus
- Peripheral artery disease (PAD)
- Other indication independent wound formation in the foot and ankle area
- Postoperative and trauma
- Gait unsteadiness (vertigo, visual impairment)



Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Width (cm)	Product Code
2.5 - 3.5	35 - 36	23.8	9	216R1/36
4 - 5	37 - 38	25.1	9.2	216R1/38
5.5 - 6.5	39 - 40	26.4	9.5	216R1/40
7 - 8	41 - 42	27.7	9.7	216R1/42
9 - 9.5	43 - 44	29	10	216R1/44
10.5 - 11	45 - 46	30.3	10.2	216R1/46
12 - 13	47 - 48	31.6	10.5	216R1/48
13.5 - 14.5	49 - 50	33	11.2	216R1/50





Foot pressure measurement -Left: Balance shoe Right: Pluto S

Forefoot Protective Cap			
Shoe Size (UK)	Shoe Size (EU)	Product Code	
2.5 - 3.5	35 - 36	216R12/36	
4 - 5	37 - 38	216R12/38	
5.5 - 6.5	39 - 40	216R12/40	
7 - 8	41 - 42	216R12/42	
9 - 9.5	43 - 44	216R12/44	
10.5 - 11	45 - 46	216R12/46	
12 - 13	47 - 48	216R12/48	
13.5 - 14.5	49 - 50	216R12/50	



Cubical Cushioning Insole			
Shoe Size (UK)	Shoe Size (EU)	Product Code	
2.5 - 3.5	35 - 36	216R11/36	
4 - 5	37 - 38	216R11/38	
5.5 - 6.5	39 - 40	216R11/40	
7 - 8	41 - 42	216R11/42	
9 - 9.5	43 - 44	216R11/44	
10.5 - 11	45 - 46	216R11/46	
12 - 13	47 - 48	216R11/48	
13.5 - 14.5	49 - 50	216R11/50	

Sirius Therapeutic Shoes




Features and Benefits

- Compensates for the imbalance between the foot fitted with a relief shoe and the healthy side, thus avoiding negative effects on the knee or hip joint and on the spine due to the pelvic obliquity
- Anti-bacterial shoe lining and cover fabric of the insole
- Integrated anti-slip nubs

Indications for Use

- Reduced mobility when walking and standing due to the difference in height when fitted with relief shoes
- Enables and maintains an appropriate walking function
- Compensates for pelvic obliquity that is worthy of treatment



Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Width (cm)	Product Code
2.5 - 3.5	35 - 36	24	8.9	211R1/36
4 - 5	37 - 38	25.3	9.2	211R1/38
5.5 - 6.5	39 - 40	26.7	9.5	211R1/40
7 - 8	41 - 42	28	9.8	211R1/42
9 - 9.5	43 - 44	29.3	10.1	211R1/44
10.5 - 11	45 - 46	30.6	10.4	211R1/46
12 - 13	47 - 48	32	10.7	211R1/48
13.5 - 14.5	49 - 50	33.2	11	211R1/50




Cubical Cushioning Insole			
Shoe Size (UK)	Shoe Size (EU)	Product Code	
2.5 - 3.5	35 - 36	216R11/36	
4 - 5	37 - 38	216R11/38	
5.5 - 6.5	39 - 40	216R11/40	
7 - 8	41 - 42	216R11/42	
9 - 9.5	43 - 44	216R11/44	
10.5 - 11	45 - 46	216R11/46	
12 - 13	47 - 48	216R11/48	
13.5 - 14.5	49 - 50	216R11/50	



Foot pressure measurement -Left: Sirius Right: Pluto S



172



ADULT INSOLES

PAEDIATRIC



Silicone Heel Cups

Features and Benefits

- Silicone Heel Cups absorb and dissipate shock forces that may cause pain or discomfort
- Designed for relief of heel pain and light rearfoot trauma, for use during sporting activities or simply to improve comfort
- Manufactured from medical grade silicone, each heel cup incorporates calcaneal support

Indications for Use

- Heel pain
- Calcaneal spur
- · Plantar fasciitis
- Metatarsalgia
- Achilles tendonitis



3/4 Length Silicone Insoles

Features and Benefits

- Silicone insoles absorb and dissipate shock forces that may cause pain or discomfort
- Designed for relief of heel pain and light rearfoot trauma, for use during sporting activities or simply to improve comfort
- Manufactured from medical grade silicone, each insole incorporates calcaneal support

Indications for Use

- Heel pain
- Calcaneal spur
- · Plantar fasciitis
- Metatarsalgia

Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Product Code
Small	11 - 1.5	29 - 33	YT-S04S
Medium	2 - 4.5	34 - 37	YT-S04M
Large	5 - 9.5	38 - 43	YT-S04L
X-Large	10 - 12	44 - 46	YT-S04XL

Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Product Code
Small	1 - 4	33 - 37	YT-S02S
Medium	4.5 - 8.5	38 - 42	YT-S02M
Large	9 - 14	43 - 48	YT-S02L



More Information

Wear and Care Instructions



www.steepergroup.com



Full Length Silicone Insoles with Metatarsal Support

Features and Benefits

- Silicone insoles absorb and dissipate shock forces that may cause pain or discomfort
- Designed for relief of heel pain and light rearfoot trauma, for use during sporting activities or simply to improve comfort
- Manufactured from medical grade silicone, each insole incorporates calcaneal support
- The full length insole also features metatarsal support

Indications for Use

- · Heel pain
- Calcaneal spur
- · Plantar fasciitis
- Metatarsalgia

Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Product Code
Small	4 - 5	37 - 38	YT-S01S
Medium	6 - 6.5	39 - 40	YT-S01M
Large	7 - 7.5	41 - 42	YT-S01L
X-Large	8.5 - 9.5	43 - 44	YT-S01XL





Slim Gel-Fit™ Bunion Sleeve

Features and Benefits

- The Slim Gel-Fit[™] Bunion Sleeve is an elastic sleeve with a gel pad incorporated inside two layers of a comfortable and conforming fabric
- Thin design that fits most footwear and can be worn under hoisery
- Sleeve can also be worn on the hand to protect the thumb
- Washable and reuseable
- The covered gel pad shields the painful and sensitive bony protrusion of the forefoot
- Thin and compact design makes it perfect for everyday use and fits with most dress, casual or athletic shoes
- · Latex free and hypo-allergenic

Indications for Use

- Helps to relieve shoe pressure, friction and shear forces on the bunion
- May also be worn on the hand to protect the thumb from discomfort during gripping activities

Ordering Information

Size	Product Code
Small	SI0048
Medium	SI0050

ADULT INSOLES

PAEDIATRIC

ADULT INSOLES

Pair

Diaped Duosoft Therapeutic

Features and Benefits

- Constructed with a multi-density material combination of plastazote and cellular urethane, these advanced insoles offer superb resilience, cushioning and shock absorption
- The insoles are designed to fit into both traditional footwear and 'additional depth' orthopaedic shoes
- The addition of an increased metatarsal support in the anatomical design is designed to facilitate off-loading of the metatarsal head area

Indications for Use

- For individuals with diabetes, arthritis, or other conditions that result in sensitive feet
- Metatarsalgia
- Pre-ulcerated sites
- Plantar fasciitis
- Foot deformities
- Morton's neuralgia
- Forefoot accommodation
- · Skin breakdown

Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Product Code
2 - 4.5	34 - 37	SB2335
5 - 6.5	36 - 38	SB2336
7 - 8.5	39 - 42	SB2337
9 - 10.5	43 - 45	SB2338
11 - 13	46 - 48	SB2339



Diaped Duosoft Plus

Features and Benefits

- Therapeutic insoles for insensate feet
 - The dual-layer design of Duosoft Plus consists of a resilient, shock absorbing urethane base and a PORON 94* top cover
- The slow rebound quality of the material cradles the foot and allows enhanced redistribution of body weight over the foot's plantar surface
- Soft, comfortable and allows the feet to remain cool, dry and odourless
- A Vygel heel accompanies the dual-layer formation to provide optimal shock attenuation at heel strike

Indications for Use

- For individuals with diabetes, arthritis, or other conditions that result in sensitive feet
- Metatarsalgia
- Pre-ulcerated sites
- Plantar fasciitis
- · Foot deformities
- Morton's neuralgia
- · Forefoot accommodation
- · Skin breakdown

Size	Shoe Size (UK)	Product Code
А	2 - 4	SB2402
В	5 - 7	SB2404
С	8 - 9	SB2406
D	10 - 12	SB2408
E	13 - 15	SB2410

ARM & SHOULDER

ADULT INSOLES

Tred-lite Ezi-Rize **Heel Raise**

Features and Benefits

- Adjustable heel raise insoles manufactured in multiple layers of EVA (5x3mm) with simple peel away layers to achieve the required raise height
- Depth ranges from 3mm to 15mm in 3mm increments

Indications for Use

- · Back pain
- Posture related symptoms
- Arch pain
- Leg length discrepancy
- Achilles tendonitis



Trio Heel Pain

Features and Benefits

- The Trio Heel Pain orthosis is contoured to support the navicular bone of the arch, directly reducing the cause of plantar fascia strain
- Higher loading forces under the heel can also lead to inflammation of the skin and supporting structures of the foot
- Features a large cushioning protection for the softer structure under the heel bones

Indications for Use

- · Plantar fasciitis
- Plantar heel pain
- Calf strain
- · Anterior medial knee pain
- Hip instability

Ordering Information

Size	Heel Width (cm)	Product Code
Men's	7	EZI-RIZE/MALE
Women's	5.5	EZI-RIZE/FEMALE

Size	Shoe Size (UK)	Density	Product Code
Small	5 - 6.5	Medium	THP-S
Medium	7 - 8.5	Medium	THP-M
Large	9 - 10.5	Medium	THP-L

HEAD & NECK

Trio Easy Fit

Features and Benefits

- The Trio Easy Fit is intrinsically a 3/4 length medially posted orthotic with an additional posting option for triplanar control of pronation featuring a lateral cutaway design
- · This model maintains the essential features for pronation control, including an optional biplanar rearfoot post
- · A medial forefoot post can be used to replace the first metatarsal cushion if required

Indications for Use

- Excessive pronation
- Forefoot pain
- Knee pain
- Heel valgus

Note: Paediatric sizing also available. See Page 247.



Trio Slimfit

Features and Benefits

- The recessed medial heel post within the Trio Slimfit design offers control over pronation without risking irritation of the skin within the tighter fitting fashion
- The forefoot contouring provides a stabilising effect upon the forefoot joints and reduces cramping pain

Indications for Use

- Excessive pronation
- Forefoot pain
- Knee pain
- · Heel valgus

Ordering Information

Size	Shoe Size (UK)	Product Code
X-Small	3 - 4.5	TEF-XS
Small	5 - 6.5	TEF-S
Medium	7 - 8.5	TEF-M
Large	9 - 10.5	TEF-L
X-Large	11 - 12.5	TEF-XL

Ordering Information

Size	Shoe Size (UK)	Product Code	
X-Small	3 - 4.5	TSF-XS	
Small	5 - 6.5	TSF-S	
Medium	7 - 8.5	TSF-M	
Large	9 - 10.5	TSF-L	

ADULT INSOLES

186

Trio

Features and Benefits

- The Trio insole is available in ³/₄ and full length and is designed to resist excessive valgus rotation
- Inverts the calcaneus within the heel cup thereby maintaining mid tarsal joint
- Supports the proximal medial arch of the foot by providing navicular stability
- Each pair is supplied with additional 3° medial biplanar rearfoot posts and 3° lateral forefoot posts

Indications for Use

- Plantar fasciitis
- Heel pain
- · Calf strain
- · Anterior medial knee pain
- Hip instability

Note: Paediatric sizing also available. See Page 248.

Ordering Information

3/4 Length				
Size	Shoe Size	Product Code		
Size	(UK)	Medium	Firm	Extra Firm
X-Small	3 - 4.5	TSLM-XS	TSLF-XS	TSLE-XS
Small	5 - 6.5	TSLM-S	TSLF-S	TSLE-S
Medium	7 - 8.5	TSLM-M	TSLF-M	TSLE-M
Large	9 - 10.5	TSLM-L	TSLF-L	TSLE-L
X-Large	11 - 12.5	TSLM-XL	TSLF-XL	TSLE-XL
XX-Large	13 - 14.5	TSLM-XXL	TSLF-XXL	TSLE-XXL

Full Length				
G:	Shoe Size	Product Code		
Size	(UK)	Medium	Firm	Extra Firm
X-Small	3 - 4.5	TFLM-XS	TFLF-XS	TFLE-XS
Small	5 - 6.5	TFLM-S	TFLF-S	TFLE-S
Medium	7 - 8.5	TFLM-M	TFLF-M	TFLE-M
Large	9 - 10.5	TFLM-L	TFLF-L	TFLE-L
X-Large	11 - 12.5	TFLM-XL	TFLF-XL	TFLE-XL
XX-Large	13 - 14.5	TFLM-XXL	TFLF-XXL	TFLE-XXL



Equiflex

Features and Benefits

- By elevating the heel the Equiflex will restore ankle joint extension reducing the demand for compensation
- Lateral forefoot is supported at the level of the metatarsal heads to facilitate this and accommodate any additional postings that may be required
- Supplied with 0°, 3° and 5° rear posts and 6mm elevators

Indications for Use

- High arches
- Excessive pronation
- Forefoot pain
- Heel valgus

Size	Shoo Sizo (LIV)	Product Code	
Size	Shoe Size (UK)	Medium	Firm
X-Small	3.5 - 4.5	EQM-XS	EQF-XS
Small	5 - 6	EQM-S	EQF-S
Medium	6.5 - 7.5	EQM-M	EQF-M
Large	8 - 9	EQM-L	EQF-L
X-Large	9.5 - 10.5	EQM-XL	EQF-XL
XX-Large	11 - 12	EQM-XXL	EQF-XXL



Prostep

Features and Benefits

- Prostep is designed to restore the natural position of the talus by supporting the inside edge of the heel and supports the spring ligament
- Resists the valgus rotation of the rearfoot, which is associated with pronation
- Anti-bacterial top cover and forefoot and rearfoot integral posts of 0°, 3° & 5°
- Available in full or 3/4 length

Indications for Use

- Low arches
- Excessive pronation
- Forefoot pain
- Knee pain
- Heel valgus

Note: Paediatric sizing also available. See Page 249.

Ordering Information

34 Length		
Shoe Size	Product Code	
(UK)	3/4 Medium	¾ Firm
3.5 - 4.5	PSSLM-XS	PSSLF-XS
5 - 6	PSSLM-S	PSSLF-S
6.5 - 7.5	PSSLM-M	PSSLF-M
8 - 9	PSSLM-L	PSSLF-L
9.5 - 10.5	PSSLM-XL	PSSLF-XL
11 - 12	PSSLM-XXL	PSSLF-XXL
	(UK) 3.5 - 4.5 5 - 6 6.5 - 7.5 8 - 9 9.5 - 10.5	3.5 - 4.5 PSSLM-XS 5 - 6 PSSLM-S 6.5 - 7.5 PSSLM-M 8 - 9 PSSLM-L 9.5 - 10.5 PSSLM-XL

Full Length			
Size	Shoe Size	Product Code	
Size	(UK)	Full Length Medium	Full Length Firm
X-Small	3.5 - 4.5	PSFLM-XS	PSFLF-XS
Small	5 - 6	PSFLM-S	PSFLF-S
Medium	6.5 - 7.5	PSFLM-M	PSFLF-M
Large	8 - 9	PSFLM-L	PSFLF-L
X-Large	9.5 - 10.5	PSFLM-XL	PSFLF-XL
XX-Large	11 - 12	PSFLM-XXL	PSFLF-XXL



Control Tech SoftFull Length

Features and Benefits

- Available in low, moderate and high arch angles (4°, 6° & 8°)
- Plantar fascial grooves reduce arch irritation as well as provide comfort and superior shock absorption

Indications for Use

- Metatarsalgia
- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

Ordering Information

Shoe Size		Product Code	
(UK)	4° Low Arch	6° Moderate Arch	8° High arch
10 - 13	SC3086	-	-
1 - 2	SC3088	-	-
2 - 3	SC3090	SC3040	-
3 - 4	SC3091	SC3041	SC3031
5 - 6	SC3092	SC3042	SC3032
7 - 8	SC3093	SC3043	SC3033
9 - 10	SC3094	SC3044	SC3034
11 - 12	SC3095	SC3045	-
13 - 15	SC3096	SC3046	-



Control Tech Soft 3/4 Length

Features and Benefits

• The rearfoot wedging reduces the need for excessive arch height while maintaining support and comfort, providing superior shock absorption

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- · Back pain
- Postural related symptoms

Note: Paediatric sizing also available. See Page 246.



Control Tech Flex 3/4

Features and Benefits

- ¾ length insole eliminates pronation moments around the subtalar joint created by the lateral heel cup
- Supplied with optional self adhesive
- Polyprolylene integral posting over the subtalar joint axis

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

Ordering Information

Shoe Size (UK)	Product Code	Shoe Size (UK)	Product Code
2 - 3	SC3055	7 - 8	SC3058
3 - 4	SC3056	9 - 10	SC3059
5 - 6	SC3057	11 - 12	SC3060

Ordering Information

Shoe Size		Product Code		
(UK)	4° Low Arch	6° Moderate Arch	8° High arch	
3 - 4	SC3021	SC3011	SC3001	
5 - 6	SC3022	SC3012	SC3002	
7 - 8	SC3023	SC3013	SC3003	
9 - 10	SC3024	SC3014	SC3004	
11 - 12	SC3025	SC3015	SC3005	

ADULT INSOLES

Interpod Diabetic

Features and Benefits

- Inverted rearfoot improves biomechanical function and enhances the windlass effect
- The 6° moderate arch height provides adequate functionality for the majority
- PORON 96 top cover sets a new level of advanced underfoot technology
- The slow rebound cradles the foot and allows enhanced distribution of body weight over the foot's plantar surface
- · Cradle functionality makes this the ultimate stock diabetic orthosis

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- · Tendon pain
- Back pain
- Postural related symptoms

Ordering Information

Shoe Size (UK)	Product Code	
2 - 3	SC3130	
3 - 4	SC3131	
5 - 6	SC3132	
7 - 8	SC3133	

Shoe Size (UK)	Product Code
9 - 10	SC3134
11 - 12	SC3135
13	SC3136



Darco Universal Toe Alignment Splint

Features and Benefits

- Universal toe alignment splint
- Constructed of soft cotton backed elastic material that will fit into standard footwear

Indications for Use

- Designed to maintain position of soft tissues until they heal following surgery for bunions or hammer toes
- Reduces incidence of recurrence of the deformity following surgical correction

Size	Product Code
Universal	NV0338

Slimflex Simple Kit ³/₄ Length

Features and Benefits

- Brush Nylon cover allows for comfort and moisture absorption
- High density EVA (A65) construction to cover most pathologies
- Integrated first metatarsal cut-out to facilitate the windlass mechanism
- Metatarsal raise and sustanticulum tali support
- Flat base for adhering wedging
- · Deep heel cup for greater control
- Fully heat mouldable and grindable
- 6mm wedging comes as standard

Indications for Use

- Generalised foot pain
- · Heel and arch pain
- Anterior knee pain
- Plantar Fasciitis



Pair

Slimflex Green InsolesFull Length

Features and Benefits

- Ideal initial introduction to orthotic therapy, perfect for first line treatment
- Perforated under layer to increase breathability and reduce the risk of hyperhydration
- High arch profile
- Hard wearing top cover
- Exceptional value for money
- Ideal to trial different postings at budget friendly prices
- Moulded EVA with a deep heel cup which allows greater foot control and correction when used with Podotech wedges
- Heat moulding possibilities which allows for off-loading implementation
- Antimicrobial cover to reduce odour and microbial colonisation

Indications for Use

- · Generalised foot pain
- Heel and arch pain
- · Anterior knee pain
- Plantar Fasciitis

Ordering Information

Shoe Size (UK)	Product Code
3	SB2223
4	SB2224
5	SB2225
6	SB2226
7	SB2227
8	SB2228
9	SB2229
10	SB2230
11	SB2231
12	SB2232
13	SB2233
14	SB2234
15	SB2235

Ordering Information

Size	Shoe Size (UK)	Product Code
X-Small	3 - 4.5	SB4017
Small	5 - 6.5	SB4018
Medium	7 - 8.5	SB4019
Large	9 - 10.5	SB4020
X-Large	11 - 12.5	SB4021
XX-Large	13 - 14.5	SB4022

ADULT INSOLES



Vectorthotic Full Length Insoles

Features and Benefits

- This modular system brings all the advantages of custom foot orthoses without the need for casting or waiting
- Each pack contains a polypropylene composite shell, top covers and a range of medial clip-in rearfoot (2°, 4°, & 6°) and medial and lateral pre-glued forefoot (2° & 4°) posts
- The shell can be customised to the patient's requirements by clipping in the appropriate post to invert the rear foot, support the forefoot with a varus, valgus or no post as deemed suitable
- If desired, apply the moulded top cover incorporating a 2mm neoprene comfort layer, which can then be cut to the required length
- Other additions such as heel lifts or extension postings can easily be added
- Being an injection moulded span device, its strength lies in the integrity of the design and durability of the shell

Indications for Use

- Arch ache or pain
- Exercise induced kneecap pains
- Heel and ankle symptoms
- Hip and back pain on exercise
- Musculoskeletal pains linked to poor foot position and function

Ordering Information

	Product Code	
Shoe Size (UK)	Standard Top Cover (Lower arch profile)	Extra Top Cover (Higher arch profile and metatarsal support)
3 - 4	VEC05	VEC05X
5 - 6	VEC04	VEC04X
6 - 7	VEC03W	VEC03WX
7 - 8	VEC03	VEC03X
8 - 9	VEC02W	VEC02WX
9 - 10	VEC02	VEC02X
11 - 12	VEC01	VEC01X





Features and Benefits

- · Features heel cradle and midfoot contoured support for maximum comfort compared to a traditional 'arch support'
- Insole is full-length and can be quickly and easily trimmed to fit most comfortable shoes
- The brushed nylon surface allows easy passage of the foot into the shoe

Indications for Use

- Arch aches or pains from foot strain
- Pain under forefoot and toe area (metatarsalgia)
- Heel and ankle pains, including mild Plantar Fasciitis (policeman's heel) and Achilles symptoms (Achilles Tendonitis)

Size	Shoe Size (UK)	Product Code
Small	3 - 5	XLINE04
Medium	6 - 8	XLINE03
Large	9 - 11	XLINE02
X-Large	12 - 14	XLINE01
XX-Large	15 - 16	XLINE01X

X-Line Reinforced 3/4 Features and Benefits



- The deep heel cup offers reassuring and effective control of the rearfoot
- There is an extra anti-pronatory feature within the heel cup so the rearfoot is actively supported
- The midfoot saddle provides essential support and improves stability throughout gait
- Unique medial arch under-bevel
- Designed to be your go-to device for many conditions particularly where footwear accommodation is an issue



- · Arch and foot strain
- Ankle pain
- Plantar heel pain
- Posterior Tibial Tendon Dysfunction
- · Situations where full-length devices aren't a good option



X-Line 70 Insoles

Features and Benefits

- · Features heel cradle and midfoot contoured support for maximum comfort compared to a traditional 'arch
- Insole is full-length and can be quickly and easily trimmed to fit most comfort shoes
- The brushed nylon surface allows easy passage of the foot into the shoe
- Made from high-density foam giving a firm, highly durable foot support

Indications for Use

- Arch aches or pains from foot strain
- Pain under forefoot and toe area (metatarsalgias)
- Heel and ankle pains, including mild Plantar Fasciitis (policeman's heel) and Achilles symptoms (Achilles Tendonitis)

Ordering Information

Size	Shoe Size (UK)	Product Code
Small	3 - 5	XLINE3/404
Medium	6 - 8	XLINE3/403
Large	9 - 11	XLINE3/402
X-Large	12 - 14	XLINE3/401

Ordering Information

Size	Shoe Size (UK)	Product Code	
Small	3 - 5	XLINE0470	
Medium	6 - 8	XLINE0370	
Large	9 - 11	XLINE0270	
X-Large	12 - 14	XLINE0170	
XX-Large	15 - 16	XLINE01X70	

198


X-Line Pressure Perfect Insoles

Features and Benefits

- · Features special cushioning top cover material
- Heel cradle and midfoot contour supports and protects the foot by reducing heel and forefoot pressures
- Insole is full-length and can be quickly and easily trimmed to fit most comfort
- · Cushioning and support protects vulnerable feet from abnormal heel and forefoot pressure
- · Helps keep you moving in comfort for longer

Indications for Use

- · Ideal for diabetic feet
- · Pain in feet due to arthritis
- burning sensations (metatarsalgias)
- Active older feet
- Flat feet due to weight issues



- · High impact activities such as squash, badminton, netball and basketball
- Roadrunners



Size	Shoe Size (UK)	Product Code
Small	3 - 5	PRE04
Medium	6 - 8	PRE03
Large	9 - 11	PRE02
X-Large	12 - 14	PRE01
XX-Large	15 - 16	PRE01X



X-Line RIF (ReInForcement) Insoles

Features and Benefits

- Features plastic wrap ReInForcement enhancing durability and strength
- Heel cradle and midfoot contour supports the foot throughout every step providing maximum comfort compared to a traditional 'arch support'
- Insole is full-length and can be quickly and easily trimmed to fit most comfort
- Strong Cambrelle® top cover for long life and enhanced insole rigidity
- Ideal for those regularly active on their
- Fits in any shoe with a removable insock/insole, most lace up shoes, trainers and walking shoes

Indications for Use

- Arch and foot strain
- · Pain under forefoot and toe area (metatarsalgias)
- Ankle pains, including mild Plantar Fasciitis (policeman's heel) and Achilles symptoms
- Calf strains and other recurrent heel pain after running

Ordering Information

Size	Shoe Size (UK)	Product Code	
Small	3 - 5	XLINERE04	
Medium	6 - 8	XLINERE03	
Large	9 - 11	XLINERE02	
X-Large	12 - 14	XLINERE01	
XX-Large	15 - 16	XLINERE01X	

X-Line Achilles Tendonitis Insoles

Features and Benefits

- Reinforced for improved stability
- Boasts a hindfoot lift which raises the heel reducing the strain on the Achilles during the driving phase of gait
- Features a mid-foot saddle and a 'profiled arch stabilising contour', designed to reduce the speed of distal calcaneal plantarflexion of the foot during weight acceptance and prevent excess loss of navicular height

Indications for Use

• Resolves over strain and pain of the Achilles (calcaneal) tendon



X-Line TPD Insoles Firm Density

Features and Benefits

- For Tibialis Posterior Dysfunction
- Features high heel cups for heel stabilisation
- High arch support placed specifically to offload the tibialis posterior tendon
- Midfoot cradle and metatarsal support
- High plastic wrap reinforcement to strengthen insole
- Camberelle® top cover for strength and durability

Indications for Use

- Severe arch ache or pain
- Severe medial ankle and shin pain
- Any symptoms caused by large arch
- Perfect for immediate relief of Tibialis Posterior Dysfunction

Ordering Information

Size	Shoe Size (UK)	Product Code
Small	3 - 5	XLINEAT04
Medium	6 - 8	XLINEAT03
Large	9 - 11	XLINEATO2
X-Large	12 - 14	XLINEAT01
XX-Large	15 - 16	XLINEAT01X

Size	Shoe Size	Product Code	
Size	(UK)	Medium Density	Firm Density
Small	3 - 5	XLINETPD50NW04	XLINETPD04
Medium	6 - 8	XLINETPD50NW03	XLINETPD03
Large	9 - 11	XLINETPD50NW02	XLINETPD02
X-Large	12 - 14	XLINETPD50NW01	XLINETPD01
XX-Large	15 - 16	XLINETPD50NW01X	XLINETPD01X

202

X-Line DJD Insoles

Features and Benefits

- Hallux Rigidus insoles and Hallux Limitis Insoles designed to alleviate big toe joint pain or ache during and after activity
- The first insole designed for arthritis big toe. Widely used in the NHS for the treatment of Hallux Rigidus and Hallux Limitus
- Features supports and contours that reduce the amount of bend in the big toe when walking
- Forefoot rocker precisely located to reduce big toe stress
- Heel and midfoot cradled to guarantee comfort
 - The Cambrelle® top cover gives durability and an easy surface to clean

Indications for Use

- Osteoarthritis or Degenerative Joint Disease (DJD) of the big toe joint
- · Big toe pain during movement
- Deep pain from bunions
- Big toe sesamoiditis
- Golf or hiking induced big toe joint pain
- Excellent for use after big toe joint surgery including fusion



X-Line Plantar Fasciitis Insoles

Features and Benefits

- Features a heel lift which reduces the strain on the Plantar Fascia
- Arch support to take the load from the
- Heel indentation and shock dot positioned to take pressure off the most sensitive area of the heel (medial process of plantar calcaneal tubercle)
- The insole also comes with an intrinsic lateral skive to the rear foot and lateral forefoot which reduces 'sudden' eversion moments
- The 1st metatarsal recess and lesser metatarsal transverse plane support the insole to improve 1st metatarsophalangeal joint function and create a stable forefoot during heel lift

Indications for Use

- To relieve heel pain caused by Plantar Fasciitis
- · Relief of heel spur pain

Ordering Information

Size	Shoe Size (UK)	Product Code
Small	3 - 5	XLINEDJD04
Medium	6 - 8	XLINEDJD03
Large	9 - 11	XLINEDJD02
X-Large	12 - 14	XLINEDJD01
XX-Large	15 - 16	XLINEDJD01X

Ordering Information

Size	Shoe Size (UK)	Product Code
Small	3 - 5	XLINEPF04
Medium	6 - 8	XLINEPF03
Large	9 - 11	XLINEPF02
X-Large	12 - 14	XLINEPF01
XX-Large	15 - 16	XLINEPF01X

To order contact order@steepergroup.com



X-Line Additional Wedges

Features and Benefits

- Self adhesive
- Rear foot lateral wedges available in two heights, 4° and 8°
- Alleviate forefoot wedges available in 4°
- Alleviate heel raise 6mm

Indications for Use

• Increasing biomechanical correction of the pronated or supinated foot



Note: Wedges are designed for use with the corresponding size of insole, however, some adjustment may be required.

Rear Foot Wedges (Supplied in 5 pairs)			
Ci	Product Code		
Size	4 °	8°	
Small/Medium	ALW01	ALW02	
Large/X-Large	ALW03	ALW04	





Alleviate Forefoot Wedges (Supplied in 5 pairs)		Self-Adhesive 6mm Heel Raise (Supplied individually)	
Size	Product Code	Size	Product Code
Small/Medium	ALVFW02	Small/Medium	ALVHR02
Large/X-Large/XX-Large	AI VFW01	Large/X-Large/XX-Large	ALVHR01





PAEDIATRIC

HEAD & NECK

Starlight Standard Features and Benefits



- Impact absorbing foam construction
- · Cotton jersey lining
- · Soft to the skin
- · Wicks away moisture
- · Good moisture circulation
- Outer layer made from robust cotton for Iongevity
- Lightweight approximately 100g
- Machine washable
- Available in head circumferences 45cm to 62cm in 1cm increments

Indications for Use

· Simple yet effective protection for people at risk of head injury

Styles Available



Note: This product is non-returnable for hygiene reasons.

*Starlight Standard in Brown is not kept in stock, delivery times may vary.

Ordering Information

Colour	Product Code
Brown	5180-
Red	5180-1-
Blue	5180-2-

= Circumference of head (45cm to 62cm in 1cm increments)

Example: For a red Starlight Standard Helmet measuring 50cm - 5180-1-50





Starlight Quilted Soft


Features and Benefits

- The Starlight Quilted Soft is upholstered with highly efficient high-tech foam
- The foam has excellent shock absorbing properties
- The forces applied by a fall are redistributed over a larger area
- Adapts better to the shape of the head offering improved comfort
- Available in head circumferences 45cm to 62cm in 1cm increments
- Excellent moisture wicking properties

Indications for Use

 Advanced protection for people at risk of head injury

Styles Available



Note: This product is non-returnable for hygiene reasons.

Ordering Information

Colour	Product Code
Oxford Rose	5180SG-1-
Oxford Light Blue	5180SG-2-

= Circumference of head (45cm to 62cm in 1cm increments)

Example: For an Oxford Rose Starlight Quilted Soft Helmet measuring 50cm -5180SG-1-50



Starlight Secure Fabric



Features and Benefits

- Impact absorbing soft/firm foam combination construction
- · High quality cotton
- · Soft to the skin
- Lightweight approximately 100g 130g (size dependent)
- Easy to wash

Indications for Use

 Advanced protection for people at risk of head injury

Note: Starlight Secure model shown with optional forehead protector, neck cushion and chin guard.

Note: This product is non-returnable for hygiene reasons.

Starlight Secure Leather

Features and Benefits

- Impact absorbing soft/firm foam combination construction
- · Leather lined
- Soft to the skin
- Lightweight approximately 150g 240g (size dependent)
- Heat and moisture regulating
- Breathable design
- · Easy to clean

Indications for Use

 Advanced protection for people at risk of head injury

Note: This product is non-returnable for hygiene reasons.

Ordering Information

Size	Head Circ. (cm)	Product Code	Accessories	Product Code
1	47 - 49	C401101B-1	Chin Guard	C4012B-1
2	50 - 52	C401102B-1	Forehead Protector	C4013B-1
3	53 - 55	C401103B-1	Neck Cushion	C4014B-1
4	56 - 59	C401104B-1		
5	60 - 62	C401105B-1		

Is grey not your colour?

Made to Measure service is available for this product!

Other colours are available when ordering made to measure.





Ordering Information

Size	Head Circ. (cm)	Product Code	Accessories	Product Code
1	47 - 49	C401101	Chin Guard	C4012
2	50 - 52	C401102	Forehead Protector	C4013
3	53 - 55	C401103	Neck Cushion	C4014
4	56 - 59	C401104		
5	60 - 62	C401105		

Is beige not your colour?

Made to Measure service is available for this product!

Other colours are available when ordering made to measure.





ADULT INSOLES

PAEDIATRIC

Universal Wrist Brace



Features and Benefits

- Breathable felt and neoprene mix material
- Aluminium palmer splint and stays can be shaped to fit the wrist with additional hook and loop strap for enhanced fit
- Designed for easy on and off application
- Circumferential seam for improved fit around waist

Indications for Use

- Arthritis
- Sprain
- Swelling
- Carpal Tunnel

Note: For adult version, see Page 62.



Universal Palm Protector

Features and Benefits

- Fleece-lined, lightweight palm cover
- Designed to allow full finger movement whilst reducing finger contractures and skin damage/breakdown in the palmar
- Easily laundered
- Easily applied using hook and loop fastening straps

Indications for Use

- Post surgery/injury positioning
- Hypertonic contractures
- · Rheumatoid arthritis

Ordering Information

Description	Size	Product Code
Left	Paediatric	K1SWW-LT
Right	Paediatric	K1SWW-RT

Ordering Information

Description	Size	Product Code
Right	Universal	PPR
Left	Universal	PPL

Comfy Hand-Wrist Finger Orthosis

Features and Benefits

- Terrycloth
- Supplied with one cover and finger separators
- Ambidextrous
- Side wings prevent ulnar and radial deviation
- Arm cuffs easily adjust to width of arm
- Excellent resting splint to prevent trauma to joints and increases ROM
- Multi-layered for comfort
- Adjustable reinforced spine allows for easy adjustments without tools
- Easy to clean
- Keeps thumb in functional position without stretching thenar eminence
- Washable and absorbent covers help prevent skin breakdown
- Soft straps accommodate fluctuations in oedema
- Arm, wrist, fingers independently adjustable
- Resists strong flexion synergy

Indications for Use

 This popular splint provides support and positioning for weak or deformed hands

Note: For larger sizing see Page 67.

Ordering Information

Size	Wrist to Top of Middle Finger (cm)	Total Length (cm)	Product Code
Paediatric Small	6.3 - 8.9	17.8	PH-101-S
Paediatric Medium	7.6 - 11.4	20.3	PH-101-M
Paediatric Large	11.4 - 15.2	24.1	PH-101-L

Comfy Deviation HandOrthosis



Features and Benefits

- Terrycloth
- Comes with one liner
- Ambidextrous
- Progressively change wrist and finger positions and amount of wrist deviation by adjusting the internal frame and locking the frame's swivel with included Allen key
- Correct or control wrist deviation motion even when changes in the position of the wrist and hand are made
- Unique below the wrist swivel ability that allows for lateral and medial movements
- Accommodates for ulnar and radial deviation

Indications for Use

- Arthritic or neuromuscular deformities
- Flaccid, weak or spastic hand conditions

Styles Available



Ordering Information

Size	Colour	Wrist to Top of Middle Finger (cm)	Total Length (cm)	Product Code
Small	Purple	6.3 - 8.9	16.5	PDH-101-S
Medium	Turquoise	7.6 - 11.4	17.8	PDH-101-M
Large	Emerald Green	11.4 - 15.2	22.9	PDH-101-L

216



Comfy Hand ThumbOrthosis

Features and Benefits

- Available in Terrycloth and Headliner
- · Supplied with one liner
- Ambidextrous
- Paediatric sizing
- Adjustable thumb support
- Keeps thumb in functional position without stretching thenar eminence

Indications for Use

• Excellent for patients with highly adducted thumbs



Ordering Information

Size	Wrist to Top of Middle Finger (cm)	Total Length (cm)	Product Code
Paediatric Small	6.3 - 8.9	17.8	PHT-101-S
Paediatric Medium	7.6 - 11.4	20.3	PHT-101-M
Paediatric Large	11.4 - 15.2	24.1	PHT-101-L







Jobskin® Arm Gaiters Standard Lining

Features and Benefits

- An external orthosis with structural inserts to immobilise the affected joint and hold the limb in controlled extension, or the best possible position for rehabilitation
- Provides a prolonged stretch, immobilisation and contracture management to help maintain the gains achieved by manipulation, corrective surgery or reconstructive procedures
- Can be used by infants through to adulthood
- Designed to meet the individual's functional needs rather than to treat an individual condition or pathology
- Ready-to-wear

Indications for Use

- Certain skin conditions e.g. Eczema
- Vascular/Circulation problems
- Oedema
- Weight fluctuations
- Epilepsy
- · Contracture management

Note: This product is non-returnable. A pattern must be stipulated when ordering. For available patterns, see page 221.

Ordering Information

Size	Top Circ. (cm)	Bottom Circ. (cm)	Length (cm)	Product Code
X-Small	18 - 25	13 - 20	20	AG090220
X-Siliali	10 - 25	13 - 20	25	AG090225
			20	AG090320
Small	24 - 31	15 - 22	25	AG090325
			30	AG090330
Ma alia	71 70	18 - 25	30	AG090430
Medium	31 - 38		35	AG090435
			40	AG090540
Large	39 - 46	21 - 28	45	AG090545
			50	AG090550
V Lawas	46 57	24 - 31	45	AG090645
X-Large	46 - 53		50	AG090650

Jobskin® Arm Gaiters Deluxe Lining

Features and Benefits

- An external orthosis with structural inserts to immobilise the affected joint and hold the limb in controlled extension, or the best possible position for rehabilitation
- Provides a prolonged stretch, immobilisation and contracture management to help maintain the gains achieved by manipulation, corrective surgery or reconstructive procedures
- Can be used by infants through to adulthood
- Designed to meet the individual's functional needs rather than to treat an individual condition or pathology
- Ready-to-wear

Indications for Use

- Certain skin conditions e.g. Eczema
- Vascular/Circulation problems
- Oedema

w

- Weight fluctuations
- Epilepsy
- Contracture management

Note: This product is non-returnable. A pattern must be stipulated when ordering. For available patterns, see page 221.

Ordering Information

Size	Top Circ. (cm)	Bottom Circ. (cm)	Length (cm)	Product Code
X-Small	18 - 25	13 - 20	20	AG150220
V-2IIIali	10 - 25	13 - 20	25	AG150225
			20	AG150320
Small	24 - 31	15 - 22	25	AG150325
			30	AG150330
Madium	71 70	10 05	30	AG150430
Medium	31 - 38	18 - 25	35	AG150435
			40	AG150540
Large	39 - 46	21 - 28	45	AG150545
		50	AG150550	
V Lavas	46 E7	0.4 71	45	AG150645
X-Large	46 - 53	24 - 31	50	AG150650

Bumblebee Camouflage Denim Dinosaurs Hearts and Flowers

Styles Available

Mint Butterflies

Note: Gaiter patterns are accurate at the time of printing and are subject to change.

Black

PAEDIATRIC

222



Paediatric Hip Abduction Splint - Standard

Features and Benefits

- Maintains concentric reduction of hips in patients with hip dysplasia
- Replaces spica casts and harnesses for certain developmental dysplasia patients
- One-Piece, non-adjustable polyproplene orthosis holds infant in 90° hip flexion and 60° abduction
- Open design facilitates nappy changes
- Closed-cell foam lining does not retain moisture

Indications for Use

• Indicated for patients with hip dysplasia



Size	Average Age (mnths)	Weight Guide (lbs)	Product Code
X-Small	<3	6 - 12	B/3899
Small	3 - 6	13 - 16	B/3900
Medium	6 - 12	17 - 21	B/3901
Large	12	22 - 30	B/3902
X-Large	18	30+	B/3903



Paediatric Hip Abduction Splint - Adjustable

Features and Benefits

- Lightweight, paediatric orthosis offers advantage of abduction adjustment over the standard paediatric hip splint model described opposite
- Sturdy polypropylene frame set in 90° hip flexion - abduction easily adjustable 40° through 180° by aluminium slotted screw assembly
- Fits infants from 3 months (6-12lbs) to 18 months (over 30lbs)
- More total contact than Pavlik harness
- Easy application, removal facilitates nappy changes
- Closed-cell foam lining does not retain moisture

Indications for Use

- Indicated for patients with hip dysplasia
- X-Large size effective for certain cerebral palsy patients

Size	Average Age (mnths)	Weight Guide (lbs)	Product Code
X-Small	<3	6 - 12	B/3906
Small	3 - 6	13 - 16	B/3907
Medium	6 - 12	17 - 21	B/3908
Large	12	22 - 30	B/3909
X-Large	18	30+	B/3910

 $\stackrel{\wedge}{\mathbb{H}}$

Mini TLC™ Hip **Abduction** Brace

Features and Benefits

- Prefabricated
- Material: MPE
- Pelvic band adjustable for hip development and pelvic circumference with quick release buckle closure
- Thigh cuff bivalve design adjustable to thigh proximal and distal measurements
- Breathable, washable liner absorbs perspiration and reduces heat build up
- · Universal hip joint, TLC hip joint with infinite abduction/adduction and flexion/ extension settings
- Available in 6 sizes

Indications for Use

- Arthoplasty
- · Hip dislocation or potential high risk dislocation
- Hip dysplasia anterior or posterior
- Hip management/immobilisation
- Hip revision
- Hip surgery
- Postoperative
- Preoperative
- · Stabilise, align and reinforce hip

Ordering Information

Size	Waist Circ. (cm)	Proximal Thigh Circ. (cm)	Distal Thigh Circ. (cm)	Product Code
XX-Small	40.6 - 45.7	29.2 - 34.3	20.3 - 25.4	TLCMINI2S
X-Small	45.7 - 50.8	31.8 - 36.8	22.9 - 27.9	TLCMINIXS
Small	50.8 - 55.9	34.3 - 39.4	25.4 - 30.5	TLCMINISM
Medium	55.9 - 61	36.8 - 41.9	27.9 - 33	TLCMINIMD
Large	61 - 66	39.4 - 44.4	30.5 - 35.6	TLCMINILG
X-Large	66 - 71.1	41.9 - 47	33 - 38.1	TLCMINIXL





Jobskin® Leg Gaiters Standard Lining

Features and Benefits

- External orthosis with structural inserts to immobilise the affected joint and hold the limb in controlled extension, or the best possible position for rehabilitation
- Providing a prolonged stretch, immobilisation and contracture management to help maintain the gains achieved by manipulation, corrective surgery or reconstructive procedures
- For use on one or more limbs and can be used by infants through to adulthood
- Designed to meet the individual's functional needs rather than to treat an individual condition or pathology
- · Ready-to-wear



- Certain skin conditions e.g. Eczema
- Vascular/Circulation problems
- Oedema
- Weight fluctuations
- Epilepsy
- Contracture management of the lower limb

Note: This product is non-returnable. A pattern must be stipulated when ordering. For available patterns, see page 221.

Ordering Information

Size	Top Circ. (cm)	Bottom Circ. (cm)	Length (cm)	Product Code
X-Small	18 - 25	13 - 20	20	LG090220
X-Small	18 - 25	13 - 20	25	LG090225
			20	LG090320
Small	24 - 31	15 - 22	25	LG090325
			30	LG090330
Medium	31 - 38	18 - 25	30	LG090430
Medium	31 - 38		35	LG090435
			40	LG090540
Large	39 - 46	21 - 28	45	LG090545
			50	LG090550
V I	4C F7	24 - 31	45	LG090645
X-Large	46 - 53		50	LG090650

PAEDIATRIC

Jobskin® Leg Gaiters Deluxe Lining

Features and Benefits

- · External orthosis with structural inserts to immobilise the affected joint and hold the limb in controlled extension, or the best possible position for rehabilitation
- · Providing a prolonged stretch, immobilisation and contracture management to help maintain the gains achieved by manipulation, corrective surgery or reconstructive procedures
- For use on one or more limbs and can be used by infants through to adulthood
- Designed to meet the individual's functional needs rather than to treat an individual condition or pathology
- · Ready-to-wear

Indications for Use

- Certain skin conditions e.g. Eczema
- Vascular/Circulation problems
- Oedema
- Weight fluctuations
- Epilepsy
- · Contracture management of the lower limb

Note: This product is non-returnable. A pattern must be stipulated when ordering. For available patterns, see page 221.



Ordering Information

Size	Top Circ. (cm)	Bottom Circ. (cm)	Length (cm)	Product Code
X-Small	18 - 25	13 - 20	20	LG150220
A-Silidii	10 - 25	13 - 20	25	LG150225
			20	LG150320
Small	24 - 31	15 - 22	25	LG150325
			30	LG150330
Medium	31 - 38		30	LG150430
Medium	31 - 38	18 - 25	35	LG150435
			40	LG150540
Large	39 - 46	21 - 28	45	LG150545
			50	LG150550
V Lawara	46 E7	24 71	45	LG150645
X-Large	46 - 53	24 - 31	50	LG150650



Made To Measure Leg Gaiters

Features and Benefits

- Choice of padding for prolonged use and a more comfortable fit
- Choice of fabrics
- Leg immobilisation

Indications for Use


· Contracture management of the lower

Note: A pattern must be stipulated when ordering. For available patterns, see page 221.

Size	Product Code
Leg Gaiter	Please contact Customer Services on
Leg Flexion Gaiter	order@steepergroup.com to request a
Lining as Required	measurement form.

Push Aequi Junior Ankle Brace

Features and Benefits



- Stabilises the ankle joint during inversion as well as eversion and neutral and flexed positions
- The freedom of dorsal and plantar flexion is sufficiently maintained to allow normal walking movement
- Optimal functioning of the brace is obtained when it is worn with a sturdy lace-up shoe



- Acute ankle ligament lesions (both inversion and eversion trauma
- Stabilisation of conservatively or operatively treated ankle fractures
- (Residual) instability
- lesions

Note: For larger sizing see Page 99.

- Secondary prevention of ankle ligament



Size	Product Code
Junior	PUS100B



UFO™

Features and Benefits

Introducing an improved prefabricated night splint for the effective treatment of plantar fasciitis and achilles tendonitis. The UFO™ comfortably supports the affected foot in a controlled amount of dorsiflexion, applying a Low Load Prolonged Stretch (LLPS) to the soft tissues of the plantar fascia and achilles tendon. LLPS promotes gentle stretching, healing and remodelling of these tissues, effectively reducing the stress and pain caused by inflammation and tissue shortening. Such support can help break the cycle of injury, pain and re-injury due to otherwise unlimited contraction and deformity of the soft tissues at night.

- The simple variable dorsiflexion design allows the patient to quickly and easily adjust the amount of applied stretch to provide the most effective pain relief
- Adjustable bilateral dorsifelxion assist straps
- CoolFoam[™] replaceable liner wicks moisture away from the skin for enhanced hygiene

Indications for Use

- Plantar fasciitis
- · Achilles tendonitis
- Contracture management of the foot and ankle

Size	Footplate Length (cm)	Brace Height (cm)	Product Code
Small	14	16.5	B/3614
Medium	17	25	B/3615
Large	19	26.5	B/3616



Universal PRAFOPaediatric

Features and Benefits

- The PRAFO® (Pressure Relief Ankle Foot Orthosis) is a totally adjustable, customfitted AFO that can help manage many ankle/foot anomalies, offering a range of adjustable features
- Custom-contoured aluminium heel connector bar helps to control dorsiplantar flexion by allowing infinite, measurable adjustments
- Ergonomically designed foot plate allows the liner to be changed easily without tools
- Integral, securely fastened vibram walking base provides a shoe-type surface for ambulatory patients
- The design protects the patient around the vulnerable areas of the foot and heel and will remain so during ambulation
- Calf and foot segments can be custommodified to accommodate patient physiology
- Anti-rotation bar made of malleable aluminium adjusts to control unwanted extremity rotation
- Suitable for day and night
- Universal sizing to reduce stock holding

Indications for Use

- · Pressure ulcers
- Can be used to resist the development of contracture

Note: For larger sizing see Page 131, 132 & 133.

Ordering Information

Size	Product Code
Universal	PRFO-550



Paediatric A-Boot

Features and Benefits

- Washable, soft fleece-lined cover wicks moisture away from the skin to maintain skin integrity and allows heel to float, thus eliminating pressure of friction on the heel and enhancing blood circulation
- Rubber sole with non-slip surface provides a walking base so patients can walk comfortably
- The splint frame resists plantar flexion while providing a dynamic push toe extension for the correction of foot drop, foot ankle contractures and deformity
- The adjustable toe post relieves pressure on the toes and adjusts to accommodate a variation of foot sizes
- Rotation bar can be positioned to the side to control hip and leg rotation
- An extender kit is available for patients with increased calf and foot circumference and or oedema
- Universal sizing to reduce stock holding

Indications for Use

• For the patient who requires positioning but can also transfer and ambulate

Note: For larger sizing see Page 130.

Ordering Information

Size	Foot Length (cm)	Calf Circ. (cm)	Heel to Top of Back (cm)	Product Code
Universal	12.7 - 20.3	<17.9	17.8	A-BOOT-P

PAEDIATRIC



Features and Benefits

- The semi-rigid insert can be positioned to accommodate the angle of the ankle/ foot combines with a removable pull strap to initiate a therapeutic dorsiflexion stretch to help combat mild plantar flexion contracture with or without inversion/eversion
- The cover has a non-skid surface sewn into the bottom to improve grip during transfers
- Toe post to provide protection from shear forces during in bed usage



- · Plantar flexion contractures of the ankle
- To protect the heel from unwanted pressure while in bed

Note: For larger sizing see Page 127.



DynaPro® AFO

Features and Benefits

• Utilises Flex-E-Core technology for adjustment without the need for tools or heat

Leada Multi-Use Boot

Indications for Use

- · Plantar fasciitis
- · Achilles tendonitis
- Night splinting
- · Foot drop

Note: For larger sizing see Page 128.



Size	Length of Foot (cm)	Max Calf Circ. (cm)	Product Code
Child	9 - 13	30	DPAFO24C
Youth	13 - 15	36	DPAFO25Y

Size	Max Foot Length (cm)	Brace Height (cm)	Product Code
Paediatric	15	25	LMU-P

Ponseti AFO Sandals Pair



Features and Benefits

- These sandals are for use with the Ponseti AFO bar, regardless of the size of the bar (except preemie bars)
- With stiff feet, progress can easily be seen as the heels come down further in the sandal each week and their location can be seen through two holes in the back of the sandals
- The sandals attach to the bar by means of a 'Quick Clip' system which enables the sandals to easily be clipped to, and unclipped from the bar
- The footwear consists of a well-moulded rubber footplate and soft straps
- The foot is held securely to the footplate and there is no slippage
- Stock supplied in grey. Pink and blue are available to order

Indications for Use

• Treatment of CTEV using the Ponseti method

Styles Available



Note: Accessories are available to order on Pages 237 & 238.

Is grey not your colour?

Other sizes and colours are available to order on demand!



Grey sizes 0000 - 10 are kept in stock, other sizes and pink and blue options are available to order on demand. Delivery times may vary on these products.

Foot Length		Product Code	
(mm)	Grey	Pink	Blue
60 - 62	AFO-P6-0-PRW	-	-
63 - 66	AFO-P5-0-PRW	-	-
67 - 70	AFO-0000-PRW	AFO-0000-PRP	AFO-0000-PRB
71 - 75	AFO-000-PRW	AFO-000-PRP	AFO-000-PRB
76 - 80	AFO-00-PRW	AFO-00-PRP	AFO-00-PRB
81 - 87	AFO-0-PRW	AFO-0-PRP	AFO-0-PRB
88 - 92	AFO-1-PRW	AFO-1-PRP	AFO-1-PRB
93 - 102	AFO-2-PRW	AFO-2-PRP	AFO-2-PRB
103 - 112	AFO-3-PRW	AFO-3-PRP	AFO-3-PRB
113 - 122	AFO-4-PRW	AFO-4-PRP	AFO-4-PRB
123 - 132	AFO-5-PRW	AFO-5-PRP	AFO-5-PRB
133 - 142	AFO-6-PRW	AFO-6-PRP	AFO-6-PRB
143 - 152	AFO-7-PRW	-	-
153 - 162	AFO-8-PRW	-	-
163 - 172	AFO-9-PRW	-	-
173 - 182	AFO-10-PRW	-	-
183 - 192	AFO-11-PRW	-	-
193 - 202	AFO-12-PRW	-	-
	(mm) 60 - 62 63 - 66 67 - 70 71 - 75 76 - 80 81 - 87 88 - 92 93 - 102 103 - 112 113 - 122 123 - 132 133 - 142 143 - 152 153 - 162 163 - 172 173 - 182 183 - 192	(mm) Grey 60 - 62 AFO-P6-O-PRW 63 - 66 AFO-P5-O-PRW 67 - 70 AFO-0000-PRW 71 - 75 AFO-000-PRW 76 - 80 AFO-00-PRW 81 - 87 AFO-0-PRW 88 - 92 AFO-1-PRW 93 - 102 AFO-2-PRW 103 - 112 AFO-3-PRW 113 - 122 AFO-4-PRW 123 - 132 AFO-5-PRW 143 - 152 AFO-6-PRW 143 - 152 AFO-7-PRW 153 - 162 AFO-8-PRW 163 - 172 AFO-9-PRW 173 - 182 AFO-10-PRW 183 - 192 AFO-11-PRW	(mm) Grey Pink 60 - 62 AFO-P6-0-PRW - 63 - 66 AFO-P5-0-PRW - 67 - 70 AFO-0000-PRW AFO-0000-PRP 71 - 75 AFO-000-PRW AFO-000-PRP 76 - 80 AFO-00-PRW AFO-00-PRP 81 - 87 AFO-0-PRW AFO-0-PRP 88 - 92 AFO-1-PRW AFO-1-PRP 93 - 102 AFO-2-PRW AFO-2-PRP 103 - 112 AFO-3-PRW AFO-3-PRP 113 - 122 AFO-4-PRW AFO-4-PRP 123 - 132 AFO-5-PRW AFO-5-PRP 133 - 142 AFO-6-PRW AFO-6-PRP 143 - 152 AFO-7-PRW - 153 - 162 AFO-8-PRW - 163 - 172 AFO-9-PRW - 173 - 182 AFO-10-PRW - 183 - 192 AFO-11-PRW -



Ponseti AFOSandals Single

Features and Benefits

- These sandals are for use with the Ponseti AFO bar, regardless of the size of the bar (except preemie bars)
- Please note sandals are provided with the buckles on the inside
- With stiff feet, progress can easily be seen as the heels come down further in the sandal each week and their location can be seen through two holes in the back of the sandals
- The sandals attach to the bar by means of a 'Quick Clip' system which enables the sandals to easily be clipped to, and unclipped from the bar
- The footwear consists of a well-moulded rubber footplate and soft straps
- The foot is held securely to the footplate and there is no slippage

Styles Available

Grey	Pink	Blue

Note: Accessories are available to order on Pages 237 & 238.

Indications for Use

 Treatment of CTEV using the Ponseti method

Ordering Information

Pink			
Size	Foot Length (mm)	Product Code	
Size		Left	Right
0000	67 - 70	AFO-0000-LP	AFO-0000-RP
000	71 - 75	AFO-000-LP	AFO-000-RP
00	76 - 80	AFO-00-LP	AFO-00-RP
0	81 - 87	AFO-0-LP	AFO-0-RP
1	88 - 92	AFO-1-LP	AFO-1-RP
2	93 - 102	AFO-2-LP	AFO-2-RP
3	103 - 112	AFO-3-LP	AFO-3-RP
4	113 - 122	AFO-4-LP	AFO-4-RP
5	123 - 132	AFO-5-LP	AFO-5-RP
6	133 - 142	AFO-6-LP	AFO-6-RP

Ordering Information

Grey				
Cina	Fact Longth (mm)	Product Code		
Size	Foot Length (mm)	Left	Right	
0000	67 - 70	AFO-0000-LW	AFO-0000-RW	
000	71 - 75	AFO-000-LW	AFO-000-RW	
00	76 - 80	AFO-00-LW	AFO-00-RW	
0	81 - 87	AFO-0-LW	AFO-0-RW	
1	88 - 92	AFO-1-LW	AFO-1-RW	
2	93 - 102	AFO-2-LW	AFO-2-RW	
3	103 - 112	AFO-3-LW	AFO-3-RW	
4	113 - 122	AFO-4-LW	AFO-4-RW	
5	123 - 132	AFO-5-LW	AFO-5-RW	
6	133 - 142	AFO-6-LW	AFO-6-RW	
7	143 - 152	AFO-7-LW	AFO-7-RW	
8	153 - 162	AFO-8-LW	AFO-8-RW	
9	163 - 172	AFO-9-LW	AFO-9-RW	
10	173 - 182	AFO-10-LW	AFO-10-RW	
11	183 - 192	AFO-11-LW	AFO-11-RW	
12	193 - 202	AFO-12-LW	AFO-12-RW	

Blue			
Size	Foot Length (mm)	Product Code	
Size	Foot Length (IIIII)	Left	Right
0000	67 - 70	AFO-0000-LB	AFO-0000-RB
000	71 - 75	AFO-000-LB	AFO-000-RB
00	76 - 80	AFO-00-LB	AFO-00-RB
0	81 - 87	AFO-0-LB	AFO-0-RB
1	88 - 92	AFO-1-LB	AFO-1-RB
2	93 - 102	AFO-2-LB	AFO-2-RB
3	103 - 112	AFO-3-LB	AFO-3-RB
4	113 - 122	AFO-4-LB	AFO-4-RB
5	123 - 132	AFO-5-LB	AFO-5-RB
6	133 - 142	AFO-6-LB	AFO-6-RB

PAEDIATRIC

On Demand

Ponseti AFO Flexion Stop Sandals

Features and Benefits

- These sandals are for use with the Ponseti AFO bar, regardless of the size of the bar (except preemie bars)
- The sandals are used in conjunction with the standard Ponseti AFO 'Quick Clip' system
- As with the standard Ponseti AFO, the PFS holds the foot out at an angle to maintain the correction of the clubfoot
- The sandals must be worn with good fitting socks to prevent any sores or blisters
- For ease of use the sandals are supplied with the buckles on the inside

Indications for Use

 Treatment of CTEV using the Ponseti method

Note: Accessories are available to order on Pages 237 & 238.

Ordering Information

Foot Size Length (mm)		Product Code					
		Left	Right	Pair			
00	76 - 80	AFO-00-LW-PFS	AFO-00-RW-PFS	AFO-00-PRW-PFS			
0	81 - 87	AFO-0-LW-PFS	AFO-0-RW-PFS	AFO-0-PRW-PFS			
1	88 - 92	AFO-1-LW-PFS	AFO-1-RW-PFS	AFO-1-PRW-PFS			
2	93 - 102	AFO-2-LW-PFS	AFO-2-RW-PFS	AFO-2-PRW-PFS			
3	103 - 112	AFO-3-LW-PFS	AFO-3-RW-PFS	AFO-3-PRW-PFS			
4	113 - 122	AFO-4-LW-PFS	AFO-4-RW-PFS	AFO-4-PRW-PFS			
5	123 - 132	AFO-5-LW-PFS	AFO-5-RW-PFS	AFO-5-PRW-PFS			
6	133 - 142	AFO-6-LW-PFS	AFO-6-RW-PFS	AFO-6-PRW-PFS			
7	143 - 152	AFO-7-LW-PFS	AFO-7-RW-PFS	AFO-7-PRW-PFS			
8	153 - 162	AFO-8-LW-PFS	AFO-8-RW-PFS	AFO-8-PRW-PFS			
9	163 - 172	AFO-9-LW-PFS	AFO-9-RW-PFS	-			
10	173 - 182	AFO-10-LW-PFS	AFO-10-RW-PFS	-			
11	183 - 192	AFO-11-LW-PFS	AFO-11-RW-PFS	-			
12	193 - 202	AFO-12-LW-PFS	AFO-12-RW-PFS	-			



Ponseti AFO Mitchell Ponseti Bar

Features and Benefits

- Abduction bar for managing the prescribed width and orientation of the corrected clubfoot
- Patented 'Quick Clip' system allows for easy attachment and removal
- Facilitates quick and easy clothing changes, placing the child in a high chair or shopping trolley
- Can be cut down to accommodate less than 7 inches
- Prepackaged at 60° abduction (may be varied using a Phillips screwdriver)

Indications for Use

Treatment of CTEV using the Ponseti method

Extra Short Bar 10° 17.3 25 BAR-PO	N-XSB10
Extra Short Bar 15° 17.3 25 BAR-PO	N-XSB15
Short Bar 10° 20 30 BAR-PO	N-SB10
Short Bar 15° 20 30 BAR-PO	N-SB15
Long Bar 10° 23.5 37.8 BAR-PO	N-LB10
Long Bar 15° 23.5 37.8 BAR-PO	N-LB15

N N

Ponseti AFO DBar Europe



Features and Benefits

- Independent Dorsi Assist Springs
- Adjustable bar width
- External rotation in 10° increments
- Reinforced polycarbonate plastic
- Minimalistic design to reduce weight

Indications for Use

• Treatment of CTEV using the Ponseti method

Description	Size (inch)	Product Code
DBar Europe Mitchell Spring Assist Bar	6 - 9	69MDSA
DBar Europe Mitchell Adjustable Spring Assist Bar	9 - 14	914MDASA
DBar Europe Mitchell Spring Assist Bar	9 - 14	914MDSA



Ponseti AFOPressure Saddle

Features and Benefits

- Slim design
- Designed to locate the heel within the AFO sandal

Indications for Use

• Difficulty locating patients heel within the AFO sandal

Size	Product Code
Small	PS-SM
Large	PS-LG



Calzamedi Footwear



Features and Benefits

- Extended stiffeners on both medial & lateral sides provide improved stability of the foot
- Lightweight rubber soles contribute to a more energy efficient gait
- Hook and loop fastening or lace styles help children develop independent dressing skills
- Low opening style makes footwear easy to take on and off
- Extra padding around heel stiffeners to improve comfort
- Fully breathable lining reduces skin irritation
- Adaptable sole provides enhanced function
- · Leather uppers are durable and breathable
- Padded collars ensure a comfortable fit
- Supplied with removable anatomically shaped footbeds



Black Leather

Indications for Use

- Flexible flat feet
- Hypermobility
- Hindfoot instability
- Sensitive feet
- · Children who need extra support during gait



Blue Nubuck

Note: Please see footwear size chart on the next page (page 244-245).



Pink Nubuck

Ordering Information

Fastening Type		Product Code		
Lace		-4200-0101-		
Hook and Loop		-4201-0101-		
Pink Nubuck		Blue Nubuck		
Fastening Type	Product Code	Fastening Type	Product Code	
Lace	-4202-0203-	Lace	-4202-0202-	
Hook and Loop	-4203-0203-	Hook and Loop	-4203-0202-	

See footwear sizing chart on following pages for more sizing information.

Example: To order a pair of black leather lace up shoes in size 37 with a width of 7.7cm (4 width), the code would be 37-4200-0101-4.



Calzamedi Footwear Size Chart

				Wid	lth 4	Wid	lth 6
Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Foot Length (cm)	Width Circ. (cm)	Joint Circ. (cm)	Width Circ. (cm)	Joint Circ. (cm)
1	16	10.3	9.2 - 9.7	12.7	4.8	13.7	5.2
1.5	17	11	9.8 - 10.3	13.1	5	14.1	5.4
2.5	18	11.7	10.4 - 10.9	13.5	5.1	14.5	5.5
3	19	12.3	11 - 11.5	13.9	5.2	14.9	5.6
4	20	13	11.6 - 12.1	14.3	5.3	15.3	5.7
5	21	13.7	12.2 - 12.8	14.7	5.4	15.7	5.8
5.5	22	14.3	12.9 - 13.3	15.1	5.5	16.1	5.9
6	23	15	13.4 - 14	15.5	5.6	16.5	6
7	24	15.7	14.1 - 14.7	15.9	5.7	16.9	6.1
8	25	16.3	14.8 - 15.2	16.3	5.8	17.3	6.2
8.5	26	17	15.3 - 15.9	16.7	6	17.7	6.4
9.5	27	17.7	16 - 16.6	17.1	6.2	18.1	6.6
10	28	18.3	16.7 - 17.2	17.5	6.4	18.5	6.8
11	29	19	17.3 - 17.9	17.9	6.6	18.9	7
12	30	19.7	18 - 18.5	18.3	6.8	19.3	7.2
12.5	31	20.3	18.6 - 19.1	18.7	6.9	19.7	7.3
13.5	32	21	19.2 - 19.8	19.1	7	20.1	7.4
1	33	21.7	19.9 - 20.5	19.5	7.1	20.5	7.5
2	34	22.3	20.6 - 21.1	19.9	7.3	20.9	7.7
2.5	35	23	21.2 - 21.8	20.3	7.4	21.3	7.8
3.5	36	23.7	21.9 - 22.5	20.7	7.6	21.7	8
4	37	24.3	22.6 - 23.1	21.1	7.7	22.1	8.1
5	38	25	23.2 - 23.8	21.5	7.9	22.5	8.3
6	39	25.7	23.9 - 24.5	21.9	8.1	22.9	8.5
6.5	40	26.3	24.6 - 25.1	22.3	8.3	23.3	8.7

W	Width 8 Anti Varus DAFO						
Width Circ. (cm)	Joint Circ. (cm)	Width Circ. (cm)	Joint Circ. (cm)	Width Circ. (cm)	Joint Circ. (cm)	Shoe Size (UK)	Shoe Size (EU)
14.7	5.6	-	-	-	-	1	16
15.1	5.8	-	-	-	-	1.5	17
15.5	5.9	-	-	-	-	2.5	18
15.9	6	-	-	-	-	3	19
16.3	6.1	15.1	5.6	17.6	6.2	4	20
16.7	6.2	15.5	5.8	18	6.3	5	21
17.1	6.3	15.9	5.9	18.4	6.4	5.5	22
17.5	6.4	16.4	6.1	18.8	6.5	6	23
17.9	6.5	16.8	6.2	19.2	6.6	7	24
18.3	6.6	17.2	6.4	19.6	6.8	8	25
18.7	6.8	17.6	6.5	20	6.9	8.5	26
19.1	7	18	6.7	20.3	7	9.5	27
19.5	7.2	18.4	6.9	20.7	7.1	10	28
19.9	7.4	18.9	7	21.1	7.3	11	29
20.3	7.6	19.3	7.2	21.5	7.4	12	30
20.7	7.7	19.7	7.3	21.9	7.5	12.5	31
21.1	7.8	20.1	7.5	22.3	7.6	13.5	32
21.5	7.9	20.5	7.6	22.7	7.8	1	33
21.9	8.1	20.9	7.8	23.1	7.9	2	34
22.3	8.2	21.3	7.9	23.5	8	2.5	35
22.7	8.4	21.8	8.1	23.9	8.1	3.5	36
23.1	8.5	22.2	8.2	24.3	8.2	4	37
23.5	8.7	22.6	8.4	24.6	8.4	5	38
23.9	8.9	23	8.6	25	8.5	6	39
24.3	9.1	23.4	8.7	25.4	8.6	6.5	40

HEAD & NECK

Control Tech Soft 3/4 Kids



Features and Benefits

 Rearfoot wedging reduces the need for excessive arch height whilst maintaining support, comfort and providing superior shock absorption

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

Note: For larger sizing, see page 190.

Ordering Information

Product Code
SC3051
SC3052
SC3053
SC3054



Trio Easy Fit

Features and Benefits

- The Easy Fit is intrinsically a ¾ length medially posted orthotic with an additional posting option for triplanar control of pronation featuring a lateral cutaway design
- This model maintains the essential features for pronation control, including an optional biplanar rearfoot post
- A medial forefoot post can be used to replace the first metatarsal cushion if required
- Comfort is achieved without compromising natural movement and comfort

Indications for Use

- · Medium arches
- Excessive pronation
- Forefoot pain
- Knee pain
- Heel valgus

Note: For larger sizing, see page 184.

Ordering Information

Size	Shoe Size (UK)	Product Code
Kids	1 - 2.5	TEF-K

Pair

Trio

Features and Benefits

- The Trio orthotic is designed to resist excessive valgus rotation, by inverting the calcaneus within the heel cup thereby maintaining mid tarsal joint. It also supports the proximal medial arch of the foot providing navicular stability
- Each pair is supplied with additional 3° medial biplanar rearfoot posts and 3° lateral forefoot posts
- Available in ¾ and full length

Indications for Use

- · Plantar fasciitis
- Heel pain
- Calf strain
- Anterior medial knee pain
- Hip instability

Note: For larger sizing, see page 186.

Ordering Information

34 Length							
Size	Shoe Size	Product Code					
Size	(UK)	Medium	Firm	Extra Firm			
Kids Small	12.5 - 13.5	TSLM-KS	TSLF-KS	TSLE-KS			
Kids Large	1 - 2.5	TSLM-KL	TSLF-KL	TSLE-KL			

Full Length						
Size	Shoe Size	Product Code				
Size	(UK)	Medium	Firm	Extra Firm		
Kids Small	12.5 - 13.5	TFLM-KS	TFLF-KS	TFLE-KS		
Kids Large	1 - 2.5	TFLM-KL	TFLF-KL	TFLE-KL		



Prostep Paediatric

Features and Benefits

- Prostep is designed to restore the natural position of the talus by supporting the inside edge of the heel and supports the spring ligament. This resists the valgus rotation of the rearfoot, which is associated with pronation
- Anti-bacterial top cover and forefoot and rearfoot integral posts of 0°, 3° and 5°
- Full length or ¾ length available

Indications for Use

- Low arches
- Excessive pronation
- Forefoot pain
- Knee pain
- Heel valgus

Note: For larger sizing, see page 188.

Ordering Information

3/4 Length						
Size	Chan Cina (UIC)	Product Code				
Size	Shoe Size (UK)	3/4 Medium	3/4 Firm			
Kids Small	13.5 - 1.5	PSSLM-KS	PSSLF-KS			
Kids Large	2 - 3	PSSLM-KL	PSSLF-KL			

Full Length				
Shoe Size (UK)				
3110e 312e (OK)	Full Length Medium	Full Length Firm		
13.5 - 1.5	PSFLM-KS	PSFLF-KS		
2 - 3	PSFLM-KL	PSFLF-KL		
		13.5 - 1.5 PSFLM-KS		



Kids Feet In Motion

Features and Benefits

- These brightly coloured pairs consist of one blue (left) and one right (red) to help children learn whilst providing a functional solution
- Medium density EVA orthosis specifically designed for children's feet
- A flat base on the standard orthotic allows for modification with wedges and components
- Insoles are available as standard or with a pre-posted heel of 3° or 5°

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Planovalgus feet
- Postural related symptoms

EZIKID (Peddy) Insoles

Features and Benefits

- Dynamic insole for children
- Modular sandwich-like construction and flexible thermoformable wedge
- Soft sole to support the natural movement of the foot

Indications for Use

- Back pain
- Posture related symptoms
- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- · Planovalgus feet

Ordering Information

Shoe Size	Shoe Size Shoe Size		Product Code		
(UK)	(EU)	Standard	3°	5°	
4	19 - 20	SC3304	SC3330	SC3350	
5	21 - 22	SC3305	SC3331	SC3351	
6.5	23 - 24	SC3306	SC3332	SC3352	
8	25 - 26	SC3307	SC3333	SC3353	
9 - 10	27 - 28	SC3308	SC3334	SC3354	
11	29 - 30	SC3309	SC3335	SC3355	
12 - 13	31 - 32	SC3310	SC3336	SC3356	
1 - 2	33 - 34	SC3311	SC3337	SC3357	

Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Product Code
4	19 - 20	ZKR1-0640-19/20
5	21 - 22	ZKR1-0640-21/22
6.5	23 - 24	ZKR1-0640-23/24
8	25 - 26	ZKR1-0640-25/26
9 - 10	27 - 28	ZKR1-0640-27/28
11	29 - 30	ZKR1-0640-29/30
12 - 13	31 - 32	ZKR1-0640-31/32
1 - 2	33 - 34	ZKR1-0640-33/34
2 - 3	35 - 36	ZKR1-0640-35/36
4 - 5	37 - 38	ZKR1-0640-37/38
6 - 7	39 - 40	ZKR1-0640-39/40



Podo Tech Junior 3/4 Length

Features and Benefits

- Designed especially for children
- Easily ground and heat mouldable

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Product Code
4	19 - 20	SC3220
5	21 - 22	SC3222
6.5	23 - 24	SC3224
8	25 - 26	SC3226
9 - 10	27 - 28	SC3228
11	29 - 30	SC3230
12 - 13	31 - 32	SC3232
1 - 2	33 - 34	SC3234
2 - 3	35 - 36	SC3236
4 - 5	37 - 38	SC3238
6 - 7	39 - 40	SC3240
8 - 9	41 - 42	SC3242



Pedipod Full Length Insoles

Features and Benefits

PediPod insoles are uniquely designed to help children develop healthy foot and body posture as they grow. The development of poor posture while growing potentially disturbs normal walking, running and standing which can lead to recurrent back, knee and foot pain as an adult - often resulting in a lifetime of treatment.

- Anatomical support features in insoles built for children's feet
- Provides comfort and pain relief by deep heel cup and paediatric orthotic insole specific arch support
- The brushed nylon surface is comfortable and allows easy passage of the foot into the shoe
- A full-length insole that can be quickly and easily trimmed to fit most comfort children shoes

Indications for Use

- Aching arches on exercise
- Children's heel pain during and after exercise (Sever's Disease)
- Ankle, calf or leg aches on exercise
- Children's kneecap pain during or after sport or on stairs

Ordering Information

Size	Shoe Size (UK)	Product Code
С	12 - 1	PED03
D	10 - 11.5	PED04
E	8 - 9.5	PED05
F	6 - 7.5	PED06

Pair

Peapod Junior

Features and Benefits

- 1.5mm thermoplastic base offers excellent fitting compatibility in most footwear
- Heat mouldable and grindable

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms



Kraemer Heel Cups

Features and Benefits

- Children's multi-layered prefabricated foot orthosis covered in Multiform A30 shore
- Thermo-mouldable reinforced base and wedge

Indications for Use


- Ligament weakness
- Sprains/strains
- Planovalgus feet

Ordering Information

SC2502
SC2504
SC2506
SC2508
SC2510
SC2512
SC2514
SC2516

Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Product Code
4	19 - 20	ZKR1-0607-19/20
5	21 - 22	ZKR1-0607-21/22
6.5	23 - 24	ZKR1-0607-23/24
8	25 - 26	ZKR1-0607-25/26
9 - 10	27 - 28	ZKR1-0607-27/28
11	29 - 30	ZKR1-0607-29/30
12 - 13	31 - 32	ZKR1-0607-31/32
1 - 2	33 - 34	ZKR1-0607-33/34
2 - 3	35 - 36	ZKR1-0607-35/36
4 - 5	37 - 38	ZKR1-0607-37/38
6 - 7	39 - 40	ZKR1-0607-39/40



Index

Symbols		ComfySplints™ Spring Loaded Goniometer Hand Orthosis	d 7
3/4 Length Silicone Insoles	177	Control Tech Flex 3/4	19
		Control Tech Soft 3/4 Kids	24
		Control Tech Soft ¾ Length	19
A		Control Tech Soft Full Length	18
	400	CTLSO Bivalve	4
A60™ Ankle Support	102		
A-Boot Universal	130	_	
AFO Light	115	D	
AirLift™ PTTD	104	Daniel Habrania I Tara	
Airsport™	103	Darco Universal Toe	19
Airtex Transpirable Ankle Brace	96	Alignment Splint DH Active Offloading	19
All Purpose Boot	142	Post-Op Shoe	14
Ankle Sleeve Plus	96	Diaped Duosoft Plus	18
Apollo Therapeutic Shoes	160	Diaped Duosoft Therapeutic	18
Arthritic Knee Sleeve	82	Dual Maternity Belt	10
Aspen Vista® Universal Collar	24	DynaPro® AFO	127, 23
		Dyliario Aro	127, 23
В		_	
		E	
Basic Foam Collar	16	Equiflex	18
BlueROCKER® 2.0	123	EVENup Shoe Lift	13
BlueROCKER® 2½	124	EZIKID (Peddy) Insoles	25
Bluerocker® Short	122	EZIKID (Peddy) Ilisoles	25
С		F	
Calzamedi Footwear	242	Felt Boots	14
Collamed® Knee Brace	242 86	Foot-Up®	10
Comfort Night Splint	125	Formfit® Ankle Brace	
Comfy Air Hand Orthosis	66	with Figure 8	10
•		Formfit® Shoulder Brace	5
Comfy Deviation Hand Orthosis	215	Formfit® Walker	13
Comfy Hand Thumb Orthosis	217	Full Length Silicone Insoles	
Comfy Hand-Wrist	7 216	with Metatarsal Support	17
Finger Orthosis 6 Comfyprene® Contour	7, 216		
Hand Orthosis	65		
ComfySplints™ Deviation	05	Н	
Hand Thumb Orthosis	69		
ComfySplints™ Large Pan		Heads First Collar	1
Hand Thumb Orthosis	71	Hernia Support Briefs Men's	. 7

Hernia Support Briefs Women's

Index

Hernia Truss with Bilateral Pads	77	Miami J® Select Collar	21
Hyperextension Brace	38	Mini TLC™ Hip Abduction Brace	224
I		N	
Interpod Diabetic	192	Neoprene Boots	140
		Neptun B Therapeutic Shoes	168
		Neptun Therapeutic Shoes	166
J		Neurodyn Foot Lift Orthosis	110
Jobskin® Arm Gaiters		Neuro-Lux Shoulder II®	-1
Deluxe Lining	220	(Sporlastic)	51
Jobskin® Arm Gaiters			
Standard Lining	218	0	
Jobskin® Leg Gaiters			
Deluxe Lining	226	OA Ease Knee Brace	90
Jobskin® Leg Gaiters	225	Oasis LSO	34
Standard Lining	225	Oasis TLSO	42
Jupiter Therapeutic Shoes	156	Ortel Hernia Truss	76
K		P	
Kids Feet In Motion	250	Paediatric A-Boot	231
Komet Therapeutic Shoes	144	Paediatric Hip Abduction	201
Kraemer Heel Cups	255	Splint - Adjustable	223
		Paediatric Hip Abduction	
		Splint - Standard	222
L		Peapod Junior	254
Leada Multi-Use Boot 128,	233	Pedipod Full Length Insoles	253
Leaf Spring Orthosis	114	PediSupport.night	126
Luna Therapeutic Shoes	162	PeroSupport.pro	117
		PeroSupport.tec	116
		Pluto S Therapeutic Shoes	170
M		Podo Tech Junior ¾ Length	252
		Podus Boot	129
Made To Measure Leg Gaiters	227	Ponseti AFO DBar Europe	240
Malleo Dynastab® Boa®	100	Ponseti AFO Flexion Stop Sandals	238
Ankle Brace	100	Ponseti AFO	230
Mars S Therapeutic Shoes Mars Therapeutic Shoes	164 152	Mitchell Ponseti Bar	239
Merkur A Therapeutic Shoes	150	Ponseti AFO Pressure Saddle	240
Merkur Therapeutic Shoes	146	Ponseti AFO Sandals Pair	234
Miami J® Cervical Collar	22	Ponseti AFO Sandals Single	236

 $\mathsf{ComfySplints}^{\scriptscriptstyle\mathsf{TM}}\,\mathsf{Rest}\;\mathsf{Hand}$

70

Index

PRIM Dorsi Assist	108
PRIM Dorsolumbar Back Support	36
PRIM MAX Medium	
Lumbosacral Back Support	32
PRIM SPINAIR SP500 TLSO	37
PRIMSPINE Semi-Rigid	
Lumbosacral Back Support	33
Prostep	188
Prostep Paediatric	249
Push Aequi Ankle Brace	99
Push Aequi Junior Ankle Brace	228
Push AFO	111
Push Care Ankle Brace	97
Push Care Wrist Brace	59
Push Med Ankle Brace	98
Push Med Wrist Brace	60
Push Ortho Thumb Brace CMC	58

R

Rebound® Air Walker	135
Rebound® Diabetic Walker	136
Rebound® Post Op Elbow	52
Rebound® Post-Op Knee Brace	87
Rebound Foot-Up®	105
Repose Wedge	128

S

Saturn Therapeutic Shoes	158
Shaped Foam Collar	17
Short Wrist Brace	61
Silicone Heel Cups	176
Sirius Therapeutic Shoes	172
Slimflex Green Insoles	
Full Length	195
Slimflex Simple Kit ¾ Length	194
Slim Gel-Fit™ Bunion Sleeve	179
SoftPro® Functional	
Resting WHFO	68
Stabilised Knee Sleeve	83
Stabimed® Knee Brace	85

Standard Knee Brace		84
Standard Maternity Support		30
Starlight Quilted Soft	13,	211
Starlight Secure Fabric	14,	212
Starlight Secure Leather	15,	213
Starlight Standard	12,	210
Swedish Type Lightweight AF	С	112
Swedish Type Lightweight		
AFO - Carbon Reinforced		113

Т

errapin bivaive TLSO Brace	43
Terrapin LSO Bivalve Brace	35
ToeOFF® 2.0	119
ToeOFF® Flow 2½	121
ToeOFF® Short	118
Tred-lite Ezi-Rize Heel Raise	182
<u> Frio</u>	186, 248
Trio Easy Fit	184, 247
Trio Heel Pain	183
Trio Slimfit	185

229

U UFO™

Universal Abdominal Binder	28
Universal Deep Abdominal Bir	ider <mark>28</mark>
Universal Palm Protector	56, 215
Universal PRAFO 650 Kodel	131
Universal PRAFO APU	
(Adjustable Posterior Upright)	132
Universal PRAFO Dual Action	133
Universal PRAFO Paediatric	230
Universal PROGLIDE 172	19
Universal PROGLIDE 174	20
Universal Temporary	
Hernia Support	29
Universal Wrist Brace	62
Unloader® One X	92
Uranus Therapeutic Shoes	154

Index

V

VACOCAST DIADETIC	138
VACOped	137
Vectorthotic Full Length Insoles	196
VENUM Universal HYBRID TLSO	39
VENUM Universal Soft Hip	
Abduction Brace	80
VENUM Universal TLSO PPX	40
Venus Therapeutic Shoes	148

V

Wrist and Thumb Spica	63
Wrist and Thumb Spica -	
Neoprene	64

X

K-Line 70 Insoles	199
K-Line Achilles Tendonitis Insoles	202
K-Line Additional Wedges	206
K-Line DJD Insoles	204
K-Line Plantar Fasciitis Insoles	205
K-Line Pressure Perfect Insoles	200
K-Line Reinforced ¾	198
K-Line RIF	
ReInForcement) Insoles	201
K-Line Standard Insoles	197
K-Line TPD Insoles Firm Density	203
K-ROM™ Post-Op Knee Brace	88

Notes	Notes



Steeper Group Unit 3, Stourton Link Intermezzo Drive Leeds LS10 1DF United Kingdom

Tel: +44 (0) 113 207 0435 Email: orthotics@steepergroup.com

www.steepergroup.com

STPOR018 Issue 4 Mar 2024

