

Cable and Spring Operated Hands



Cable and Spring Operated Hands

Mechanical hands designed to provide a functional terminal device with a natural cosmetic appearance.

Manufactured from high-strength aerospace alloy and housed inside a flexible PVC shell, the hands can be operated by body-powered cable control or manual operation.





Cable and Spring Operated Hands

- Four sizes available from child to adult
- Voluntary opening with choice of front pull, back pull or spring close options
- Choice of three cable terminations supplied with cable operated hands
- Available with plain hand base, M12 x 1.5 or ½ x 20 TPI threaded stem
- Compatible with our range of full-length Elegance™ cosmetic gloves in either PVC or silicone. Available in 19 skin shades with micropigmentation TrueFinish™ for a more lifelike appearance

Ordering Information - Spring Operated Hands

Part No. Plain Base	Part No. M12 x 1.5 Stud	Part No. ½ x 20 TPI Stud	Size	Left / Right	Application	Wrist Diameter	TG = TrueFinish PVC TSG = TrueFinish Silicone	
							Male Glove	Female Glove
SPA675L	SPA625LEU	SPA625LNA	6 ¼	Left	Child (Approx.8-12)	45mm	310	310
SPA675R	SPA625REU	SPA625RNA	6 ¼	Right	Child (Approx.8-12)	45mm	311	311
SPA725L	SPA725LEU	SPA725LNA	7 ¼	Left	Adolescent / Adult	50mm	306	304
SPA725R	SPA725REU	SPA725RNA	7 ¼	Right	Adolescent / Adult	50mm	307	305
SPA775L	SPA775LEU	SPA775LNA	7 ¾	Left	Adult	54mm	302	312
SPA775R	SPA775REU	SPA775RNA	7 ¾	Right	Adult	54mm	303	313
SPA825L	SPA825LEU	SPA825LNA	8 ¼	Left	Adult	54mm	300	
SPA825R	SPA825REU	SPA825RNA	8 ¼	Right	Adult	54mm	301	

Ordering Information - Cable Operated Hands - Front Pull Voluntary Opening

Part No. Plain Base	Part No. M12 x 1.5 Stud	Part No. ½ x 20 TPI Stud	Size	Left / Right	Application	Wrist Diameter	TG = TrueFinish PVC TSG = TrueFinish Silicone	
							Male Glove	Female Glove
MHA625LF	MHA625LFEU	MHA625LFNA	6 ¼	Left	Child (Approx.8-12)	45mm	310	310
MHA625RF	MHA625RFEU	MHA625RFNA	6 ¼	Right	Child (Approx.8-12)	45mm	311	311
MHA725LF	MHA725LFEU	MHA725LFNA	7 ¼	Left	Adolescent / Adult	50mm	306	304
MHA725RF	MHA725RFEU	MHA725RFNA	7 ¼	Right	Adolescent / Adult	50mm	307	305
MHA775LF	MHA775LFEU	MHA775LFNA	7 ¾	Left	Adult	54mm	302	312
MHA775RF	MHA775RFEU	MHA775RFNA	7 ¾	Right	Adult	54mm	303	313
MHA825LF	MHA825LFEU	MHA825LFNA	8 ¼	Left	Adult	54mm	300	
MHA825RF	MHA825RFEU	MHA825RFNA	8 ¼	Right	Adult	54mm	301	

Ordering Information - Cable Operated Hands - Back Pull Voluntary Opening

Part No. Plain Base	Part No. M12 x 1.5 Stud	Part No. ½ x 20 TPI Stud	Size	Left / Right	Application	Wrist Diameter	TG = TrueFinish PVC TSG = TrueFinish Silicone	
							Male Glove	Female Glove
MHA625LB	MHA625LBEU	MHA625LBNA	6 ¼	Left	Child (Approx.8-12)	45mm	310	310
MHA625RB	MHA625RBEU	MHA625RBNA	6 ¼	Right	Child (Approx.8-12)	45mm	311	311
MHA725LB	MHA725LBEU	MHA725LBNA	7 ¼	Left	Adolescent / Adult	50mm	306	304
MHA725RB	MHA725RBEU	MHA725RBNA	7 ¼	Right	Adolescent / Adult	50mm	307	305
MHA775LB	MHA775LBEU	MHA775LBNA	7 ¾	Left	Adult	54mm	302	312
MHA775RB	MHA775RBEU	MHA775RBNA	7 ¾	Right	Adult	54mm	303	313
MHA825LB	MHA825LBEU	MHA825LBNA	8 ¼	Left	Adult	54mm	300	
MHA825RB	MHA825RBEU	MHA825RBNA	8 ¼	Right	Adult	54mm	301	

To order the appropriate glove, prefix the number shown in the table with TG for TrueFinish PVC or TSG for TrueFinish Silicone. Then add the preferred skin shade. For example **TSG300/E4** orders a 8¼ left adult silicone glove in TrueFinish Color E4. To order a Color swatch use part number **E24655**.

3619 Paesanos Parkway, Suite 200,
San Antonio, TX 78231
Tel: 210 481 4126
Fax: 210 481 4117
Email: sales@steeperusa.com
www.steeperusa.com