

19 TrueFinish™ Colour Shades



Elegance Extra Length Silicone and PVC Cosmetic Gloves

- Extra Length as standard in 7 styles
- Recommended for SteeperUSA, Otto Bock, Motion Control and Hosmer adult and youths, myo electric and mechanical hand ranges
- A choice between the Standard colour range or the advanced micropigmentation TrueFinish™ for added realism and depth

Available in 19 colour shades, 5 sizes and 7 styles



Elegance Silicone

By using Tufsil 3™ our innovative polymer, Elegance Silicone gloves offer the toughness associated with PVC but provide the ease of cleaning provided by silicones.

In addition, Tufsil 3™ allows additional stretch, reducing the energy expended by the hand motor each time it opens. Properties of this polymer include improved tear resistance proven by clinical evaluation and up to 15% stretch, to allow even bulbous sockets to be accommodated.

All Elegance Silicone gloves are easy to fit and remove, allowing quick access to the hardware for maintenance and adjustment. The Silicone gloves are simple to keep clean using soap and water.

Elegance PVC

Elegance PVC gloves provide a cost effective solution for electric and mechanical hands when the extra toughness of PVC is required.

Size	Approx Length	Application	Left Silicone	Right Silicone	Left PVC	Right PVC	Colour Ranges
Standard Colours							
6 ¾	460mm	Youths	SG310	SG311	CG310	CG311	EO - E18
7	520mm	Female / Youths	SG308	SG309	CG308	CG309	EO - E18
7 ¼	645mm	Youths	SG306	SG307	CG306	CG307	EO - E18
7 ½	645mm	Female	SG304	SG305	CG304	CG305	EO - E18
7 ¾	620mm	Female	SG312	SG313	CG312	CG313	EO - E18
7 ¾	645mm	Male	SG302	SG303	CG302	CG303	EO - E18
8 ¼	645mm	Male	SG300	SG301	CG300	CG301	EO - E18
TrueFinish™ Colours							
6 ¾	460mm	Youths	TSG310	TSG311	TG310	TG311	EO - E18
7	520mm	Female / Youths	TSG308	TSG309	TG308	TG309	EO - E18
7 ¼	645mm	Youths	TSG306	TSG307	TG306	TG307	EO - E18
7 ½	645mm	Female	TSG304	TSG305	TG304	TG305	EO - E18
7 ¾	620mm	Female	TSG312	TSG313	TG312	TG313	EO - E18
7 ¾	645mm	Male	TSG302	TSG303	TG302	TG303	EO - E18
8 ¼	645mm	Male	TSG300	TSG301	TG300	TG301	EO - E18

19 TrueFinish™ Colour Shades



Ordering Examples

TSG300/E5
Orders a 8 ¼ Left Silicone Glove in TrueFinish colour E5

CG300/E5
Orders a 8 ¼ Left PVC Glove in Standard Finish colour E5

Swatches are available on request:
E24654 EO - E18 Standard colours
E24655 EO - E18 TrueFinish colours

3619 Paesanos Parkway, Suite 200,
San Antonio, TX 78231
Tel: 210 481 4126
Fax: 210 481 4117
Email: sales@steeperusa.com
www.steeperusa.com